# Bure Rise & Go DB

SystemRoMedic™

Instructions for use – English

Bruksanvisning – Svenska

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Max. patient weight: 150 kg / 330 lbs

## Table of contents

English .....	3
Svenska .....	17
Dansk .....	31
Suomi .....	45
Français .....	59
Lietuvių .....	73
čeština .....	87
<b>EMC</b> .....	<b>101</b>
<b>Parts List</b> .....	<b>104</b>
<b>Product Overview</b> .....	<b>106</b>

# Bure Rise & Go DB

## SystemRoMedic™

### Instructions for use - English

The Bure Rise & Go Double is a conventional, electrically operated Walker that is enhanced with a power rise function. Not only does it act as a Walker, it also makes getting about in everyday life easier thanks to its smart power rise function.

The Bure Rise & Go can be used by anyone who regularly uses a Walker. The user can remain active while the power rise function is in use.

The product consists of a Walker and its associated harness and straps. The harness is available in different sizes (S, M, L, XL) and is split at the rear to make remaining in the harness possible while using the WC. The harness is secured to the two attachment straps (one each side), which are in turn attached to the Walker.

When the patient stands up the electric power rise function is employed and the smart attachment reduces the distance between the patient and the Walker until the patient is standing upright.

The Walker can be used to advantage with a platform to facilitate patient transportation.

## 1. Intended use

Walkers are indoor, movable and portable assistive aids, which together with approved accessories are intended to be used for transfer from sitting to standing positions, and walking training. Walkers are to be used by persons with reduced muscular strength, who shall be able to stand and perform simple leg movements. Walkers are to be used for leaning on when walking and standing during early mobilization training. Walkers and accessories are intended for indoor use only in Home Healthcare environment and Professional Healthcare facility environment.

## 2. Contraindications, Precautions and Warnings

<b>Contraindications</b>
The device may not be used by patients above the maximum weight indicated on the label of the device.
<b>Precautions</b>
Check that the device is used on flat surface.
Check the device is correctly mounted/assembled before its first use.
Check the device after every folding after any transport.
Check height adjustment and leg spreading motions and inspect the actuators' or gas pistons' full ranges.
Activate the brakes whenever the device is not in use.
Check the function of the brakes before use.
Inspect the device regularly to detect any signs of damage. If there are any visible signs of damage, contact Direct Healthcare Group Customer Service.
In case of an electrical model of the device, check that the hand control does not show signs of wear before every use.
Check that the hand control markings are in accordance with the lifting functions.
Check the battery status of an electrically operated model, before using the device.
It is important to never leave the user alone during the transfer.
Warranty applies only if repairs or alterations are made by personnel who are authorized by Direct Healthcare Group.
Ensure that there are no obstacles or people in the way of the device, when moving.
Handle batteries with care. Do not drop.

Use only batteries and cables that are intended for the device as instructed by the manufacturer.
Verify that potential accessories that might be used are properly attached to the device before use.
Low speed is recommended when moving with the device.
Take care not to drive the device over thresholds with high speed or force. Approach the threshold with lower speed and communicate the coming threshold to the user.
Keep the instructions for use together with the product so that it can be used as a reference if necessary.
Beware of moving parts during usage.
Perform a risk assessment on the user before using the Walker to ensure they can stand and perform simple movements.

<b>Warnings</b>
The caregiver must be able to read and understand the Manual/IFU of the product.
It is important to use only approved accessories to prevent unintended detachment of components and subsequently a fall that may lead to patient injury.
Use careful and gentle maneuvers when moving the device.
Perform maintenance/service of the device, according to the instructions in the manual/IFU, at least once every 12 months.
Accessories must be properly attached. They should be tested in relation to the user's needs and functional ability.
Special care must be taken when using strong electrical power sources such as diathermy and the like so that diathermy cables are not placed on or near the device. In case of doubt consult with a Direct Healthcare Group representative.
When lifting the device, only use the denoted handholds. At least two people are required for lifting the device.
Do not leave a care recipient unattended during a transfer.
Never move the lift by pulling on the actuator or gas piston.
The device must not be lowered into water.
The device must not be left or stored in a damp or humid environment.
The device must not be cleaned using steam.
The device must not be charged in a wet room.
The device must not be used outdoors, only indoors.
To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
No modification of this equipment is allowed
Do not modify this equipment without authorisation of the manufacturer
If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of this equipment.

**Before use:**

- Inspect the packaging for any damage.
- Check for correct product being delivered.
- Check to ensure that the device is free from damage.
- Check if all the parts/components are included in the packaging.
- Check if all approved accessories are included in the package.
- Check if Product Quality Approval document is included in the packaging. Save this document for future contact with manufacturer.
- Make certain the device is properly assembled.
- Check lifting function and bottom-frame adjustment function.
- Check if it is possible to move the device in all directions.

### 3. Standing Up Using the Bure Rise & Go

Because standing up or transporting patients always entails risks we recommend that two people are always present to assist the user at such times.

1. Fit the harness as low as possible across the hips and bottom with the broad upper part around the hips and the narrow lower part around the bottom. Attach the strap with the Velcro fastener and tighten it using the steplessly adjustable safety belt. Note that the lower part of the harness attached around the patient's bottom can be opened as required, e.g. for visits to the WC.

2. Place the user's legs at a 90° angle with their feet on the floor or platform. Adjust the shin supports to fit the user's shins; remember to avoid pressure on the patellae and patellar tendons. Refer to page 5 for shin support adjustment.

3. Thread both attachment belts through the Walker's frame and the green locking devices. Make sure the belts are threaded completely through the locking devices. NB: Check especially that the belts are undamaged. Check that belts are not worn or frayed and that the holder and straps are NOT bent or crooked. If the holder and belts are crooked (pic 3c-red line) it is a sign that the product has been overloaded, in which case the frame must be replaced before further use.

4. Connect the power rise harness and attachment belts (that are secured to the Walker) and tighten them using the loops; the belts must be taught on both sides. (see fig. 4a)

Note: Slide the loop back once the straps have been tightened – this prevents tangling (see fig. 4b).

1

2

3a

3b

3c

4a

4b

5. Swing the forearm supports to the side; the user must now grasp the handles WITHOUT support for the fore-arms.
6. Lock at least two of the Walker's castors (any two) and fold down the stirrup.
7. Perform the standing manoeuvre and check that the harness pulls evenly on both sides; if not, adjust immediately. During the standing manoeuvre one carer must stand in front of the user while applying their weight to the stirrup and maintaining eye contact with the patient. The other carer may stand next to the patient.
8. In most cases the patient will not achieve a fully upright standing position; encourage the patient to achieve a final upright posture independently or assist the patient by tightening the adjustable harness straps.
9. Lower the Walker to its individually adjusted height and swing the forearm supports back into position to provide support for the user. If a platform is not being used, swing the shin supports away for increased walking room: refer to page 5, Shin support adjustment. The harness may remain attached for enhanced safety, e.g. during transportation using the platform.

## 4. Using Bure Walker with Gait training Kit 56-388 - KIT

The walking training harness is an accessory used for increased safety when walking in the product Bure Rise & Go. The product's both straps and fastening hooks minimize risks if users suffer from dizziness or similar issues.

Step 1: Applying the Gait training kit on the patient

The walking training harness is to be used together with Belt MULTI in various sizes. The Gait training kit contains 4 parts made of soft materials: 2 straps with soft padding and 2 straps with fastening hooks.

- Straps with padding are to be applied between the patient's legs and locked in the intended lock.
- Straps with fastening hook hangs loose until patient is secured in Bure Rise & Go (see later step 3)


Step 2: Secure attachment modules to the Bure Rise & Go Double

- The attachment modules (2 pcs) is pushed in from the inside on the left and right sides of the armrests under the cushions. We recommend that the modules are placed in the middle of the armrests.
- Lock the module with the supplied screw. Make sure that the module is secured.


Step 3: Securing a patient in Bure Rise & Go Double

If there is a need to further secure the patient, hook the carabiner hooks on the green bands (right and left side) in the fastening module. Tighten the fastening strap so that it is as short as possible.

Pre-caution:

- Green strap is threaded through the buckles from the inside / lower part (A) - up through the buckle & upper part (B).
- Loose straps must be outside the tensioned strap during use


## 5. Raising/Lowering the Frame

The desired height is adjusted with Arrow <UP> or Arrow <DOWN>. Adjust bottom frame using <WIDER> arrow or <NARROWER> arrow.

Handset from Plus also have charge indicator / battery (low, medium, high) and Service indicator (see further under "21. Hand Control Plus").

## 6. Adjusting the Handles

To adjust handles to the desired position; undo knob, adjust handle, tighten knob.

## 7. Left and Right Adjustment of Armrests

The armrests can be adjusted sideways for optimal support: the whole arm can be moved in or out and locked with the knob.

## 8. Locking/Unlocking Castors

Depress the foot brake to set it and lock the castor. To unlock, press down on the protruding part of the brake.

## 9. Adjusting the Shin Supports

a - knob for height adjustment.

b - knob for locking shin supports in forward position (when not in use). Supports it.

## 10. Control Box Wiring

Remove the control box cover to make the connections as described below. Use a screwdriver or similar; refer to the special labels on the control box.

NOTE: Make sure that plugs are connected to the correct sockets otherwise the actuators can be damaged.

1. Socket for the height adjustment actuator.
2. Socket for the width adjustment actuator.
3. Hand control socket.
4. Battery terminal socket.
5. Plug socket (NOTE: Plugs must ALWAYS be connected to the control box for IP classification to be valid).

## 11. Battery Charging

- Plug needs to connect to 230V wall socket. Charging starts automatically.
  - LED / battery flashes green.
- The battery is of the LiOn type and insensitive to how long / how often you charge the battery. Full charge on discharged battery takes about 24 hours.
- We recommend that the battery plug be removed from the control box if the product will be out of use for an extended period (more than one week) or if the battery is switched off by means of the ON/OFF switch. This eliminates the risk of disabling the battery as the control box always uses a small quantity of electricity even when the Walker is not in use.

## 12. LED on Control Box and Battery

Control box – green light when power is available from battery/mains socket.

Battery:

- The LED flashes green during charging.
- ON/OFF: To use the product, press ON.
- The green lamp will begin to flash.
- Use the OFF button (hold down for 5 s) to switch off the battery if the product will not be used for an extended period.
- The LED will show a constant green light when the battery is fully charged.
- When the battery is in normal state/in operation, the LED flashes green at REGULAR intervals, namely 0.5 s ON, 4 s OFF.
- When the battery requires charging (remaining charge <25%) the LED will flash orange/red at REGULAR intervals, namely 0.5 s ON, 4 s OFF.

NOTE: The battery has an in-built device that automatically switches it off when remaining charge is <20%. The LED is then extinguished. The battery must now be fully charged before the electronics can be used again.

## 13. Maintenance

Check that everything is clean and working before using.

Recommended at least once per year

Check that	Action
<p><b>General</b>                      The Walker feels firm/free from play.                      The Walker does not rattle during manoeuvres.                      The Walker is level and all of its castors are in contact with the floor.                      The Walker is not dirty.</p>	<p>Clean with lukewarm soap solution or alcohol-based cleaning agent (no petroleum products).</p>
<p><b>Armrests</b>                      Armrests are intact and clean.                      Armrest widening is functional.</p>	<p>Clean; replace if damaged.                      Replace lock knob.</p>
<p><b>Handle</b>                      Grips are not damaged/dirty.                      Handle adjustment is functional.</p>	<p>Clean with lukewarm soap solution (no petroleum products).                      Replace handle grips.                      Fit lock screws/plates as required.</p>
<p><b>Battery</b>                      Check that the charging cable is connected to the control box.                      Check that the hand control is connected to the control box.                      Connect or fit new hand control as required.                      Check that the battery, actuators and control box are not loose.</p>	<p>Fit new charging cable; this must always be connected to the control box (IP65 protection).                      Connect or fit new hand control as required.                      Refer to the wiring diagram in the control box wiring section.                      Tighten; replace with new fasteners as required.</p>
<p><b>Height adjustment</b>                      Raising and lowering are functional.                      The Walker feels stable at maximum settings.                      There is no play between the vertical frame and bottom frame.                      Height adjustment clamps lock.                      The electric actuator fitting must be free from play.                      Height adjustment via hand control is functional. The electric actuator must run smoothly at a constant speed.</p>	<p>Fit new guide sleeves (in frame) or end plugs (in chrome tubes) as required.                      Tighten all lock bolts.                      Replace clamps.                      Tighten the attachment fitting concerned; replace bolts/lock nuts as necessary.                      Check that the battery/hand control/ actuator are connected in accordance as in the control box wiring section.                      Charge the battery. Contact customer services.                      The part list shows the design and which fasteners must be checked.</p>
<p><b>Frame parts</b>                      There is no mechanical damage.                      There are no scratches.                      End plugs/lower frame fitted</p>	<p>If damage is present, contact DHG customer services.                      Touch up as necessary.                      Fit new end plugs.</p>
<p><b>Castors/brakes</b>                      Castors roll easily/tread not damaged.                      The castors are firmly fastened to the lower frame.                      Castor brakes functional on all castors.                      Handbrake functional.</p>	<p>Clean or replace castors (NOTE! Castors are always fixed to the lower frame with thread locker or lock nuts. The castors are sealed and we do not recommend their dismantling; instead, replace the whole castor assembly. Tighten the castor bolt and fit locking nut or use thread locker (depends on model).                      Replace castor assembly.                      Adjust the brake or fit new castors.</p>
<p><b>Harness</b>                      Belts, seams and Velcro fasteners are whole/intact.                      The Walker lockplate (green) is functional.</p>	<p>Replace harness. Available in three sizes (S, M, &amp; L)                      Replace lockplate.</p>
<p><b>Holder, belts</b>                      Holder, belts not bent (see Figure s3).                      Harness and attachment belts must be replaced every three years.</p>	<p>A bent holder is a sign that the product has been overloaded. The frame must be replaced before further use. The action is then "Check harness labeling. If the harness production date is older than 3 years, the harness and belts must be scrapped.</p>


DHG can provide upon request necessary technical information in order to assist maintenance and repair of the device.

## 14. General Maintenance Advice & Cleaning of Walker

- Before use, check that the holder and belts are not bent (see Figure s3).
- Before use, check that the harness is not more than 3 years old (see labeling).
- Clean the walker armrests and handles after every use. Wipe using a mild soap solution/surface disinfectant. However, the product is not corrosion proof which means washing by machine may reduce its lifetime.
- Replace damaged parts immediately; must be done by experienced personnel.
- The product is only intended for indoor use.
- With environments with high risk of infections use the walker with PU-cushions (Accessory number 57-366-L, 57-366-R and 57-365)

## 16. Storage and transportation

If the walker is not to be used for some time or e.g., during transport, we recommend that the emergency stop button be pressed in. The walker should be transported and stored in  $-10^{\circ}\text{C}$  to  $+50^{\circ}\text{C}$  and in normal humidity, 20% to 80% non-condensing. The air pressure should be between 700 and 1060 hPa. Leftmost symbol indicates storage and transportation. Let the walker reach room temperature before the batteries are charged or the walker is used. The walker should not be stored so that it is exposed to dust, or so that the battery or gas piston is exposed to direct sunlight.

## Operation


The operating environment should be  $5^{\circ}\text{C}$  to  $40^{\circ}\text{C}$ , relative humidity 20% to 80% non-condensing, and atmospheric pressure 700 to 1060 hPa. See also Technical Information below. Leftmost symbol indicates operating condition.

## 17. Technical Specifications

Dimensions apply to Walkers with 125mm castors.

Part No.	Product width	Max user weight	Length	Height	Weight	Armrest/ internal dimensions
56-309	66-106 cm	150 kg	101 cm	95-135 cm	40 kg	36-65 cm

## 18. Recycling/Disposal

When the product is no longer usable it should be recycled in accordance with legislation and regulations in the country concerned. All electrical parts including batteries must be removed and recycled as electrical components. (Ask your local recycling station for further information about how the different materials in the product (metals, plastics, electronics) should be recycled).

## 19. Warranty

For further information or enquiries please contact:

Direct Healthcare Group Sverige AB  
 Torshamnsgatan 35,  
 SE-164 40 Kista, Sweden

Tel: +46 (0)8-557 62 200  
 info.se@directhealthcaregroup.com

[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)

## 20. Electrical System - Technical Data

Actuator	Lift	Widening
Producer	TiMOTION	TiMOTION
Voltage	24V	24V
Push	2000N/1500N	1000N
Pull	2000N/1500N	1000N
Self-locking push	1000N	1000N
Current	6.6 (+/-1,0)A	2.6 (+/-0,4)A
Colour	Black	Aluminium
Stroke	400mm	55mm
Mounting dim	605mm	167mm
Encapsulation	IP66	UIP

IEC60601-1 3rd revision

## Control box + Battery

Producer	TiMOTION
Voltage	24V
Capacity	2.0AH
Encapsulation	IP66
Colour	BLACK
Electrical safety	ClassII

IEC60601-1 3rd revision

## 21. Hand Control Plus

\*\*1. Service indicator – flashes red after 11 months' operation.

Constant light (red) after 12 months' operation.

The indicator is a reminder to review the maintenance instruction (see Maintenance Instructions). This should be done at least once every 12 months.

Reset: hold down buttons <2a> + <2b> for at least 30 seconds.

\*\*The indicator will go out around 30 seconds after the hand control is activated (battery conservation function).


## 22. Expected lifetime:

The life expectancy of the product is 8 years.

## 23. Accessories

Part number	Description
56-339-4-0	Manual brake, 2 handles
56-339-4-1	Manual brake, 1 handle
56-333-DB	Standing plate, Short
56-333/SW-DB	Standing plate Short, SW, Bure DB, Rise DB
56-383-DB	Standing plate, Long
56-383/SW-DB	Standing plate Long, SW, Bure DB, Rise DB
56-335-3	Holder for Dripstand
56-336	Oxygen tube holder
56-337	Side support
56-338	Dripstand holder, Complete
56-387	Wheel w brake 75 mm
56-358	Wheel w brake 100 mm
56-361	Wheel w brake 125 mm
57-365	Armrestpad center, PUR
57-366-L/R	Armrestpad left/right, PUR
57-367	Kneepad, PUR
58-374-2	Hand Control (TH7)
58-375-X2	Hand Control PLUS, DB (TH10)
58-326 (-1, -2, -3)	Coiled Cable/EUR (-UK, -US, -AUS)
58-411	Hand Control TIMOTION
56-384-MULTI-S	Belt MULTI S, Rise&Go kit
56-384-MULTI-M	Belt MULTI M, Rise&Go kit
56-384-MULTI-L	Belt MULTI L, Rise&Go kit
56-384-MULTI-XL	Belt MULTI XL, Rise&Go kit
56-388-KIT	Kit/Gait training for Belt Mult

## 24. Table of Symbols

	This product complies with the requirements of the Medical Device Regulation 2017/745		Visual Inspection
	Medical Device		Read the manual
	Caution		Product Code
	Manufacturer information		Batch Code
	Type BF, according to the degree of protection against electric shock		The device is intended for indoor use
	WEEE Symbol May not be discarded in domestic waste		Class II Equipment
	Duty Cycle: 2 min in active (ON) mode. 18 min in rest (OFF) mode.		Wash at 75 degrees
	Do not iron		Do not dry clean
	Do not bleach		

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# Electromagnetic Compatibility (EMC)

## Electromagnetic emissions and test levels

The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Emission test	Standards	Compliance
RF emissions	CISPR 11	Group 1
RF emissions	CISPR 11	Class B
Harmonic current emissions	IEC 61000-3-2	Class A
Voltage fluctuations and Flicker emissions	IEC 61000-3-3	Complies

## Electromagnetic Immunity

Acceptance criteria for the EMC pass criteria is unintentional movement above  $\pm 10$ mm is not allowed. The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Acceptance criteria for the EMC pass criteria is unintentional movement above  $\pm 10$ mm is not allowed. The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Radiated Fields in close proximity, Immunity test and compliance. Standard: IEC 61000-4-39			
Dwell time: 3 sec			
Test levels (A/m)	Modulation Pulse modulation	Mod. Frequency (kHz)	Test Frequency
8	CW	CW	30 kHz
65	50% Duty cycle	2.1	134.2 kHz
7.5	50% Duty cycle	50	13.56 MHz

Immunity test	Basic EMC standard or test method	IEC 60601-1-2 Edition 4 Test levels and compliance	
		Professional healthcare facility environment†	Home healthcare environment
Surge	IEC 61000-4-5	± 0.5 kV, ± 1 kV	
Voltage dips	IEC 61000-4-11	0% U <sub>T</sub> ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
		0% U <sub>T</sub> ; 1 cycle and 70% U <sub>T</sub> ; 25/30 cycles. Single phase: at 0°	
Voltage Interruptions	IEC 61000-4-11	0% U <sub>T</sub> ; 250 / 300 cycle	
Conducted disturbances induced by RF fields	IEC 61000-4-6	6V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	6V in ISM and amateur bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
		Dwell time: 3 sec Frequency step size: 1%	
Radiated RF Electro-magnetic fields	IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz
		Dwell time: 3 sec Frequency step size: 1%	
Electrical fast transients / burst	IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	
RATED Power frequency magnetic field	IEC 61000-4-8	30 A/m 50 Hz or 60Hz	
Electrostatic Discharge (ESD)	IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV Number of discharges: 10 for each polarity	


Radiated Radio – Frequency, Proximity fields from wireless communication equipment & compliance. (Dwell time: 3 sec.)  
 Standard: IEC 61000-4-3

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test levels (V/m)
Spot Frequencies	385	380 to 390	TETRA 400	Pulse modulation 18 Hz	27
	450	430 to 470	FMRS 460, FRS 460	FM ± 5 kHz Deviation 1 kHz sine	28
	710	704 to 787	LTE Band 13, 17	Pulse modulation 217 Hz	9
	745				
	780				
	810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
	870				
	930				
	1720	1700 to 1990	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 217 Hz	28
	1845				
	1970				
	2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
	5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
	5500				
5785					
Other Identified frequencies	433	-	-	Pulse modulation 2 Hz	3
	-	860-960	-	Pulse modulation 2 Hz	54
	2450	-	-	Pulse modulation 2 Hz	54



## Parts List

No.	Art No.	Description	No.	Art No.	Description
1	58-309	Leg support vinyl connection	25	58-389	Battery box 4.0 Ah TBB5/Li-On
2	58-312 57-367	Leg support pad vinyl Rise&Go Leg support pad PU Rise&Go	26	58-377	Actuator Timotion 2000N/5200rpm Bure
3	57-887	Screw K6S M6x16	27	56-353	Bure Handle assy
4	58-306-L 57-366-L	Armrest Outer pad Left Rise&Go Armrest Outer pad PU Left Rise&Go	28	58-338	Bushing Vertical suport frame
5	58-306-R 57-366-R	Armrest Outer pad Right Rise&Go Armrest Outer pad PU Right Rise&Go	29	56-308	Cable holder Bure charger
6	58-305 57-365	Armrest Center Pad Rise&Go Armrest Center Pad PU Rise&Go	30	57-639	Screw K6S M6x8
7	58-351	Middle Frame Rise&Go	31	57-805	Washer BRB 6,4x12x1,6
8	-	Widening assembly Bure Double, R&G, see seperate view	32	57-817	Lock nut CL6 M10
9	58-311	Telescope inner pipe	33	57-825	Washer 10,5x18x2 FZB
10	58-363	Armrest frame center Part Rise&Go	34	57-848	Nordlock 10,7x21x2,5
11	58-374-X2	Hand control TiMotion TH7	35	57-787	Washer 10x25x1,5
12	-	Axis middle frame (Included in 58-351)	36	57-885	Screw MF6S M6x10
13	57-310	Middle frame bushing	37	57-883	Screw MF6S M6x40
14	58-355-L	Bottom Beam Left Rise&Go DB	38	57-879	Cap nut M6
15	58-355-R	Bottom Beam Right Rise&Go DB	39	57-878	Screw K6S M6x25
16	57-049 57-319 57-051	Castor wheel 125mm Brake Castor wheel 100mm Brake Castor wheel 75mm Brake	40	57-811	Serrated lock washer M16
17	57-775	Dome cover wheel	41	57-816	Screw MC6S M10x40 FZB
18	57-776	Dome cover wheel top	44	57-806	Washer 8.4x16x1.5
19	58-354-1	Wheel attachment	45	57-769	Washer 12,5x35x2 FZB
20	58-360	Cover Wheel attachment	46	57-766	Lock nut CL6 M12 FZB
21	57-041-3	Holder remote control TiMotion	47	58-322	Band washer green Rise & Go
22	57-369	Battery clamp Alu	48	58-304-L	Armrest Hinge Left Rise&Go
23	57-370	Battery bracket Alu	49	58-304-R	Armrest Hinge Right Rise&Go
24	58-373	Control box 220V Digital Li-On TiMotion Bure	50	57-880	Cap nut M8


No.	Art No.	Description	No.	Art No.	Description
51	57-867	Screw M8x90 MC6S FZB	63	57-859	Screw M6SF M8x50
52	57-720	Bushing Leg support bracket	64	57-733	Washer Nylon 30x8x7 Rise&Go
53	58-310-1	Leg support hanger	65	58-321	Bushing Armrest Rise&Go
54	58-308-1-L	Leg Support Bracket 14mm LEFT Rise&Go	66	57-968	Screw MC6S M10x55 FZB
55	58-329-PVC 80	Clip Green PE Soft Rise&Go	67	57-839	Star knob handle
56	58-308-1-R	Leg Support Bracket 14mm RIGHT Rise&Go	68	57-710	Star knob handle
57	58-302	Bushing PA6 Black	69	57-917	Wing nut M8
58	57-734	Washer 22x6,5x8	70	57-801	Lock nut M8 FZB
59	57-713	Roll Strap guide Rise&Go	71	57-864	Screw K6S M6x20 FZB
60	58-323	Lock bushing leg support	72	57-812	Tapping screw RXS 2,9x4,5
61	57-857	Screw K6S M4x10 FZB	73	58-326	Spiral Cord EUR - TiMotion, Bure
62	57-858	Washer BRB 8x4,3x0,8	74	57-881	Screw K6S DIN 7380 8.8 M8 x 30 FZB

No.	Art No.	Description
8		Widening assembly Bure Double, R&G

No.	Art No.	Description	No.	Art No.	Description
1	58-356	Housing Turn motor	14	57-826	Washer 8,4x29x1,5
2	58-352	Crosection bottom Bure Double	15	57-760	PVC Spacer 8x12x6,4
3	58-378	Actuator TA16 TIMotion	16	57-765	Screw MF6S M8x25
4	58-364	Turning wheel	17	58-302	Crosection bottom bushing
5	57-763	Ball joint angled M8	18	57-818	Screw K6S M6x10
6	58-381	Threaded steel rod M8	19	57-778	Washer Nylon 05-8-12-10
7	58-380-R	Turning tube right	20	57-774	Washer Nylon 13x34x1,5
8	58-380-L	Turning tube left	21	57-876	Nut M6M M8
9	57-801	Lock nut M8	22	57-882	Screw K6S 10.9 M8x55
10	57-806	Washer 8.4x16x1.5	23	57-767	Washer Starloc 8x15,3
11	58-358	Cover Housing Turn motor	24	57-771	Screw RXS 3,5x6,5
12	57-777	Washer nylon 05-8-12-3	25	57-875	Set screw SK6S M8x6
13	57-836	Spring lock washer M8			

# Product Overview

## REF

Article number	Product name
56-309	Bure Rise&Go DB Timotion, 125mm
56-309-AUS	Bure Rise&Go DB Timotion, 125mm, AUS
56-309-B/0-PU	Bure Rise&Go DB Timotion, 125mm, handbrakes, PU
56-309-PLUS	Bure Rise&Go DB Timotion, 125mm, Plus
56-309-PLUS-PU	Bure Rise&Go DB Timotion, 125mm, PU, Plus
56-309-PLUS-PU-RS	Bure R&G DB Timotion, PU, Plus, Locked direction
56-309-PU	Bure Rise&Go DB Timotion, 125mm, PU
56-309-PU-AUS	Bure Rise&Go DB Timotion, 125mm, PU, AUS
56-309-UK	Bure Rise&Go DB Timotion, 125 mm, UK
56-309-PU-UK	Bure Rise&Go DB Timotion, 125 mm, PU, UK
56-309/100	Bure Rise&Go DB Timotion, 100mm
56-309/75	Bure Rise&Go DB Timotion, 75mm
56-309/75-AUS	Bure Rise&Go DB Timotion, 75mm, AUS
56-309/75-B/0	Bure Rise&Go DB Timotion, 75mm, handbrakes
56-309/75-B/0-PLUS	Bure Rise&Go DB Timotion, 75mm, handbrakes, Plus
56-309/75-B/0-PS-PU-R	Bure Rise&Go DB Timotion, 75mm, handbrakes, PU, Plus
56-309/75-B/0-PU	Bure Rise&Go DB Timotion, 75mm, handbrakes, PU
56-309/75-Plus	Bure Rise&Go DB Timotion, 75mm, Plus
56-309/75-PLUS-PU	Bure Rise&Go DB Timotion, 75mm, PU, Plus
56-309/75-PLUS-PU-UK	Bure Rise&Go DB Timotion, 75mm, PU, Plus, UK
56-309/75-PU	Bure Rise&Go DB Timotion, 75mm, PU
56-309/75-PU-AUS	Bure Rise&Go DB Timotion, 75mm, PU, AUS
56-309/75-PU-UK	Bure Rise&Go DB Timotion, 75mm, PU, UK
56-309/75-PU-US	Bure Rise&Go DB Timotion, 75mm, PU, US
56-309/75-UK	Bure Rise&Go DB Timotion, 75mm, UK
56-309/75-US-TRG	Bure Rise&Go DB Timotion, 75mm, direction lock, US

Main configuration is following: 125 mm wheels, no handbrakes, cushions made of synthetic leather, regular hand control, no special features. Following variations are possible:

- 125/100/75 mm – indicates diameter of wheels
- Handbrakes – indicates that the model has a handbrake on both handles
- 1 handbrake – indicates that the model has only 1 handbrake on one of the handles
- UK, AUS, US – models designed specifically for these countries
- Timotion- means that the model includes electronics from Timotion for adjustment of the height
- Plus – means that the model includes a Plus hand control
- PU – means that the model includes cushions made of Polyurethane
- RS – means that the model has rotation stop.

# SystemRoMedic™

## Simple solutions for great results

SystemRoMedic® is the name of Direct Healthcare Group's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic® is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

## Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic® is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic® offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

## Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading at/from our website;  
[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)

## For questions about the product and its use

Please contact your local Direct Healthcare Group and SystemRoMedic® representative. A complete list of all our partners with their contact details can be found on our website; [www.directhealthcaregroup.com](http://www.directhealthcaregroup.com).

Moving Health Forward

 Type BF Applied Parts

Medical Device Class I. The product complies with the requirements of the Medical Devices Regulation 2017/745.

Direct Healthcare Group Sverige AB  
Torshamnsgatan 35,  
SE-164 40 Kista, Sweden

Tel: +46 (0)8-557 62 200  
[info@directhealthcaregroup.com](mailto:info@directhealthcaregroup.com)  
[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)