Instructions for use – English

Bruksanvisning – Svenska

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**Max Patient Weight: 150 kg / 330 lbs**

## Table of contents

|                        |     |
|------------------------|-----|
| English .....          | 3   |
| Svenska .....          | 15  |
| Norsk .....            | 27  |
| Danske .....           | 39  |
| Suomi .....            | 51  |
| Deutsch .....          | 63  |
| Netherlands .....      | 75  |
| Français .....         | 87  |
| EMC .....              | 99  |
| Dimensions .....       | 102 |
| Parts list .....       | 104 |
| Product overview ..... | 106 |

Bure Space is a walker that has been specifically designed for the orthopaedics departments and situations where patients require extra leg room to get up closer to wheelchairs or similar. Bure Space is available with electrically operated or hydraulic height adjustment.

## 1. Intended use

Walkers are indoor, movable and portable assistive aids, which together with approved accessories are intended to be used for transfer from sitting to standing positions, and walking training. Walkers are to be used by persons with reduced muscular strength, who shall be able to stand and perform simple leg movements. Walkers are to be used for leaning on when walking and standing during early mobilisation training. Walkers and accessories are intended for indoor use only in home healthcare environment and professional healthcare facility environment.

The Bure Space is intended to be used in orthopaedics department and situations where patients require extra leg room to get up closer to wheelchairs or similar.

## 2. Contraindications, Precautions, Warnings

### Contraindications

- The device may not be used by patients above the maximum weight indicated on the label of the device.

### Precautions

- Check that the device is used on a dry and flat surface.
- Check the device is correctly mounted/assembled before its first use.
- Check the device after every folding/disassembly after any transport.
- Check height adjustment and leg spreading motions and inspect the actuators or gas pistons full ranges.
- Activate the brakes whenever the device is not in use.
- Check the function of the brakes before use.
- Inspect the device regularly to detect any signs of damage. If there are any visible signs of damage, contact Direct Healthcare Group Customer Service.
- Check that the hand control does not show signs of wear before every use.
- Check that the hand control markings are in accordance with the lifting functions.
- Check the battery status of an electrically operated model before using the device.
- Warranty applies only if repairs or alterations are made by personnel who are authorized by Direct Healthcare Group.
- Ensure that there are no obstacles or people in the way of the device, when moving.
- Handle batteries with care. Do not drop.
- Use only batteries and cables that are intended for the device as instructed by the manufacturer.
- Verify that potential accessories that might be used are properly attached to the device before use.
- Low speed is recommended when moving with the device.
- Take care not to drive the device over thresholds with high speed or force. Approach the threshold with lower speed and communicate the coming threshold to the user.
- In case of an electrical model of the device, check that the hand control does not show signs of wear before every use.
- Keep the instructions for use together with the product so that it can be used as a reference if necessary.

- Beware of moving parts during usage.
- Perform a risk assessment on the user before using the Walker to ensure they can stand and perform simple movements.

### Warnings

- The device may not be used by infants, children or adults who are 146 cm or shorter.
- The caregiver must be able to read and understand the Manual/IFU of the product.
- Caregiver: Never use the device if you have not had a proper training on how to use the device including the accessories.
- Maintenance/Service staff: Never maintain the device if you have not had a proper training on how to use and maintain the device including the accessories.
- It is important to use only approved accessories to prevent unintended detachment of components and subsequently a fall that may lead to patient injury.
- Use careful and gentle maneuvers when moving the device.
- Perform maintenance/service of the device, according to the instructions in the IFU (Instruction for use), at least once every 12 months.
- No parts of the device shall be serviced, while the device is in use
- Accessories must be properly fitted and tested in relation to the user's needs and functional ability.
- Special care must be taken when using strong electrical power sources such as diathermy and the like so that diathermy cables are not placed on or near the device. In case of doubt consult with a Direct Healthcare Group representative.
- When lifting the device, only use the denoted handholds. At least two people are required for lifting the device.
- Do not leave a patient unattended when the device is being used.
- The device can not be used by patient alone.
- Never move the walker by pulling on the actuator device or gas piston.
- The device must not be lowered into water.
- The device must not be cleaned using steam.
- The device must not be left or stored in a damp or humid environment.
- The device must not be charged in a wet room.
- The device must not be used outdoors, only indoors.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective grounding.
- No modification of this equipment is allowed
- Do not modify this equipment without authorisation of the manufacturer.
- If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of this equipment.
- Do not use the device while it is charging.
- Do not use the device if any of the cables shows signs of wear and tear.
- Do not put the device in direct sunlight as it may get hot.
- Do not run over the cables with the walker or other devices.

**Always read the user manual**

Always make sure that you have the correct version of the manual. The most recent version of all manuals are available for downloading at/from our website: [www.directhealthcaregroup.com](http://www.directhealthcaregroup.com).

**Before first use:**

- Inspect the packaging for any damage.
- Check for correct product being delivered.
- Check to ensure that the device is free from damage.
- Check if all the parts/components are included in the package.
- Check if all approved accessories are included in the package.
- Check if Product Quality Approval document is included in the packaging. Save this document for future contact with manufacturer.
- Make certain the device is properly assembled.
- Check lifting function and bottom-frame adjustment function.
- Check if it is possible to move the device in all directions.




### **3. Using Bure Walker**

**Standing up and walk training**

1. Ask the user to sit on the edge of the bed or chair. Bring out the walking table and place it in front of the user. Lock the brakes on the wheels.
2. Ask the user to place their arms on the armrests. If necessary, ask the user to move forward on the seat base so that the user reaches the handles.
3. Ask the user to lean forward. To start raising the walking table, press the UP button on the controller. Raise the walking table to a desired height.
4. Unlock the brakes. Adjust the height of the walking table to a desired position
5. Assist the user while walk training.

## 4. Using Bure Walker with Gait Training Kit 56-388

- The support harness has four blue straps with a plastic buckle allowing for adjustments in length, having steel carabiner clips at the ends to be fasten to the walker. The carabiner clips on to the metal mounts placed under the table of the walker.
- Apply the walker support harness by clipping on the two front carabiners to the front mounters, one on each side of the table.
- Place the user standing up in the walker. Pull back the harness in between the legs and clip on the carabiners to the back mounters, one on each side of the table. Adjust the harness to the user by pulling in the back strap loops.
- **Observe!** Intended use of the trainer support harness is to practise walking. **It is not** a sitting or lifting aid.

### Standing Plate

Attach standing plate by hooks on the main frame of the walker.

## Raising/Lowering the Frame

Adjust the height with the aid of the remote control UP and DOWN arrows.

Handset Plus has charge indicator / battery (low, medium, high) and Service indicator (see further under "Hand Control Plus").

## Adjusting the Handles

To adjust handles to the desired position; loosen screw, adjust handle, tighten screw.

## Left and Right Adjustment of Armrests

The armrests can be adjusted sideways for optimal support; loosen screw, adjust armrest, tighten screw.

NOTE: The rear attachments on the frame/armrest are interchangeable (left with right and vice versa) to provide different widths. The locating pins on the cushion's underside/rear can also be moved. These adjustment capabilities make customisation possible.

## Using the Parking Brakes on Castor Wheels

When parking the device, activate the parking brake on the lockable castor wheels. Lock the wheels by pressing the green pedal brake down with your foot. To release the brake, press the green protruding part with your foot. The caregiver should assist the patient in activating the parking brake in case the patient cannot activate them him or herself.

## Using Hand Brakes

**Running brakes:** Pull the brake levers towards you to reduce speed.

You will brake the rear wheel that is located on the same side as the brake lever.

If you apply only one of the brakes, the walker will turn to that side. When the lever is released, the braking will stop. Do not use the brakes continuously while the walker is in motion.

**Parking brakes:** Push the brake levers forwards from you to activate the parking brakes.

Be aware that the lever has some resistance before it stops and locks the rear wheel.

Pull the brake levers back towards you to release the parking brakes.

## Hand Control Plus

\*\*1. Service indicator – flashes red after 11 months' operation.

Constant light (red) after 12 months' operation.

The indicator is a reminder to review the maintenance instruction (see Maintenance Instructions). This should be done at least once every 12 months.

Reset: hold down buttons <2a> + <2b> for at least 30 seconds. \*\*The indicator will go out around 30 seconds after the hand control is activated (battery conservation function).

## Battery Charging

Battery charging must take place when the device is not in use.

Charging starts automatically when power plug is connected to 230V wall socket. Green light on the control box will turn on, once the power plug is connected.

- The LED flashes green during charging.
- ON/OFF: To use the product, press ON.
- The green lamp will begin to flash.
- Use the OFF button (hold down for 5 s) to switch off the battery if the product will not be used for an extended period.
- The LED will show a constant green light when the battery is fully charged.
- When the battery is in normal state/in operation, the LED flashes green at REGULAR intervals, namely 0.5 s ON, 4 s OFF.
- When the battery requires charging (remaining charge <25%) the LED will flash orange/red at REGULAR intervals, namely 0.5 s ON, 4 s OFF.

NOTE: The battery has an in-built device that automatically switches it off when remaining charge is <20%. The LED is then extinguished. The battery must now be fully charged before the electronics can be used again.

## Control box

The control box cover can be opened by using a screwdriver or similar; refer to the instructions which are sitting on the control box.

NOTE: Make sure that plugs are connected to the correct sockets otherwise the actuators can be damaged.

1. Connector for actuator up / down.
2. Hand control connector.
3. Battery outlet.
4. Connector for mains cable. NOTE: The mains cable must ALWAYS be connected to the control box for IP protection to be valid.

## Trouble shooting

If the lifting function cannot be activated, check the following:

- That all cables are properly and securely connected, Pull out the contact and plug it in again firmly.
- That battery charging is not in progress.
- That the battery is charged.

If the device is not working properly, contact your distributor.

If the device makes unusual noises:

- Try to determine the source of the sound. Take the device out of operation and contact your distributor.

## Accessories

| ART NR.         | DESCRIPTION                             |
|-----------------|-----------------------------------------|
| 58-328-1        | Heel strap Standing plate               |
| 56-336-4        | Oxygen holder                           |
| 56-338          | Drip rod complete with attachment       |
| 56-337          | Side support, adjustable                |
| 56-337-PU       | Side support, adjustable PUR            |
| 56-339-1-0      | Handbrake kit/2 handles                 |
| 56-339-1-1      | Handbrake kit/One-handed brake          |
| 56-352          | Lower Leg Support                       |
| 56-388          | training belt                           |
| 57-375-X1       | Hand control PLUS (TH10) incl.holder    |
| 56-333-SPACE    | Standing plate Short Bure Space         |
| 56-333/SW-SPACE | Standing plate Short, SW, Bure Space    |
| 56-383-SPACE    | Standing plate Long Bure Space          |
| 56-383/SW-SPACE | Standing plate Long, SW, Bure Space     |
| 56-305          | Touch-up paint Bure, Silver 8           |
| 56-306          | Touch-up paint Bure, Dark grey RAL 7021 |

## 5. Maintenance

The device must undergo thorough inspection at least once per year. Inspection must be performed by authorised personnel and in accordance with Direct Healthcare Group's instructions.

Repairs and maintenance may only be done by authorised personnel using original spare parts.


| Check that                                                                                                                                                                                                                                                                                                                                                                                     | Action if a problem is noticed                                                                                                                                                                                                                                                                                                                                                                                                 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>General</b><br>The Walker feels firm/free from play.<br>The Walker does not rattle during manoeuvres.<br>The Walker is level and all of its castors are in contact with the floor.<br>The Walker is not dirty.                                                                                                                                                                              | Clean with lukewarm soap solution or alcohol-based cleaning agent (no petroleum products).                                                                                                                                                                                                                                                                                                                                     |
| <b>Armrests</b><br>Armrests are intact and clean.<br>Armrest widening is functional.                                                                                                                                                                                                                                                                                                           | Clean; replace if damaged.<br>Replace lock knob.<br>Fit lock screws/plates as required.                                                                                                                                                                                                                                                                                                                                        |
| <b>Handle</b><br>Grips are not damaged/dirty.<br>Handle adjustment is functional.                                                                                                                                                                                                                                                                                                              | Clean with lukewarm soap solution (no petroleum products).<br>Replace handle grips.<br>Fit lock screws/plates as required.                                                                                                                                                                                                                                                                                                     |
| <b>Electronics</b><br>Check that the charging cable is connected to the control box.<br>Check that the hand control is connected to the control box.<br>Check that the battery, actuators and control box are not loose.                                                                                                                                                                       | If any part needs to be replaced refer to the wiring diagram in the control box wiring section.<br>Fit new charging cable; this must always be connected to the control box.<br>Connect or fit new hand control as required.<br>Tighten; replace with new fasteners as required.                                                                                                                                               |
| <b>Height adjustment</b><br>Raising and lowering are functional.<br>The Walker feels stable at maximum settings.<br>There is no play between the vertical frame and bottom frame.<br>Height adjustment clamps lock.<br>The electric actuator fitting must be free from play.<br>Height adjustment via hand control is functional. The electric actuator must run smoothly at a constant speed. | Fit new guide sleeves (in frame) or end plugs (in chrome tubes) as required.<br>Tighten all lock bolts.<br>Replace clamps.<br>Tighten the attachment fitting concerned; replace bolts/lock nuts as necessary.<br>Check that the battery/hand control/ actuator are connected in accordance as in the control box wiring section.<br>Charge the battery.<br>The part list shows the design and which fasteners must be checked. |
| <b>Frame parts</b><br>There is no mechanical damage.<br>There are no scratches.<br>End plugs/lower frame fitted                                                                                                                                                                                                                                                                                | If mechanical damage is present, contact DHG Customer Services.<br>Touch up as necessary.<br>Fit new end plugs.                                                                                                                                                                                                                                                                                                                |
| <b>Castors/brakes</b><br>Castors roll easily/tread not damaged.<br>The castors are firmly fastened to the lower frame.<br>Castor brakes functional on all castors.<br>Handbrake functional.                                                                                                                                                                                                    | Clean or replace castors.<br>NOTE! Castors are always fixed to the lower frame with thread locker or lock nuts. The castors are sealed and we do not recommend their dismantling; instead, replace the whole castor assembly.<br>Tighten the castor bolt and fit locking nut or use thread locker (depends on model).<br>Replace castor assembly.<br>Adjust the brake or fit new castors.                                      |

Power cord may only be replaced by DHG service personnel or technical personnel trained by DHG.

DHG can provide upon request necessary technical information in order to assist maintenance and repair of the device.

## Storage and transportation

The walker should be transported and stored in 5 °C to 50 °C and in normal humidity, 20% to 80% non-condensing. The air pressure should be between 700 and 1060 hPa. Leftmost symbol indicates storage and transportation. Let the walker reach room temperature before the batteries are charged or the walker is used. The walker should not be stored so that it is exposed to dust, or so that the battery or gas piston is exposed to direct sunlight.

## Operation

The operating environment should be 5 °C to 40°C, relative humidity 20% to 80% non-condensing, and atmospheric pressure 700 to 1060 hPa. See also Technical specification below. Leftmost symbol indicates operating condition.

## Cleaning and disinfection

After every use:

Clean the walker armrests and handles. Wipe using a mild soap solution/surface disinfectant. However, the product is not corrosion proof.

Weekly and possibly more often depending on frequency of use:

If necessary, clean the whole walker with a cloth with warm water or a soap solution and check that the castors are free from dirt and hair. The electronic components (battery, control box, hand control, actuators, cables) should be cleaned with a damp cloth only. To avoid degreasing of the piston rods, the actuators should be retracted to minimum stroke and without load before cleaning. Ensure that the device is dried thoroughly after cleaning.

If disinfection is needed:

70% ethanol, 45% isopropanol or similar should be used.

Do not steam clean due to risk of corrosion.

Do not use cleaning agents containing phenol or chlorine, as this could damage the materials.

With environments with high risk of infections use the walker with PU-cushions (Accessory number 57-361-L and 57-361-R)

## General Maintenance for gait training Sling:

Read the product label.

Do not use rinsing agent.

To prolong product life, avoid tumble-drying.

## Recycling/Disposal

When the product is no longer usable it should be recycled in accordance with legislation and regulations in the country concerned. All electrical parts including batteries must be removed and recycled as electrical components. Ask your local recycling station for further information about how the different materials in the product (metals, plastics, electronics) should be recycled. Used batteries can also be returned to Direct Healthcare Group or a Direct Healthcare Group distributor for recycling.

## 6. Technical Specifications

|                                               |                                                                                                                                                                                                                                                                                            |
|-----------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Control box</b>                            | TiMotion, Voltage 220V-240VAC, 50Hz, 1250mA/24VDC, 70VA, IP66, Class II (Doubled insulated)                                                                                                                                                                                                |
| <b>Batteries</b>                              | TiMotion, 24 V - 4A, 96Wh, Rechargeable Li-Ion Battery, IP66, Cable length 600mm, Duty cycle: 10%, Max. 2min. ON/18 min./OFF                                                                                                                                                               |
| <b>Actuator</b>                               | 12VDC, Max 8,0A, Push: 2000N, Self-locking capacity (min): 1000N, Stroke 400mm +/- 2 mm, IP66.<br>Duty cycle: 10%, Max. 2min.ON/18min. OFF                                                                                                                                                 |
| <b>Hand Control</b>                           | Cable length straight 300mm + coiled 400 mm, IP66                                                                                                                                                                                                                                          |
| <b>IP class</b>                               | IPX5                                                                                                                                                                                                                                                                                       |
| <b>Material</b>                               | Steel                                                                                                                                                                                                                                                                                      |
| <b>Castors</b>                                | 125/100/75mm                                                                                                                                                                                                                                                                               |
| <b>Weight</b>                                 | 20 kg                                                                                                                                                                                                                                                                                      |
| <b>Expected Service life</b>                  | 8 years                                                                                                                                                                                                                                                                                    |
| <b>Max patient weight</b>                     | 150 kg / 330 lbs                                                                                                                                                                                                                                                                           |
| <b>Operating environment</b>                  | 5 °C to 40°C, relative humidity 20% to 80% non-condensing, and atmospheric pressure 700 to 1060 hPa.                                                                                                                                                                                       |
| <b>Storage and transportation environment</b> | 5 °C to 50 °C, relative humidity 20% to 80% non-condensing, and atmospheric. Environment pressure 700 to 1060 hPa.                                                                                                                                                                         |
| <b>Height</b>                                 | 125 mm wheels model minimum 900 mm and maximum 1300 mm<br>100 mm wheels model minimum 880 mm and maximum 1280 mm<br>75mm wheels model minimum 850 mm and maximum 1250 mm                                                                                                                   |
| <b>Length</b>                                 | 125mm wheels 850mm<br>100mm wheels 835mm<br>75mm wheels 830mm                                                                                                                                                                                                                              |
| <b>Width</b>                                  | 775 mm                                                                                                                                                                                                                                                                                     |
| <b>Turning Width</b>                          | 125 mm wheels: 1158 mm, 100 mm wheels: 1130 mm, 75 mm wheels: 1130 mm                                                                                                                                                                                                                      |
| <b>Materials in contact with human skin</b>   | Materials that user and operator/caregiver are exposed to during usage of the product for prolonged time:<br>Poroflex (synthetic leather cushions) P <sup>x</sup><br>Polyurethane (on PU cushions) P<br>Dryflex (handles) P<br>PA6 + 20% GF (brake lever) P<br>Hand control C <sup>x</sup> |
| <b>Applied parts:</b>                         | Armrest P<br>Handles P<br>Knee pad P<br>Foot rest support P<br>Foot rest support clamps C<br>Handles adjusting knobs C<br>Hand controller and its holder. C<br>Handles to adjust the armrest width. C<br>Hand control C<br>Charger C<br>Charger cable C<br>Footplate                       |

<sup>x</sup> P - patient, <sup>x</sup> C - caregiver

7. Product label

**Bure Space El, Timotion**  
**SystemRoMedic**

**REF** 56-304

UDI (Unique Device Identification) can be found within the barcode human readable text.  
It contains.

(01)07350013304019(11)220301(21)610231377

**LOT** 2022-03-R

(01) EAN code

(11) Date of Manufacture. The date is in the format

YYMMDD

Y = Year, M = Month and D = Date.

150 kg / 330 lbs

230 VAC, 50/60 Hz

IPX5

775mm

900-1300 mm

**Direct Healthcare Group Sverige AB**

Torshamnsgatan 35, 16440 Kista, SWEDEN

[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)

Made in Sweden

Rev .01


(21) Serial Number displayed as five digit

code.

Product label

## 8. Table of Symbols

|                                                                                    |                                                                                       |                                                                                    |                                       |
|------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|---------------------------------------|
|   | This product complies with the requirements of the Medical Device Regulation 2017/745 |   | Visual Inspection                     |
|   | Medical Device                                                                        |   | Read the manual                       |
|   | Caution                                                                               |   | Product Code                          |
|   | Legal manufacturer<br>Direct Healthcare Group<br>Sverige AB                           |   | Batch Code                            |
|   | Type BF applied parts, according to the degree of protection against electric shock   |   | The device is intended for indoor use |
|   | WEEE Symbol<br>May not be discarded in domestic waste                                 |   | Class II Equipment                    |
|   | Duty Cycle:<br>1 min in active (ON) mode.<br>19 min in rest (OFF) mode.               |   | Maximum patient weight                |
|   | Walker Height                                                                         |   | Walker Width                          |
|   | Do not step on the device                                                             |   | Refer to instruction manual (IFU)     |
|  |                                                                                       | Weight (mass) of the device, the Maximum Patient Weight, and total sum. All in kg. |                                       |

## 9. How to report a serious incident

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the MHRA or another competent authority of the country in which the user and/or patient is established.

UK

info@directhealthcaregroup.com

T: +44 (0) 800 043 0881

F: +44 (0) 845 459 9832

Other:

info.export@directhealthcaregroup.com

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T: +46 8-557 62 200

# Electromagnetic Compatibility (EMC)

## Electromagnetic emissions and test levels

The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

| Emission test                              | Standards     | Compliance |
|--------------------------------------------|---------------|------------|
| RF emissions                               | CISPR 11      | Group 1    |
| RF emissions                               | CISPR 11      | Class B    |
| Harmonic current emissions                 | IEC 61000-3-2 | Class A    |
| Voltage fluctuations and Flicker emissions | IEC 61000-3-3 | Complies   |

## Electromagnetic Immunity

Acceptance criteria for the EMC pass criteria is unintentional movement above  $\pm 10\text{mm}$  is not allowed. The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Acceptance criteria for the EMC pass criteria is unintentional movement above  $\pm 10\text{mm}$  is not allowed. The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.



| Radiated Fields in close proximity, Immunity test and compliance.<br>Standard: IEC 61000-4-39 |                                |                      |                |
|-----------------------------------------------------------------------------------------------|--------------------------------|----------------------|----------------|
| Dwell time: 3 sec                                                                             |                                |                      |                |
| Test levels (A/m)                                                                             | Modulation<br>Pulse modulation | Mod. Frequency (kHz) | Test Frequency |
| 8                                                                                             | CW                             | CW                   | 30 kHz         |
| 65                                                                                            | 50% Duty cycle                 | 2.1                  | 134.2 kHz      |
| 7.5                                                                                           | 50% Duty cycle                 | 50                   | 13.56 MHz      |

| Immunity test                               | Basic EMC standard or test method | IEC 60601-1-2 Edition 4 Test levels and compliance                                                   |                                                                            |
|---------------------------------------------|-----------------------------------|------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------|
|                                             |                                   | Professional healthcare facility environment†                                                        | Home healthcare environment                                                |
| Surge                                       | IEC 61000-4-5                     | ± 0.5 kV, ± 1 kV                                                                                     |                                                                            |
| Voltage dips                                | IEC 61000-4-11                    | 0% U <sub>T</sub> ; 0.5 cycle<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°                    |                                                                            |
|                                             |                                   | 0% U <sub>T</sub> ; 1 cycle<br>and<br>70% U <sub>T</sub> ; 25/30 cycles.<br>Single phase: at 0°      |                                                                            |
| Voltage Interruptions                       | IEC 61000-4-11                    | 0% U <sub>T</sub> ; 250 / 300 cycle                                                                  |                                                                            |
| Conducted disturbances induced by RF fields | IEC 61000-4-6                     | 6V in ISM bands between 0.15 MHz and 80 MHz<br>80% AM at 1 kHz                                       | 6V in ISM and amateur bands between 0.15 MHz and 80 MHz<br>80% AM at 1 kHz |
|                                             |                                   | Dwell time: 3 sec<br>Frequency step size: 1%                                                         |                                                                            |
| Radiated RF Electro-magnetic fields         | IEC 61000-4-3                     | 3 V/m<br>80 MHz – 2.7 GHz<br>80% AM at 1 kHz                                                         | 10 V/m<br>80 MHz – 2.7 GHz<br>80% AM at 1 kHz                              |
|                                             |                                   | Dwell time: 3 sec<br>Frequency step size: 1%                                                         |                                                                            |
| Electrical fast transients / burst          | IEC 61000-4-4                     | ± 2 kV<br>100 kHz repetition frequency                                                               |                                                                            |
| RATED Power frequency magnetic field        | IEC 61000-4-8                     | 30 A/m<br>50 Hz or 60Hz                                                                              |                                                                            |
| Electrostatic Discharge (ESD)               | IEC 61000-4-2                     | Contact: ± 8 kV<br>Air: ± 2 kV, ±4 kV, ± 8 kV, ± 15 kV<br>Number of discharges: 10 for each polarity |                                                                            |

Radiated Radio – Frequency, Proximity fields from wireless communication equipment & compliance. (Dwell time: 3 sec.)  
Standard: IEC 61000-4-3

|                              | Test Frequency (MHz) | Band (MHz)      | Service                                                            | Modulation                            | Immunity test levels (V/m) |
|------------------------------|----------------------|-----------------|--------------------------------------------------------------------|---------------------------------------|----------------------------|
| Spot Frequencies             | 385                  | 380 to 390      | TETRA 400                                                          | Pulse modulation<br>18 Hz             | 27                         |
|                              | 450                  | 430 to 470      | FMRS 460, FRS 460                                                  | FM<br>± 5 kHz Deviation<br>1 kHz sine | 28                         |
|                              | 710                  | 704 to 787      | LTE Band 13, 17                                                    | Pulse modulation<br>217 Hz            | 9                          |
|                              | 745                  |                 |                                                                    |                                       |                            |
|                              | 780                  |                 |                                                                    |                                       |                            |
|                              | 810                  | 800 to 960      | GSM 800/900,<br>TETRA 800,<br>iDEN 820,<br>CDMA 850, LTE<br>Band 5 | Pulse modulation<br>18 Hz             | 28                         |
|                              | 870                  |                 |                                                                    |                                       |                            |
|                              | 930                  |                 |                                                                    |                                       |                            |
|                              | 1720                 | 1700 to<br>1990 | GSM 800/900,<br>TETRA 800,<br>iDEN 820,<br>CDMA 850, LTE<br>Band 5 | Pulse modulation<br>217 Hz            | 28                         |
|                              | 1845                 |                 |                                                                    |                                       |                            |
|                              | 1970                 |                 |                                                                    |                                       |                            |
|                              | 2450                 | 2400 to<br>2570 | Bluetooth,<br>WLAN, 802.11<br>b/g/n, RFID<br>2450, LTE<br>Band 7   | Pulse modulation<br>217 Hz            | 28                         |
|                              | 5240                 | 5100 to<br>5800 | WLAN 802.11<br>a/n                                                 | Pulse modulation<br>217 Hz            | 9                          |
|                              | 5500                 |                 |                                                                    |                                       |                            |
|                              | 5785                 |                 |                                                                    |                                       |                            |
| Other Identified frequencies | 433                  | -               | -                                                                  | Pulse modulation 2 Hz                 | 3                          |
|                              | -                    | 860-960         | -                                                                  | Pulse modulation 2 Hz                 | 54                         |
|                              | 2450                 | -               | -                                                                  | Pulse modulation 2 Hz                 | 54                         |


Dimensions

| Wheel diameter | A: Turning width | B: Armrest height | B: PU Armrest height | C: Wheel area diameter | D: Length |
|----------------|------------------|-------------------|----------------------|------------------------|-----------|
| 125mm          | 1158             | 900-1300          | 870-1270             | 218                    | 850       |
| 100mm          | 1130             | 880-1280          | 850-1250             | 190                    | 835       |
| 75mm           | 1130             | 850-1250          | 820-1220             | 188                    | 830       |

Parts

| No. | ART No.                    | DESCRIPTION                                                                     | No. | ART No.              | DESCRIPTION                                           |
|-----|----------------------------|---------------------------------------------------------------------------------|-----|----------------------|-------------------------------------------------------|
| 1   | 58-373/Lion                | Control box 220V Digital Li-On<br>TiMotion Bure                                 | 24  | 57-812               | Tapping screw RXS 2,9x4,5                             |
| 2   | 58-377                     | Actuator TiMotion 2000N/5200rpm,<br>Bure                                        | 25  | 57-811               | Serrated lock washer M16                              |
| 3   | 57-356                     | Armrest frame                                                                   | 26  | 57-925               | Endcover pipe 39mm black                              |
| 4   | 58-389                     | Battery box 4.0 Ah TBB5/Li-On                                                   | 27  | 56-308               | Cable holder Bure charger                             |
| 5   | 57-356-1                   | Bure Battery plate                                                              | 28  | 57-808               | Pin bolt M6x25                                        |
| 6   | 58-374-X1                  | Hand control TiMotion TH7                                                       | 29  | 57-825               | Washer 10,5x18x2                                      |
| 7   | 57-373                     | Bottom frame Bure Space                                                         | 30  | 57-828               | Pop rivet 4,8x10                                      |
| 8   | 57-355-1                   | Middle frame                                                                    | 31  | 57-792               | Washer Cushion                                        |
| 9   | 58-335                     | Telescopic Inner pipe                                                           | 32  | 57-848               | Nordlock 10,7x21x2,5                                  |
| 10  | 57-747                     | Plug dome plastic 25dia                                                         | 33  | 57-817               | Lock nut CL6 M10                                      |
| 11  | 57-049<br>57-319<br>57-051 | Castor wheel 125mm brake<br>Castor wheel 100mm brake<br>Castor wheel 75mm brake | 34  | 57-918               | Washer 6,5x15x10                                      |
| 12  | 57-307                     | Protective cover plastic                                                        | 35  | 57-904               | Washer 6.4x12.5x2 nylon                               |
| 13  | 57-310                     | Middle frame bushing                                                            | 36  | 57-877               | Washer 10,5x25x3                                      |
| 14  | 57-041-2                   | Holder remote control, Bure                                                     | 37  | 57-357-L<br>57-361-L | Armrest Cushion Low Left<br>Armrest Cushion PU Left   |
| 15  | 57-818                     | Screw K6S M6x10                                                                 | 38  | 57-357-R<br>57-361-R | Armrest Cushion Low Right<br>Armrest Cushion PU Right |
| 16  | 57-816                     | Screw MC6S M10x40                                                               | 39  | 57-919               | Star knob handle                                      |
| 17  | 57-803                     | Lock washer M10 6.1x11.8x1.6                                                    | 40  | 57-820               | Washer BRB 5.3x10x1                                   |
| 18  | 57-816                     | Screw MC6S M8x35                                                                | 41  | 57-821               | Screw MC6S M5x20                                      |
| 19  | 57-805                     | Washer BRB 6,4x12x1,6                                                           | 42  | 57-850               | Cap nut M5                                            |
| 20  | 56-353                     | Bure Handle assy (pair)                                                         | 43  | 57-790               | Plug Black                                            |
| 21  | 57-917                     | Wing nut M8                                                                     | 44  | 58-331               | Cover actuator rod                                    |
| 22  | 57-806                     | Washer 8.4x16x1.5                                                               | 45  | 58-326               | Spiral Cord EUR - TiMotion,<br>Bure                   |
| 23  | 57-810                     | Screw MC6S M8x45                                                                |     |                      |                                                       |

# Product Overview

**REF**

| Article number     | Product name                                      |
|--------------------|---------------------------------------------------|
| 56-304             | Bure Space Timotion, 125 mm                       |
| 56-304-PU          | Bure Space Timotion, 125mm, PU                    |
| 56-304-PLUS-PU     | Bure Space Timotion, 125mm, PU, Plus              |
| 56-304-PLUS        | Bure Space Timotion, 125mm, Plus                  |
| 56-304-B/1-PL-PU   | Bure Space Timotion, 125mm, 1 handbrake, PU, Plus |
| 56-304-B/0-PU      | Bure Space Timotion, 125mm, handbrakes, PU        |
| 56-304/100         | Bure Space Timotion, 100mm                        |
| 56-304/100-PU-PLUS | Bure Space Timotion, 100mm, PU, Plus              |
| 56-304/75          | Bure Space Timotion, 75mm                         |
| 56-304/75-PU       | Bure Space Timotion, 75mm, PU                     |
| 56-304/75-PLUS     | Bure Space Timotion, 75mm, Plus                   |
| 56-304/75-PLUS-PU  | Bure Space Timotion, 75mm, PU, Plus               |
| 56-304/75-B/1      | Bure Space Timotion, 75mm, 1 handbrake            |
| 56-304/75-B/0      | Bure Space Timotion, 75mm, handbrakes             |

Main configuration is following: 125 mm wheels, no handbrakes, cushions made of synthetic leather, regular hand control, no special features. Following variations are possible:

- 125/100/75 mm – indicates diameter of wheels
- Handbrakes – indicates that the model has a handbrake on both handles
- 1 handbrake – indicates that the model has only 1 handbrake on one of the handles
- UK, AUS, US – models designed specifically for these countries
- Timotion- means that the model includes electronics from Timotion for adjustment of the height
- Plus – means that the model includes a Plus hand control
- PU – means that the model includes cushions made of Polyurethane
- RS – means that the model has rotation stop.

## Simple solutions for great results

SystemRoMedic® is the name of Direct Healthcare Group's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic® is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

### Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic® is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic® offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

## Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading at/from our website;  
[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)

## For questions about the product and its use

Please contact your local Direct Healthcare Group and SystemRoMedic® representative. A complete list of all our partners with their contact details can be found on our website; [www.directhealthcaregroup.com](http://www.directhealthcaregroup.com).

Medical Device Class I. The product complies with the requirements of the Medical Devices Regulation 2017/745.

Direct Healthcare Group Sverige AB  
Torshamnsgatan 35,  
SE-164 40 Kista, Sweden

Tel: +46 (0)8-557 62 200  
[info@directhealthcaregroup.com](mailto:info@directhealthcaregroup.com)  
[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)