**Max User Weight: 150 kg / 330 lbs**

## Table of contents

English .....	4
Svenska .....	19

## Product Overview

### REF

56-304  
56-304-B/0-FKN\_R  
56-304-B/0-FKH\_R-PU  
56-304-B/0-FKH\_R-PL-PU-R  
56-304-B/0-PU  
56-304-B/1-PL-PU  
56-304-PLUS  
56-304-PLUS-PU  
56-304-PLUS-RED  
56-304-PU  
56-304-RED  
56-304/100  
56-304/75  
56-304/75-B/0  
56-304/75-B/0-FKH\_R-PS-PU  
56-304/75-B/0-FKN\_R-PU  
56-304/75-B/1  
56-304/75-PLUS  
56-304/75-PLUS-PU  
56-304/75-PU  
56-307  
56-307/100  
56-307/75  
56-307-B/1-PU  
56-307-B/1-PU-RED  
56-307-PU  
56-307/LOW

Bure Space is a walker that has been specifically designed for the orthopaedics departments and situations where patients require extra leg room to get up closer to wheelchairs or similar. Bure Space is available with electrically operated or hydraulic height adjustment.

## 1. Intended use

Walkers are indoor, movable and portable assistive aids, which together with approved accessories are intended to be used for transfer from sitting to standing positions, and walking training. Walkers are to be used by persons with reduced muscular strength, who shall be able to stand and perform simple leg movements. Walkers are to be used for leaning on when walking and standing during early mobilisation training. Walkers and accessories are intended for indoor use only in Home Healthcare environment and Professional Healthcare facility environment.

The Bure Space is intended to be used in orthopaedics department and situations where patients require extra leg room to get up closer to wheelchairs or similar.

## 2. Contraindications, Precautions, Warnings

### Contraindications

- The device may not be used by patients above the maximum weight indicated on the label of the device.

### Precautions

- Check that the device is used on a dry and flat surface.
- Check the device is correctly mounted/assembled before its first use.
- Check the device after every folding/disassembly after any transport.
- Check height adjustment and leg spreading motions and inspect the actuators or gas pistons full ranges.
- Activate the brakes whenever the device is not in use.
- Check the function of the brakes before use.
- Inspect the device regularly to detect any signs of damage. If there are any visible signs of damage, contact Direct Healthcare Group Customer Service.
- Check that the hand control does not show signs of wear before every use.
- Check that the hand control markings are in accordance with the lifting functions.
- Check the battery status of an electrically operated model before using the device.
- Warranty applies only if repairs or alterations are made by personnel who are authorized by Direct Healthcare Group.
- Ensure that there are no obstacles or people in the way of the device, when moving.
- Handle batteries with care. Do not drop.
- Use only batteries and cables that are intended for the device as instructed by the manufacturer.
- Verify that potential accessories that might be used are properly attached to the device before use.
- Low speed is recommended when moving with the device.
- Take care not to drive the device over thresholds with high speed or force. Approach the threshold with lower speed and communicate the coming threshold to the user.
- In case of an electrical model of the device, check that the hand control does not show signs of wear before every use
- Keep the instructions for use together with the product so that it can be used as a reference if necessary.

- Beware of moving parts during usage.
- Perform a risk assessment on the user before using the Walker to ensure they can stand and perform simple movements.

### Warnings

- The device may not be used by infants, children or adults who are 146 cm or shorter.
- The caregiver must be able to read and understand the Manual/IFU of the product.
- Caregiver: Never use the device if you have not had a proper training on how to use the device including the accessories.
- Maintenance/Service staff: Never maintain the device if you have not had a proper training on how to use and maintain the device including the accessories.
- It is important to use only approved accessories to prevent unintended detachment of components and subsequently a fall that may lead to patient injury.
- Use careful and gentle maneuvers when moving the device.
- Perform maintenance/service of the device, according to the instructions in the IFU (Instruction for use), at least once every 12 months.
- No parts of the device shall be serviced, while the device is in use
- Accessories must be properly fitted and tested in relation to the user's needs and functional ability.
- Special care must be taken when using strong electrical power sources such as diathermy and the like so that diathermy cables are not placed on or near the device. In case of doubt consult with a Direct Healthcare Group representative.
- When lifting the device, only use the denoted handholds. At least two people are required for lifting the device.
- Do not leave a patient unattended when the device is being used.
- The device can not be used by patient alone.
- Never move the walker by pulling on the actuator device or gas piston.
- The device must not be lowered into water.
- The device must not be cleaned using steam.
- The device must not be left or stored in a damp or humid environment.
- The device must not be charged in a wet room.
- The device must not be used outdoors, only indoors.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective grounding.
- No modification of this equipment is allowed
- Do not modify this equipment without authorisation of the manufacturer.
- If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of this equipment.
- Do not use the device while it is charging.
- Do not use the device if any of the cables shows signs of wear and tear.
- Do not put the device in direct sunlight as it may get hot.
- Do not run over the cables with the walker or other devices.

**Always read the user manual**

Always make sure that you have the correct version of the manual. The most recent version of all manuals are available for downloading at/from our website: [www.directhealthcaregroup.com](http://www.directhealthcaregroup.com).

#### **Before first use:**

- Inspect the packaging for any damage.
- Check for correct product being delivered.
- Check to ensure that the device is free from damage.
- Check if all the parts/components are included in the package.
- Check if all approved accessories are included in the package.
- Check if Product Quality Approval document is included in the packaging. Save this document for future contact with manufacturer.
- Make certain the device is properly assembled.
- Check lifting function and bottom-frame adjustment function.
- Check if it is possible to move the device in all directions.

### **3. Using Bure Walker**

#### **Standing up and walk training**

1. Ask the user to sit on the edge of the bed or chair. Bring out the walking table and place it in front of the user. Lock the brakes on the wheels.
2. Ask the user to place their arms on the armrests. If necessary, ask the user to move forward on the seat base so that the user reaches the handles.
3. Ask the user to lean forward. To start raising the walking table, press the UP button on the controller. Raise the walking table to a desired height.
4. Unlock the brakes. Adjust the height of the walking table to a desired position
5. Assist the user while walk training.

With gas piston:

1. Ask the user to sit on the edge of the bed or chair. Bring out the walking table and place it in front of the user. Lock the brakes on the wheels.
2. Raise the walking table by pressing the handle on the gas piston.
3. Ask the user to lean forward, grasp the armrests and stand up.
4. Adjust the height of the walking table to a desired position

## 4. Using Bure Walker with Gait Training Kit 56-388

- The gait support harness has four blue straps with a plastic buckle allowing for adjustments in length, having steel carabiner clips at the ends to be fasten to the gait walker. The carabiner clips on to the metal mounts placed under the table of the Gait walker.
- Apply the gait walker support harness by clipping on the two front carabiners to the front mounters, one on each side of the table.
- Place the user standing up in the Gait walker. Pull back the harness in between the legs and clip on the carabiners to the back mounters, one on each side of the table. Adjust the harness to the user by pulling in the back strap loops.
- **Observe!** Intended use of the Gait trainer support harness is to practise walking. It is **not** a sitting or lifting aid.

### Standing Plate

Attach standing plate by hooks on the main frame of the walker

## Raising/Lowering the Frame

Adjust the height with the aid of the remote control UP and DOWN arrows.

Handset Plus has charge indicator / battery (low, medium, high) and Service indicator (see further under "Hand Control Plus").

## Adjusting the Handles

To adjust handles to the desired position; loosen screw, adjust handle, tighten screw.

## Left and Right Adjustment of Armrests

The armrests can be adjusted sideways for optimal support; loosen screw, adjust armrest, tighten screw.

NOTE: The rear attachments on the frame/armrest are interchangeable (left with right and vice versa) to provide different widths. The locating pins on the cushion's underside/rear can also be moved. These adjustment capabilities make customisation possible.

## Locking/Unlocking Castors wheels

The device is equipped with lockable / braked castor wheels. Lock the wheels by pressing the green pedal brake down with your foot. To release the break, press the green protruding part with your foot.

## Hand Control Plus

\*\*1. Service indicator – flashes red after 11 months' operation.

Constant light (red) after 12 months' operation.

The indicator is a reminder to review the maintenance instruction (see Maintenance Instructions). This should be done at least once every 12 months.

Reset: hold down buttons <2a> + <2b> for at least 30 seconds. \*\*The indicator will go out around 30 seconds after the hand control is activated (battery conservation function).

## Battery Charging

Battery charging must take place when the device is not in use.

Charging starts automatically when power plug is connected to 230V wall socket. Green light on the control box will turn on, once the power plug is connected.

- The LED flashes green during charging.
- ON/OFF: To use the product, press ON.
- The green lamp will begin to flash.
- Use the OFF button (hold down for 5 s) to switch off the battery if the product will not be used for an extended period.
- The LED will show a constant green light when the battery is fully charged.
- When the battery is in normal state/in operation, the LED flashes green at REGULAR intervals, namely 0.5 s ON, 4 s OFF.
- When the battery requires charging (remaining charge <25%) the LED will flash orange/red at REGULAR intervals, namely 0.5 s ON, 4 s OFF.

NOTE: The battery has an in-built device that automatically switches it off when remaining charge is <20%. The LED is then extinguished. The battery must now be fully charged before the electronics can be used again.

## Control box

The control box cover can be opened by using a screwdriver or similar; refer to the instructions which are sitting on the control box.

NOTE: Make sure that plugs are connected to the correct sockets otherwise the actuators can be damaged.

1. Connector for actuator up / down.
2. Hand control connector.
3. Battery outlet.
4. Connector for mains cable. NOTE: The mains cable must

ALWAYS be connected to the control box for IP protection to be valid.

## Trouble shooting

If the lifting function cannot be activated, check the following:

- That all cables are properly and securely connected, Pull out the contact and plug it in again firmly.
- That battery charging is not in progress.
- That the battery is charged.

If the device is not working properly, contact your distributor.

If the device makes unusual noises:

- Try to determine the source of the sound. Take the device out of operation and contact your distributor.

## Accessories

Art. Nr.	Description
56-339-1-0	Handbrake kit/2 handles
56-339-1-1	Handbrake kit/One-handed brake
56-336-4	Oxygen holder
56-337	Side support, adjustable
56-338	Drip rod complete with attachment
56-389	Soft basket with attachment
56-352	Lower Leg Support

## 5. Maintenance

The device must undergo thorough inspection at least once per year. Inspection must be performed by authorised personnel and in accordance with Direct Healthcare Group's instructions.

Repairs and maintenance may only be done by authorised personnel using original spare parts.


Check that	Action if a problem is noticed
<p><b>General</b>            The Walker feels firm/free from play.            The Walker does not rattle during manoeuvres.            The Walker is level and all of its castors are in contact with the floor.            The Walker is not dirty.</p>	<p>Clean with lukewarm soap solution or alcohol-based cleaning agent (no petroleum products).</p>
<p><b>Armrests</b>            Armrests are intact and clean.            Armrest widening is functional.</p>	<p>Clean; replace if damaged.            Replace lock knob.            Fit lock screws/plates as required.</p>
<p><b>Handle</b>            Grips are not damaged/dirty.            Handle adjustment is functional.</p>	<p>Clean with lukewarm soap solution (no petroleum products).            Replace handle grips.            Fit lock screws/plates as required.</p>
<p><b>Electronics</b>            Check that the charging cable is connected to the control box.            Check that the hand control is connected to the control box.            Check that the battery, actuators and control box are not loose.</p>	<p>If any part needs to be replaced refer to the wiring diagram in the control box wiring section.            Fit new charging cable; this must always be connected to the control box.            Connect or fit new hand control as required.            Tighten; replace with new fasteners as required.</p>
<p><b>Height adjustment</b>            Raising and lowering are functional.            The Walker feels stable at maximum settings.            There is no play between the vertical frame and bottom frame.            Height adjustment clamps lock.            The electric actuator fitting must be free from play.            Height adjustment via hand control is functional. The electric actuator must run smoothly at a constant speed.</p>	<p>Fit new guide sleeves (in frame) or end plugs (in chrome tubes) as required.            Tighten all lock bolts.            Replace clamps.            Tighten the attachment fitting concerned; replace bolts/lock nuts as necessary.            Check that the battery/hand control/ actuator are connected in accordance as in the control box wiring section.            Charge the battery.            The part list shows the design and which fasteners must be checked.</p>
<p><b>Frame parts</b>            There is no mechanical damage.            There are no scratches.            End plugs/lower frame fitted</p>	<p>If mechanical damage is present, contact DHG Customer Services.            Touch up as necessary.            Fit new end plugs.</p>
<p><b>Castors/brakes</b>            Castors roll easily/tread not damaged.            The castors are firmly fastened to the lower frame.            Castor brakes functional on all castors.            Handbrake functional.</p>	<p>Clean or replace castors.            NOTE! Castors are always fixed to the lower frame with thread locker or lock nuts. The castors are sealed and we do not recommend their dismantling; instead, replace the whole castor assembly.            Tighten the castor bolt and fit locking nut or use thread locker (depends on model).            Replace castor assembly.            Adjust the brake or fit new castors.</p>

Power cord may only be replaced by DHG service personnel or technical personnel trained by DHG.

DHG can provide upon request necessary technical information in order to assist maintenance and repair of the device.


## Storage and transportation

The walker should be transported and stored in 5 °C to 50 °C and in normal humidity, 20% to 80% non-condensing. The air pressure should be between 700 and 1060 hPa. Leftmost symbol indicates storage and transportation. Let the walker reach room temperature before the batteries are charged or the walker is used. The walker should not be stored so that it is exposed to dust, or so that the battery or gas piston is exposed to direct sunlight.

## Operation

The operating environment should be 5 °C to 40°C, relative humidity 20% to 80% non-condensing, and atmospheric pressure 700 to 1060 hPa. See also Technical specification below. Leftmost symbol indicates operating condition.

## Cleaning and disinfection

After every use:

Clean the walker armrests and handles. Wipe using a mild soap solution/surface disinfectant. However, the product is not corrosion proof.

Weekly and possibly more often depending on frequency of use:

If necessary, clean the whole walker with a cloth with warm water or a soap solution and check that the castors are free from dirt and hair. The electronic components (battery, control box, hand control, actuators, cables) should be cleaned with a damp cloth only. To avoid degreasing of the piston rods, the actuators should be retracted to minimum stroke and without load before cleaning. Ensure that the device is dried thoroughly after cleaning.

If disinfection is needed:

70% ethanol, 45% isopropanol or similar should be used.

Do not steam clean due to risk of corrosion.

Do not use cleaning agents containing phenol or chlorine, as this could damage the materials.

With environments with high risk of infections use the walker with PU-cushions (Accessory number 57-361-L and 57-361-R)

## Recycling/Disposal


When the product is no longer usable it should be recycled in accordance with legislation and regulations in the country concerned. All electrical parts including batteries must be removed and recycled as electrical components. Ask your local recycling station for further information about how the different materials in the product (metals, plastics, electronics) should be recycled. Used batteries can also be returned to Direct Healthcare Group or a Direct Healthcare Group distributor for recycling.

## 6. Technical Specifications

<b>Control box</b>	TiMotion, Voltage 220V-240VAC, 50Hz, 1250mA/24VDC, 70VA, IP66, Class II (Doubled insulated)
<b>Batteries</b>	TiMotion, 24 V - 4A, 96Wh, Rechargeable Li-Ion Battery, IP66, Cable length 600mm, Duty cycle: 10%, Max. 2min. ON/18 min./OFF
<b>Actuator</b>	12VDC, Max 8,0A, Push: 2000N, Self-locking capacity (min): 1000N, Stroke 400mm +/- 2 mm, IP66. Duty cycle: 10%, Max. 2min.ON/18min. OFF
<b>Hand Control</b>	Cable length straight 300mm + coiled 400 mm, IP66
<b>IP class</b>	IPX5
<b>Material</b>	Steel
<b>Castors</b>	125/100/75mm
<b>Weight</b>	20 kg
<b>Expected Service life</b>	8 years
<b>Max patient weight</b>	240 kg / 529 lbs
<b>Operating environment</b>	5 °C to 40°C, relative humidity 20% to 80% non-condensing, and atmospheric pressure 700 to 1060 hPa.
<b>Storage and transportation environment</b>	5 °C to 50 °C, relative humidity 20% to 80% non-condensing, and atmospheric. Environment pressure 700 to 1060 hPa.
<b>Height</b>	125 mm wheels model minimum 900 mm and maximum 1300 mm 100 mm wheels model minimum 880 mm and maximum 1280 mm 75mm wheels model minimum 850 mm and maximum 1250 mm
<b>Length</b>	850 mm
<b>Width</b>	775 mm
<b>Turning Width</b>	125 mm wheels: 1158 mm, 100 mm wheels: 1130 mm, 75 mm wheels: 1130 mm
<b>Materials in contact with human skin</b>	Materials that user and operator/caregiver are exposed to during usage of the product for prolonged time: Poroflex (synthetic leather cushions) P <sup>x</sup> Polyurethane (on PU cushions) P Dryflex (handles) P PA6 + 20% GF (brake lever) P Hand control C <sup>x</sup>
<b>Applied parts:</b>	Armrest P Handles P Knee pad P Foot rest support P Foot rest support clamps C Handles adjusting knobs C Hand controller and its holder. C Handles to adjust the armrest width. C Hand control C Charger C Charger cable C Footplate


<sup>x</sup> P - patient, <sup>x</sup> C - caregiver

## 7. Parts

NO.	PART NO.	DESCRIPTION	NO.	PART NO.	DESCRIPTION
1	58-373	Control box 220V Digital Li-On TiMotion Bure	16	57-810	Screw MC6S M8x45
2	58-377	Actuator T-Motion 2000N/5200rpm, Bure	17	57-811	Serrated lock washer M16
3	58-389	Battery box 4.0 Ah TBB5/Li-On	18	57-917	Wing nut
4	58-335	Telescopic Inner pipe	19	57-925	Endcover pipe 39mm black
5	57-421	Telescope inner pipe Bushing	20	56-308	Cable holder Bure charger
6	57-422	Telescope inner pipe Tube	21	57-808	Pin bolt M6x25 8.8
7	57-747	Plug dome plastic 25dia	22	57-792	Washer Cushion
8	57-049	Castor wheel 125mm brake	23	57-848	Nordlock 10,7x21x2,5
9	57-307	Protective cover plastic	24	57-918	Washer 6,5x15x10
10	57-310	Middle frame bushing	25	57-357-L	Armrest Cushion Low Left
11	57-803	Locking washer M10 6.1x11.8x1.6	26	57-357-R	Armrest Cushion Low Right
12	57-804	Screw MC6S M8x35	27	57-919	Star knob handle
13	57-805	Washer BRB 6,4x12x1,6	28	58-331	Cover actuator rod
14	56-353	Bure Handle assy	29	58-374-X1	Hand control TiMotion TH10
15	57-806	Washer 8.4x16x1.5	30	58-326	Spiral Cord EUR - T-Motion, Bure

OPTIONS		
NO.	PART NO.	DESCRIPTION
1	57-361-L	PU Arm support cushion Bure left
2	57-361-R	PU Arm support cushion Bure right
3	57-051	Wheel w/o brake 75mm
4	56-387	Wheel with brake 75 mm
5	57-319	Wheel with brake 100mm
6	56-357	Wheel w/o brake 100mm
7	56-360	Wheel w/o brake 125mm
8	56-339-1-0	Handbrake kit/2 handles
8	56-339-1-1	Handbrake kit/One-handed brake
9	Multiple options	Footplate Bure Ortho/Space
	56-305	Touch-up paint Bure, Silver 8
	58-306	Touch-up paint Bure, Dark grey RAL 7021


## 8. Product label

# Bure Space EI, Timotion SystemRoMedic®

**REF** 56-304

**LOT** 2022-03-R

(01)07350013304019(11)220301(21)610231377

 150 kg / 330 lbs    230 VAC, 50/60 Hz    IPX5      

 775mm

 900-1300 mm

**Direct Healthcare Group Sverige AB**

Torshamnsgatan 35, 16440 Kista, SWEDEN

[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)

Made in Sweden

Rev .01

UDI (Unique Device Identification) can be found within the barcode human readable text. It contains:-




(01) EAN code

(11) Date of Manufacture. The date is in the format YYMMDD Y = Year, M = Month and D = Date.

(21) Serial Number displayed as five digit code.

## 9. Table of Symbols

	This product complies with the requirements of the Medical Device Regulation 2017/745		Visual Inspection
	Medical Device		Read the manual
	Caution		Product Code
	Legal manufacturer Direct Healthcare Group Sverige AB		Batch Code
	Type BF applied parts, according to the degree of protection against electric shock		The device is intended for indoor use
	WEEE Symbol May not be discarded in domestic waste		Class II Equipment
	Duty Cycle: 1 min in active (ON) mode. 19 min in rest (OFF) mode.		Maximum patient weight
	Walker Height		Walker Width
	Do not step on the device		Refer to instruction manual (IFU)
		Mass Weight (mass) of the device, the Maximum Patient Weight, and total sum. All in kg.	

## 10. How to report a serious incident

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the MHRA or another competent authority of the country in which the user and/or patient is established.

UK

info@directhealthcaregroup.com


T: +44 (0) 800 043 0881

F: +44 (0) 845 459 9832

Other:

info.export@directhealthcaregroup.com

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
## 8. Produktetikett

**Bure Space EI, Timotion**  
System**RoMedic**

**REF 56-304**

**LOT 2022-03-R**

**(01)07350013304019(11)220301(21)610231377**

 **150 kg / 330 lbs**   **230 VAC, 50/60 Hz**   **IPX5**     

 **775mm**

 **900-1300 mm**

 **Direct Healthcare Group Sverige AB**  
Torshamnsgatan 35, 16440 Kista, SWEDEN  
[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)  
Made in Sweden

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

SE:

info.se@directhealthcaregroup.com

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

# Electromagnetic Compatibility (EMC)

## Electromagnetic emissions and test levels

The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Emission test	Standards	Compliance
RF emissions	CISPR 11	Group 1
RF emissions	CISPR 11	Class B
Harmonic current emissions	IEC 61000-3-2	Class A
Voltage fluctuations and Flicker emissions	IEC 61000-3-3	Complies

## Electromagnetic Immunity

Acceptance criteria for the EMC pass criteria is unintentional movement above  $\pm 10\text{mm}$  is not allowed. The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Acceptance criteria for the EMC pass criteria is unintentional movement above  $\pm 10\text{mm}$  is not allowed. The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Radiated Fields in close proximity, Immunity test and compliance. Standard: IEC 61000-4-39			
Dwell time: 3 sec			
Test levels (A/m)	Modulation Pulse modulation	Mod. Frequency (kHz)	Test Frequency
8	CW	CW	30 kHz
65	50% Duty cycle	2.1	134.2 kHz
7.5	50% Duty cycle	50	13.56 MHz

Immunity test	Basic EMC standard or test method	IEC 60601-1-2 Edition 4 Test levels and compliance	
		Professional healthcare facility environment†	Home healthcare environment
Surge	IEC 61000-4-5	± 0.5 kV, ± 1 kV	
Voltage dips	IEC 61000-4-11	0% U <sub>T</sub> ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
		0% U <sub>T</sub> ; 1 cycle and 70% U <sub>T</sub> ; 25/30 cycles. Single phase: at 0°	
Voltage Interruptions	IEC 61000-4-11	0% U <sub>T</sub> ; 250 / 300 cycle	
Conducted disturbances induced by RF fields	IEC 61000-4-6	6V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	6V in ISM and amateur bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
		Dwell time: 3 sec Frequency step size: 1%	
Radiated RF Electro-magnetic fields	IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz
		Dwell time: 3 sec Frequency step size: 1%	
Electrical fast transients / burst	IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	
RATED Power frequency magnetic field	IEC 61000-4-8	30 A/m 50 Hz or 60Hz	
Electrostatic Discharge (ESD)	IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ±4 kV, ± 8 kV, ± 15 kV Number of discharges: 10 for each polarity	

Radiated Radio – Frequency, Proximity fields from wireless communication equipment & compliance. (Dwell time: 3 sec.)  
 Standard: IEC 61000-4-3

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test levels (V/m)
Spot Frequencies	385	380 to 390	TETRA 400	Pulse modulation 18 Hz	27
	450	430 to 470	FMRS 460, FRS 460	FM ± 5 kHz Deviation 1 kHz sine	28
	710	704 to 787	LTE Band 13, 17	Pulse modulation 217 Hz	9
	745				
	780				
	810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
	870				
	930				
	1720	1700 to 1990	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 217 Hz	28
	1845				
	1970				
	2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
	5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
	5500				
	5785				
Other Identified frequencies	433	-	-	Pulse modulation 2 Hz	3
	-	860-960	-	Pulse modulation 2 Hz	54
	2450	-	-	Pulse modulation 2 Hz	54

## Simple solutions for great results

SystemRoMedic® is the name of Direct Healthcare Group's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic® is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

## Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic® is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic® offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

## Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading at/from our website;  
[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)

## For questions about the product and its use

Please contact your local Direct Healthcare Group and SystemRoMedic® representative. A complete list of all our partners with their contact details can be found on our website; [www.directhealthcaregroup.com](http://www.directhealthcaregroup.com).

Moving Health Forward

Direct Healthcare Group Sverige AB  
Torshamnsgatan 35,  
SE-164 40 Kista, Sweden

Tel: +46 (0)8-557 62 200  
[info@directhealthcaregroup.com](mailto:info@directhealthcaregroup.com)  
[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)

Medical Device Class I. The product complies with the requirements of the Medical Devices Regulation 2017/745.