Bure Standard/S Bure Low/S

SystemRoMedic™

Instructions for use – English
Bruksanvisning – Svenska
Brukermanual – Norsk
Brugsvejledning – Dansk
Käyttöohje – Suomi
Gebrauchsanweisung – Deutsch
Handleiding – Nederlands
Manuel d'utilisation – Français
Istruzioni per l'uso - Italiano
Naudojimo instrukcija – lietuvių

Max Patient Weight: 150 kg / 330 lbs

Table of contents

English	4
Svenska	19
Norsk	34
Dansk	49
Suomi	64
Deutsch	79
Netherlands	95
Français	111
Italiano	127
lietuvių	143
Dimensions	159
Parts List	160
Product Overview	164

Bure Standard/S

Bure Low/S

SystemRoMedic

Instructions for use - English

Bure Standard is a walker which is an excellent aid to patients in situations like standing up and walking. Bure Standard is available with manual, electrical or hydraulic operating height adjustment.

1. Intended use

Walkers are to be used as indoor, movable and portable assistive aids, which together with approved accessories are intended to be used for transfer from sitting to standing positions, and walking training. Walkers are to be used by persons with reduced muscular strength, who shall be able to stand and perform simple leg movements. Walkers are to be used for leaning on when walking and standing during early mobilisation training. Walkers and accessories are intended for indoor use only in Home Healthcare environment and Professional Healthcare facility environment.

2. Contraindications, Precautions, Warnings

Contraindications

- The device may not be used by patients above the maximum weight indicated on the label of the device.

Precautions

- Check that the device is used on a dry and flat surface.
- Check the device is correctly mounted/assembled before its first use.
- Check the device after every folding/disassembly after any transport.
- Check height adjustment and leg spreading motions and inspect the actuators or gas pistons full ranges.
- Activate the brakes whenever the device is not in use.
- Check the function of the brakes before use.
- Inspect the device regularly to detect any signs of damage. If there are any visible signs of damage, contact Direct Healthcare Group Customer Service.
- In case of an electrical model of the device, check that the hand control does not show signs of wear before every use.
- In case of an electrical model of the device, check that the hand control markings are in accordance with the lifting functions.
- Check the battery status of an electrically operated model before using the device.
- Warranty applies only if repairs or alterations are made by personnel who are authorized by Direct Healthcare Group.
- Ensure that there are no obstacles or people in the way of the device, when moving.
- Handle batteries with care. Do not drop.
- Use only batteries and cables that are intended for the device as instructed by the manufacturer.
- Verify that potential accessories that might be used are properly attached to the device before use.
- Low speed is recommended when moving with the device.
- Take care not to drive the device over thresholds with high speed or force. Approach the threshold with lower speed and communicate the coming threshold to the patient.
- Keep the instructions for use together with the product so that it can be used as a reference if necessary.
- Beware of moving parts during usage.
- Perform a risk assessment on the patient before using the Walker to ensure they can stand and perform simple movements. Make sure the user can stand and make simple movements with legs.

Warnings

- Excluding Bure Low, the device may not be used by infants, children or adults who are 146 cm or shorter.
- The caregiver must be able to read and understand the manual/Instructions for use of the product.
- Caregiver: Never use the device if you have not had a proper training on how to use the device including the accessories.
- Maintenance/Service staff: Never maintain the device if you have not had a proper training on how to use and maintain the device including the accessories.
- It is important to use only approved accessories to prevent unintended detachment of components and subsequently a fall that may lead to patient injury.
- Use careful and gentle maneuvers when moving the device.
- Perform maintenance/service of the device, according to the instructions in the Instruction for use, at least once every 12 months.
- No parts of the device shall be serviced, while the device is in use
- Accessories must be properly fitted and tested in relation to the patient's needs and functional ability.
- Special care must be taken when using strong electrical power sources such as diathermy and the like so that diathermy cables are not placed on or near the device. In case of doubt consult with a Direct Healthcare Group representative.
- When lifting the device, only use the denoted handholds. At least two people are required for lifting the device.
- Do not leave a patient unattended when the device is being used.
- The device can not be used by patient alone.
- Never move the walker by pulling on the actuator device or gas piston.
- The device must not be lowered into water.
- The device must not be cleaned using steam or washed in a washing tunnel.
- The device must not be left or stored in a damp or humid environment.
- The device must not be charged in a wet room.
- The device must not be used outdoors, only indoors.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective grounding.
- Do not modify this equipment without authorisation of the manufacturer.
- If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of this equipment.
- Do not use the device while it is charging.
- Do not use the device if any of the cables shows signs of wear and tear.
- Do not put the device in direct sunlight as it may get hot.
- Do not run over the cables with the walker or other devices.

Always read the user manual

Always make sure that you have the correct version of the manual. The most recent version of all manuals are available for downloading at/from our website: www.directhealthcaregroup.com.

Before use:

- Inspect the packaging for any damage.
- Check for correct product being delivered.
- Check to ensure that the device is free from damage.
- Check if all the parts/components are included in the packaging.
- Check if all approved accessories are included in the package.
- Check if Product Quality Approval document is included in the packaging. Save this document for future contact with manufacturer.
- Make certain the device is properly assembled.
- Check lifting function and bottom-frame adjustment function.
- Check if it is possible to move the device in all directions.
- Check the condition and function of the sling/belt regularly. Always inspect the product after laundering. The occurrence of the following must be checked: frays, loose stitches, damage to the fabric, edging, handles, broken buckles. If any of these points apply, the product should be discarded.

3. Using Bure Walker

Standing up and walk training, with electrical motor:


1. Ask the patient to sit on the edge of the bed or chair. Bring out the walking table and place it in front of the patient. Lock the brakes on the wheels.
2. Ask the patient to place their arms on the armrests. If necessary, ask the patient to move forward on the seat base so that the patient reaches the handles.
3. Ask the patient to lean forward. To start raising the walking table, press the UP button on the controller. Raise the walking table to a desired height.
4. Unlock the brakes. Adjust the height of the walking table to a desired position
5. Assist the patient while walk training.

With gas actuator or manual:

1. Ask the patient to sit on the edge of the bed or chair. Bring out the walking table and place it in front of the patient. Lock the brakes on the wheels.
2. Raise the walking table by:
 - a. Pressing the handle on the gas piston.
 - b. Opening handle, adjusting the height, and closing the handle.
3. Ask the patient to lean forward, grasp the armrests and stand up.
4. Adjust the height of the walking table to a desired position
5. Assist the patient while walking.

Using Bure Walker with Gait Training Kit 56-388:

- The support harness has four blue straps with a plastic buckle allowing for adjustments in length, having steel carabiner clips at the ends to be fasten to the gait walker. The carabiner clips on to the metal mounts placed under the table of the Gait walker.
- Apply the walker support harness by clipping on the two front carabiners to the front mounters, one on each side of the table.
- Place the patient standing up in the walker. Pull back the harness in between the legs and clip on the carabiners to the back mounters, one on each side of the table. Adjust the harness to the patient by pulling in the back strap loops.
- **Observe!** Intended use of the trainer support harness is to practise walking. It is **not** a sitting or lifting aid.

Standing Plate

Attach standing plate by hooks on the main frame of the walker.

Raising/Lowering the Frame, Bure electrical actuator

Adjust the height with the aid of the hand control UP and DOWN arrows.

Hand Control has charge indicator / battery (low, medium, high) and Service indicator.


Service indicator – flashes red after 11 months' operation.

Constant light (red) after 12 months' operation. This indicates need of yearly maintenance.

Raising/Lowering the Frame - Bure gas actuator


Release the gas strut by moving the handle upward (the gas strut will raise the Walker). Let go of the handle at the desired height (the gas strut a locked).

Release the gas strut by moving the handle upward. Use your body weight to press the frame down. Let go of the handle at the desired height (the gas strut will lock).

Adjusting the Handles

To adjust handles to the desired position; loosen screw, adjust handle, tighten screw.

Adjustment of Armrests

The armrests can be adjusted sideways for optimal support; loosen screw, adjust armrest, tighten screw.

NOTE: The rear attachments on the frame/armrest are interchangeable (left with right and vice versa) to provide different widths. The locating pins on the cushion's underside/rear can also be moved. These adjustment capabilities make customisation possible.

Using the Parking Brakes on Castor Wheels:

When parking the device, activate the parking brake on the lockable castor wheels. Lock the wheels by pressing the green pedal brake down with your foot. To release the brake, press the green protruding part with your foot. The caregiver should assist the patient in activating the parking brake in case the patient cannot activate them him or herself.

Using Hand Brakes:

Running brakes: Pull the brake levers towards you to reduce speed. You will brake the rear wheel that is located on the same side as the brake lever. If you apply only one of the brakes, the walker will turn to that side. When the lever is released, the braking will stop. Do not use the brakes continuously while the walker is in motion.

Parking brakes: Push the brake levers forwards from you to activate the parking brakes. Be aware that the lever has some resistance before it stops and locks the rear wheel. Pull the brake levers back towards you to release the parking brakes.

Battery Charging

Battery charging must take place when the device is not in use.

Charging starts automatically when power plug is connected to 230V wall socket.

- LED lights YELLOW during charging (at the same time as the charging indicator on the hand control flashes yellow).
- The LED goes out when the battery is fully charged (the same applies to the charge indicator on the hand control).


Control box

Control box - LED lights up yellow during operation, is off when the table is not in use. When the control box is connected to the mains, the LED lights up green. When driving / connected to the mains, the LED lights up yellow. The door on the control box can be opened with a thumb grip.

1. Hand control socket.
2. Battery outlet.

(NOTE: Plugs must always be connected to the control box for IP classification to be valid.)

3. Socket for actuator up / down.
4. Socket for plug.

Trouble shooting

If the lifting function cannot be activated, check the following:

- That all cables are properly and securely connected. Pull out the contact and plug it in again firmly.
- That battery charging is not in progress.
- That the battery is charged.

If the device is not working properly, contact your distributor.

If the device makes unusual noises:

- Try to determine the source of the sound. Take the device out of operation and contact your distributor.


4. Maintenance

- Before maintenance, clean walker according to cleaning instructions.
- The device must undergo thorough inspection at least once per year. Inspection must be performed by authorised personnel and in accordance with Direct Healthcare Group's instructions.
- Repairs and maintenance may only be done by authorised personnel using original spare parts.
- Power cord may only be replaced by DHG service personnel or technical personnel trained by DHG.
- DHG can provide, upon request necessary technical information in order to assist maintenance and repair of the device.

Check that	Action if a problem is noticed
Frame	
There is no mechanical damage.	If mechanical damage is present, contact DHG Customer Services.
There are no scratches.	Touch up as necessary.
End plugs/lower frame fitted and not damaged.	Fit new end plugs.
Armrest	
Armrests are intact	Replace if damaged.
Armrest widening is functional.	Replace lock knob, see parts section.
Handle	
Grips are not damaged.	Replace handle grips.
Handle adjustment is functional.	Adjust lock screws/plates as required.
Electronics	
Check that the charging cable is connected to the control box.	Fit new charging cable; this must always be connected to the control box.
Check that the hand control is connected to the control box.	Connect or fit new hand control as required. Tighten; replace with new fasteners as required.
Check that the battery, actuators and control box are not loose.	If any part needs to be replaced refer to the wiring diagram in the control box wiring section.
Height adjustment	
Raising and lowering are functional.	Fit new guide sleeves (in frame) or end plugs (in chrome tubes) as required.
The walker feels stable at maximum height and width.	Tighten all nuts and bolts.
There is no movement between the vertical frame and bottom frame.	Tighten all nuts and bolts.
Height adjustment clamps lock.	Replace clamps.
The electric actuator fitting must be free from movement.	Tighten all nuts and bolts; replace bolts and nuts if necessary.
Height adjustment via hand control is functional. The electric actuator must run smoothly at a constant speed.	Check that the battery/hand control/ actuator are connected according to the control box wiring section. Charge the battery. See parts section
Castors and brakes	
Castors roll easily, tread not damaged.	Replace castors.
The castors are firmly fastened to the lower frame.	Tighten the castor bolt and nut. Use thread locker.
Castor brakes functional on all castors. Handbrake functional.	Replace castors or handbrake assembly. Clean and adjust the brake. The castors are sealed and it is not recommended to dismantle them; instead, replace the whole castor assembly.


Storage and transportation

The walker should be transported and stored in 5°C to 50°C and in normal humidity, 20% to 80% non-condensing. The air pressure should be between 700 and 1060hPa. Leftmost symbol indicates storage and transportation. Let the walker reach room temperature before the batteries are charged or the walker is used. The walker should not be stored so that it is exposed to dust, or so that the battery or gas piston is exposed to direct sunlight.

Operation

The operating environment should be 5°C to 40°C, relative humidity 20% to 80% non-condensing, and atmospheric pressure 700 to 1060 hPa. See also Technical specification below. Leftmost symbol indicates operating condition.

Cleaning and disinfection

After every use:

Clean the walker armrests and handles. Wipe using a mild soap solution/surface disinfectant. The product is not corrosion proof.

Weekly and possibly more often depending on frequency of use:

If necessary, clean the whole walker with a cloth with warm water or a soap solution and check that the castors are free from dirt and hair. The electronic components (battery, control box, hand control, actuators, cables) should be cleaned with a damp cloth only. To avoid degreasing of the piston rods, the actuators should be retracted to minimum stroke and without load before cleaning. Ensure that the device is dried thoroughly after cleaning.

If disinfection is needed:

70% ethanol, 45% isopropanol or similar should be used.

Do not steam clean due to risk of corrosion.

Do not use cleaning agents containing phenol or chlorine, as this could damage the materials.

In environments with high risk of infections use the walker with PU-cushions (Accessory number 57-361-L and 57-361-R)

General Maintenance for gait training Sling:

Read the product label.

Do not use rinsing agent.

To prolong product life, avoid tumble-drying.

Recycling/Disposal

When the product is no longer usable it should be recycled in accordance with legislation and regulations in the country concerned. All electrical parts including batteries must be removed and recycled as electrical components. Ask your local recycling station for further information about how the different materials in the product (metals, plastics, electronics, textiles) should be recycled. Used batteries can also be returned to Direct Healthcare Group or a Direct Healthcare Group distributor for recycling.

5. Technical Specifications

Control box	Linak, voltage 100V-240V, Max 1.5A, 50-60 Hz, IPX6, Class II (Double insulated)
Batteries	Linak, 25,2V/2,1Ah, sealed,Li-Ion, 300mA IPX6, Cable length 200mm
Actuator	Linak, DC 24 V, 4,3A, IPX6. Duty cycle: Max 10% or 2 minutes continuous use followed by 18 minutes not in use, Push: 1500N. Stroke 400mm +/-2 mm
Hand Control	Cable length coiled 600 mm, IPX6
Weight	STD- 24 kg STD S- 23 kg STD LOW- 15 kg
IP class	IPX5
Material	Steel
Castors	125/100/75mm
Dimensions	For dimensions, see "Dimensions" page
Expected Service life	Walkers: 8 years Textile products: 2 years under normal use. Two wash cycles a week are considered as normal use. However, the total lifetime is dependent on the use, maintenance and storage of the product.
Max patient weight	150 kg / 330 lbs
Operating environment	5 °C to 40°C, relative humidity 20% to 80% non-condensing, and atmospheric pressure 700 to 1060 hPa.
Storage and transportation environment	5 °C to 50 °C, relative humidity 20% to 80% non-condensing, and atmospheric. Environment pressure 700 to 1060 hPa.
Materials in contact with human skin	Materials that patient and operator/caregiver are exposed to during usage of the product for prolonged time: Meditap B-365 (synthetic leather cushions) P ^x Polyurethane (on PU cushions) P Dryflex (handles) P PA6 + 20% GF (brake lever) P Hand control C ^x
Applied parts:	Armrest P Handles P Knee pad P Foot rest support P Foot rest support clamps C Handles adjusting knobs C Hand controller and its holder. C Handles to adjust the armrest width. C Hand control C Charger C Charger cable C Footplate P

^x P - patient, ^x C - caregiver

Accessories

Art. No.	DESCRIPTION
58-328-1	Heel strap Standing plate
56-336-4	Oxygen holder
56-338	Drip rod complete with attachment
56-337	Side support, adjustable
56-337-PU	Side support, adjustable PU
56-352	Lower Leg Support
56-333-STD	Standing plate Short, Bure Std, S, Rise
56-333/SW-STD	Standing plate Short, SW, Bure Std, S, Rise
56-383-STD	Standing plate Long, Bure Std, S, Rise
56-383/SW-STD	Standing plate Long, SW, Bure Std, S, Rise
56-339-1-0	Handbrake kit, 2 handles, Bure
56-339-1-1	Handbrake kit, 1 handle, Bure
56-339/75-1-0	Handbrake kit, 75mm wheel, 2 handles
56-339/75-1-1	Handbrake kit, 75mm wheel, 1 handle
57-368	Guide handle standard
56-388	Gait training kit
56-389	Soft basket with attachment, Bure Standard/S
56-305	Touch-up paint Bure, Silver 8
56-306	Touch-up paint Bure, Dark grey RAL 7021

6. Product label

Bure STD Linak, 125 mm SystemRoMedic

REF 56-312L

(01)07331769044879(11)220301(21)610231377

MD CE

LOT 2022-03-R

UDI (Unique Device Identification) can be found within the barcode human readable text. It contains.

(01) EAN code

(11) Date of Manufacture. The date is in the format YYMMDD

Y = Year, M = Month and D = Date.

(21) Serial Number displayed as five digit code.

150 kg / 330 lbs

U In: 220-240 VAC, 50Hz

I In: Max 1,5A

775 mm


IPX5

Direct Healthcare Group Sverige AB
Torshamnsgatan 35, 16440 Kista, SWEDEN
www.directhealthcaregroup.com




Made in Sweden

Rev.02

910-1310 mm

7. Table of Symbols

	This product complies with the requirements of the Medical Device Regulation MDR 2017/745		Visual Inspection
	Medical Device		Read the manual
	Caution		Product Code
	Legal manufacturer Direct Healthcare Group Sverige AB		Batch Code
	Type BF applied parts, according to the degree of protection against electric shock		The device is intended for indoor use
	WEEE Symbol May not be discarded in domestic waste		Class II Equipment
	Duty Cycle: 2 min in active (ON) mode. 18 min in rest (OFF) mode.		Maximum patient weight
	Walker Height		Walker Width
	Do not step on the device		Refer to instruction manual (IFU)
	 $+$  $=$ 	Weight (mass) of the device, the Maximum Patient Weight, and total sum. All in kg.	

8. How to report a serious incident

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the MHRA or another competent authority of the country in which the user and/or patient is established.

UK

info@directhealthcaregroup.com

T: +44 (0) 800 043 0881

F: +44 (0) 845 459 9832

Other:

info.export@directhealthcaregroup.com

Electromagnetic Compatibility (EMC)

Electromagnetic emissions and test levels

The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Emission test	Standards	Compliance
RF emissions	CISPR 11	Group 1
RF emissions	CISPR 11	Class B
Harmonic current emissions	IEC 61000-3-2	Class A
Voltage fluctuations and Flicker emissions	IEC 61000-3-3	Complies

Electromagnetic Immunity

Acceptance criteria for the EMC pass criteria is unintentional movement above $\pm 10\text{mm}$ is not allowed. The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Radiated Fields in close proximity, Immunity test and compliance. Standard: IEC 61000-4-39			
Dwell time: 3 sec			
Test levels (A/m)	Modulation Pulse modulation	Mod. Frequency (kHz)	Test Frequency
8	CW	CW	30 kHz
65	50% Duty cycle	2.1	134.2 kHz
7.5	50% Duty cycle	50	13.56 MHz

Immunity test	Basic EMC standard or test method	IEC 60601-1-2 Edition 4 Test levels and compliance	
		Professional healthcare facility environment	Home healthcare environment
Surge	IEC 61000-4-5	± 0.5 kV, ± 1 kV	
Voltage dips	IEC 61000-4-11	0% U _T ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
		0% U _T ; 1 cycle and 70% U _T ; 25/30 cycles. Single phase: at 0°	
Voltage Interruptions	IEC 61000-4-11	0% U _T ; 250 / 300 cycle	
Conducted disturbances induced by RF fields	IEC 61000-4-6	6V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	6V in ISM and amateur bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
		Dwell time: 3 sec Frequency step size: 1%	
Radiated RF Electromagnetic fields	IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz
		Dwell time: 3 sec Frequency step size: 1%	
Electrical fast transients / burst	IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	
RATED Power frequency magnetic field	IEC 61000-4-8	30 A/m 50 Hz or 60Hz	
Electrostatic Discharge (ESD)	IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ±4 kV, ± 8 kV, ± 15 kV Number of discharges: 10 for each polarity	

Radiated Radio – Frequency, Proximity fields from wireless communication equipment & compliance. (Dwell time: 3 sec.) Standard: IEC 61000-4-3

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test levels (V/m)
Spot Frequencies	385	380 to 390	TETRA 400	Pulse modulation 18 Hz	27
	450	430 to 470	FMRS 460, FRS 460	FM ± 5 kHz Deviation 1 kHz sine	28
	710	704 to 787	LTE Band 13, 17	Pulse modulation 217 Hz	9
	745				
	780				
	810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
	870				
	930				
	1720	1700 to 1990	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 217 Hz	28
	1845				
	1970				
	2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
	5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
	5500				
5785					
Other Identified frequencies	433	-	-	Pulse modulation 2 Hz	3
	-	860-960	-	Pulse modulation 2 Hz	54
	2450	-	-	Pulse modulation 2 Hz	54

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



Distinguishing parts for other series


Smal/Narrow 56-302 series		
No.	Art No.	Description
1	57-300-2	Bottom frame Bure Narrow
2	57-302-3	Armrest frame Bure Narrow


Low 56-321 series		
No.	Art No.	Description
1	57-607-LOW	Telescope inner pipe LOW
2	57-301-LOW	Middle frame Bure STD LOW
3	58-401	Actuator Linak LA31 Low 1500N

56-322 Low Gas		
No.	Art No.	Description
1	57-607-LOW	Telescope inner pipe LOW
2	57-301-LOW	Middle frame Bure STD LOW
3	57-007-LOW	Gas spring 150N Bure Low

Manual 56-310 series		
No.	Art No.	Description
1	57-301-1	Middle frame Bure STD manual
2	57-352	Locking lever
3	-	Clip height adjustment Bure Manual
4	57-718	Screw K6S M8x40

Gas 56-311-series		
No.	Art No.	Description
1	57-329	Throttle handle Bure Gas
2	57-007	Gas spring 150N Bure
3	57-008-3	Gas spring release housing
4	57-806	Washer 8.4x16x1.5
5	57-807	Screw MC6S M8x35
6	57-880	Cap nut CL6 M8
7	57-846	Screw MC6S M10x35
8	57-825	Washer 10,5x18x2
9	57-817	Lock nut CL6 M10

Product Overview

REF

Article number	Product name
56-300	Bure S Manual, 125 mm
56-300-B/0	Bure S Manual, 125 mm, handbrakes
56-300-B/1	Bure S Manual, 125 mm, 1 handbrake
56-300/100	Bure S Manual, 100mm
56-300/75	Bure S Manual, 75mm
56-300/75-DIR	Bure S Manual, 75mm, locked direction
56-300/75-B/0	Bure S Manual, 75mm, handbrakes
56-300/75-B/0-PU	Bure S Manual, 75mm, handbrakes, PU
56-300/75-B/1-PU	Bure S Manual, 75mm, 1 handbrake, PU
56-301	Bure S Gas, 125 mm
56-301-B/0	Bure S Gas, 125 mm, handbrakes
56-301-PU	Bure S Gas, 125mm, PU
56-301/75	Bure S Gas, 75mm
56-301/75-PU	Bure S Gas, 75mm, PU
56-301/75-RS	Bure S Gas, 75mm, locked direction
56-310	Bure Std Manual, 125mm
56-310-B/0	Bure Std Manual, 125mm, handbrakes
56-310/100	Bure Std Manual, 100mm
56-310/75	Bure Std Manual, 75mm
56-310/75-B/0	Bure Std Manual, 75mm, handbrakes
56-310-PARSONS-KD	Bure Std Manual, 125mm, KnockDown
56-311	Bure Std Gas, 125mm
56-311-PARSONS-KD	Bure Std Gas, 125mm, KnockDown, 2 direction locks
56-311-B/1	Bure Std Gas, 125mm, 1 handbrake
56-311-BRAKE/0-PU	Bure Std Gas, 125mm, handbrakes, PU
56-311-PU	Bure Std Gas, 125mm, PU
56-311-RS	Bure Std Gas, 125mm, 1 direction lock
56-311/100	Bure Std Gas, 100mm
56-311/100-PU	Bure Std Gas, 100mm, PU
56-311/75	Bure Std Gas, 75mm
56-311/75-B/0	Bure Std Gas, 75mm, handbrakes
56-311/75-B/1	Bure Std Gas, 75mm, 1 handbrake
56-311/75-PU	Bure Std Gas, 75mm, PU

Article number	Product name
56-322/75	Bure Low Gas, 75mm
56-322-S	Bure Low S Gas, 125mm
56-322-S-B/0	Bure Low S Gas, 125mm, handbrakes
56-322-S/75	Bure Low S Gas, 75mm
56-302L	Bure S Linak, 125 mm
56-302L-PU	Bure S Linak, 125mm, PU
56-302L-B/0	Bure S Linak, 125mm, handbrakes
56-302L-B/0-PU	Bure S Linak, 125mm, handbrakes, PU
56-302L-B/1	Bure S Linak, 125mm, 1 handbrake
56-302L-B/1-PU	Bure S Linak, 125mm, 1 handbrake, PU
56-302L-UK	Bure S Linak, 125 mm, UK
56-302L-PU-UK	Bure S Linak, 125mm, PU, UK
56-302L-B/0-UK	Bure S Linak, 125mm, handbrakes, UK
56-302L-B/0-PU-UK	Bure S Linak, 125mm, handbrakes, PU, UK
56-302L-B/1-UK	Bure S Linak, 125mm, 1 handbrake, UK
56-302L-B/1-PU-UK	Bure S Linak, 125mm, 1 handbrake, PU, UK
56-302L-US	Bure S, Linak, 125mm, US
56-302L-AUS	Bure S, Linak, 125mm, AUS
56-302L/75	Bure S, Linak, 75mm
56-302L/75-PU	Bure S Linak, 75mm, PU
56-302L/75-B/0	Bure S Linak, 75mm, handbrakes
56-302L/75-B/0-PU	Bure S Linak, 75mm, handbrakes, PU
56-302L/75-B/1	Bure S Linak, 75mm, 1 handbrake
56-302L/75-B/1-PU	Bure S Linak, 75mm, 1 handbrake, PU
56-302L/75-UK	Bure S Linak, 75mm, UK
56-302L/75-PU-UK	Bure S Linak, 75mm, PU, UK
56-302L/75-B/0-UK	Bure S Linak, 75mm, handbrakes, UK
56-302L/75-B/0-PU-UK	Bure S Linak, 75mm, handbrakes, PU, UK
56-302L/75-B/1-UK	Bure S Linak, 75mm, 1 handbrake, UK
56-302L/75-B/1-PU-UK	Bure S Linak, 75mm, 1 handbrake, PU, UK
56-312L	Bure STD Linak, 125 mm
56-312L-PU	Bure STD Linak, 125mm, PU
56-312L-B/0	Bure STD Linak, 125mm, handbrakes
56-312L-B/0-PU	Bure STD Linak, 125mm, handbrakes, PU
56-312L-B/1	Bure STD Linak, 125mm, 1 handbrake

Article number	Product name
56-312L-B/1-PU	Bure STD Linak, 125mm, 1 handbrake, PU
56-312L-UK	Bure STD Linak, 125mm, UK
56-312L-PU-UK	Bure STD Linak, 125mm, PU, UK
56-312L-B/0-UK	Bure STD Linak, 125mm, handbrakes, UK
56-312L-B/0-PU-UK	Bure STD Linak, 125mm, handbrakes, PU, UK
56-312L-B/1-UK	Bure STD Linak, 125mm, 1 handbrake, UK
56-312L-B/1-PU-UK	Bure STD Linak, 125mm, 1 handbrake, PU, UK
56-312L-CAN-KD	Bure STD Linak, 125mm, Knockdown CAN
56-312L-KD	Bure STD Linak, 125mm, Knockdown
56-312L-KD-AUS	Bure STD Linak, 125mm, Knockdown Aus
56-312L-US	Bure STD Linak, 125mm, US
56-312L-AUS	Bure STD Linak, 125mm, AUS
56-312L/75	Bure STD Linak, 75mm
56-312L/75-PU	Bure STD Linak, 75mm, PU
56-312L/75-B/0	Bure STD Linak, 75mm, handbrakes
56-312L/75-B/0-PU	Bure STD Linak, 75mm, handbrakes, PU
56-312L/75-B/1	Bure STD Linak, 75mm, 1 handbrake
56-312L/75-B/1-PU	Bure STD Linak, 75mm, 1 handbrake, PU
56-312L/75-UK	Bure STD Linak, 75mm, UK
56-312L/75-PU-UK	Bure STD Linak, 75mm, PU, UK
56-312L/75-B/0-UK	Bure STD Linak, 75mm, handbrakes, UK
56-312L/75-B/0-PU-UK	Bure STD Linak, 75mm, handbrakes, PU, UK
56-312L/75-B/1-UK	Bure STD Linak, 75mm, 1 handbrake, UK
56-312L/75-B/1-PU-UK	Bure STD Linak, 75mm, 1 handbrake, PU, UK
56-312L/75-AUS	Bure STD Linak 75mm, AUS
56-321L/75	Bure Low Linak, 75mm
56-321L/75-PU	Bure Low Linak, 75mm, PU
56-321L/75-B/0	Bure Low Linak, 75mm, handbrakes
56-321L/75-B/0-PU	Bure Low Linak, 75mm, handbrakes, PU
56-321L/75-B/1	Bure Low Linak, 75mm, 1 handbrake
56-321L/75-B/1-PU	Bure Low Linak, 75mm, 1 handbrake, PU
56-321L/75-UK	Bure Low Linak, 75mm, UK
56-321L/75-PU-UK	Bure Low Linak, 75mm, PU, UK
56-321L/75-B/0-UK	Bure Low Linak, 75mm, handbrakes, UK
56-321L/75-B/0-PU-UK	Bure Low Linak, 75mm, handbrakes, PU, UK
56-321L/75-B/1-UK	Bure Low Linak, 75mm, 1 handbrake, UK
56-321L/75-B/1-PU-UK	Bure Low Linak, 75mm, 1 handbrake, PU, UK

Article number	Product name
56-321L-S/75	Bure Low S Linak, 75mm
56-321L-S/75-PU	Bure Low S Linak, 75mm, PU
56-321L-S/75-B/0	Bure Low S Linak, 75mm, handbrakes
56-321L-S/75-B/0-PU	Bure Low S Linak, 75mm, handbrakes, PU
56-321L-S/75-B/1	Bure Low S Linak, 75mm, 1 handbrake
56-321L-S/75-B/1-PU	Bure Low S Linak, 75mm, 1 handbrake, PU
56-321L-S/75-UK	Bure Low S Linak, 75mm, UK
56-321L-S/75-PU-UK	Bure Low S Linak, 75mm, PU, UK
56-321L-S/75-B/0-UK	Bure Low S Linak, 75mm, handbrakes, UK
56-321L-S/75-B/0-PU-UK	Bure Low S Linak, 75mm, handbrakes, PU, UK
56-321L-S/75-B/1-UK	Bure Low S Linak, 75mm, 1 handbrake, UK
56-321L-S/75-B/1-PU-UK	Bure Low S Linak, 75mm, 1 handbrake, PU, UK

Main configuration is following: 125 mm wheels, no handbrakes, cushions made of synthetic leather, regular hand control, no special features. Following variations are possible:

- 125/100/75 mm – indicates diameter of wheels
- Handbrakes – indicates that the model has a handbrake on both handles
- 1 handbrake – indicates that the model has only 1 handbrake on one of the handles
- UK, AUS, US – models designed specifically for these countries
- Linak- means that the model includes electronics from Linak for adjust of the height
- Gas – means that the model includes hydraulic system for adjustment of the height
- Manual – means that the height of the model is adjusted manually
- Plus – means that the model includes a Plus hand control
- Low – means that the model has lower middle frame for shorter patients
- S – means that the model is more narrow than the standard model
- PU – means that the model includes cushions made of Polyurethane
- DIR/RS – Means that the model has rotation stop

Simple solutions for great results

SystemRoMedic® is the name of Direct Healthcare Group's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic® is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic® is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic® offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading at/from our website;
www.directhealthcaregroup.com

For questions about the product and its use

Please contact your local Direct Healthcare Group and SystemRoMedic® representative. A complete list of all our partners with their contact details can be found on our website; www.directhealthcaregroup.com.

Moving Health Forward

CE Type BF Applied Parts

Medical Device Class I. The product complies with the requirements of the Medical Devices Regulation 2017/745.

Direct Healthcare Group Sverige AB
Torshamnsgatan 35,
SE-164 40 Kista, Sweden
Tel: +46 (0)8-557 62 200
info@directhealthcaregroup.com
www.directhealthcaregroup.com