ClassicSling

System**RoMedic**

Instructions for use – English

Bruksanvisning – Svenska

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REF

Polyester Net:

41222003

41222004

41222006

41222007

41222008

Commode Opening:

41221404

41221406

41221407

SWL: 300 kg/ 660 lbs

Table of contents

| | |
|------------------|----|
| English | 3 |
| Svenska | 5 |
| Norsk | 7 |
| Dansk | 9 |
| Deutsch | 11 |
| Nederlands | 13 |
| Français | 15 |
| Italiano | 17 |
| Polski | 19 |
| Český | 21 |
| Español | 23 |
| | |
| Symbols | 25 |

ClassicSling enables a comfortable and secure seated posture and gives the patient a feeling of security. It is designed to be left in the chair after the transfer is completed, which means that the patient sits on the sling. ClassicSling is available in two models, with or without commode opening. It provides support all the way from the hollow of the knees up to the shoulders. ClassicSling is designed for use in most common lifting situations such as transfers from bed to wheelchair, etc.

Functional inspection

Visual and mechanical inspection

Check the condition and function of the sling regularly. Always inspect the product after laundering. Check to ensure that seams and material are free from damage. Check to ensure that fabric is not worn or faded. Apply load to the device and check to ensure that clasps, handles, etc. withstand heavy load. If there are signs of wear, the product must be discarded.

Always read the manual

Always read the manuals for all assistive devices used during a transfer.

Keep the manual where it is accessible to users of the product.

Do not leave the patient unattended during a transfer situation.

To prevent discomfort and the risk of the patient sliding out of the sling, trial fit the sling carefully. First, apply the sling's back loops to the slingbar, and then apply the leg support loops. Use a slingbar that has been tested.


Note that there may be a risk of the patient tipping forward when a sling with standard (not separate) leg supports is used.

The lowest allowable safe working load always determines the safe working load of the assembled system.

Always check the safe working loads for the lift and accessories before use.

Contact your dealer if you have any questions.

ClassicSling can only be applied when the patient is fully recumbent

1. Turn the patient towards you. Place the sling so that the mid-section follows the spine and the lower edge is about 5 to 10 cm above the hollow of the patient's knees. Place the edge of the sling as far in under the patient as possible.
2. Turn the patient back over, ensuring that he or she is centred over the sling.
3. Preferably, raise the back support before lifting. Start raising the lift, when the straps become taut, stop and check that everything is ok before proceeding with the lift.
4. The sling is left in the chair.

ClassicSling can only be removed when the patient is lying on the bed/stretcher or floor

Follow the same steps as with application, but in the reverse order.

Lifetime information & Periodic inspection

The product has an expected lifetime depending on how the product is used according to the manual instructions. The service life of the product varies depending on usage frequency, washing, drying, etc. Please be sure to read the product label and follow the care instructions for this product on www.directhealthcaregroup.com

The periodic inspection of the non-rigid body support unit should be performed at the timestamps stated by the manufacturer, but at least every 6 months to be sure. More frequently inspections may be required if a non-rigid body support unit is used or cleaned more frequently than advised by the manufacturer. The periodic inspection instructions for this product are found on www.directhealthcaregroup.com

Material:

Polyester net: An airy material that can be left under the patient, if necessary. Withstands high laundry temperatures.

Care of the product:

Read the product label.

Do not use rinsing agent.
To prolong product life,
avoid tumble-drying.

Size information

Slings can often be used for both 2 and 4-point suspension. It is essential to ensure that the width of the slingbar is appropriate for the sling size and the patient's status and capacity.

Use patient-specific slings to ensure safety and hygiene.

| Size | | XS/Orange | S/Red | M/Yellow | L/Green | XL/Purple |
|------------------------|------|-----------|--------|----------|---------|-----------|
| Rec. user weight | lbs | 44-88 | 66-132 | 121-187 | 176-265 | 243-309 |
| | kg | 20-40 | 30-60 | 55-85 | 80-120 | 110-140 |
| Sling overall height | inch | 31,5 | 35,5 | 37,5 | 38,5 | 39,5 |
| | cm | 80 | 90 | 95 | 98 | 100 |
| Sling seat measurement | inch | 12 | 15 | 18 | 25 | 31,5 |
| | cm | 30 | 38 | 46 | 63 | 80 |

Contact your local distributor if you have any questions about the product and its use.

See www.directhealthcaregroup.com for a complete list of distributors. Always make sure that you have the right version of the manual. The most recent editions of manuals are available for downloading from our website, www.directhealthcaregroup.com.

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Symbols

| | | | |
|--|---|---|-------------------|
|  | This product complies with the requirements of the Medical Device Regulation 2017/745 |  | Visual Inspection |
|  | Medical Device |  | Read the manual |
|  | Caution |  | Product Code |
|  | Manufacturer information |  | Batch Code |

Simple solutions for great results

SystemRoMedic® is the name of Direct Healthcare Group's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic® is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic® is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic® offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading at/from our website;
www.directhealthcaregroup.com

For questions about the product and its use

Please contact your local Direct Healthcare Group and SystemRoMedic® representative. A complete list of all our partners with their contact details can be found on our website; www.directhealthcaregroup.com.

Moving Health Forward

Direct Healthcare Group Sverige AB
Torshamnsgatan 35,
SE-164 40 Kista, Sweden

Tel: +46 (0)8-557 62 200
info@directhealthcaregroup.com
www.directhealthcaregroup.com

Medical Device Class I. The product complies with the requirements of the Medical Devices Regulation 2017/745.