# ClassicHighBackSling

SystemRoMedic™

Instructions for use – English

Bruksanvisning – Svenska

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**REF**

41522003

41522004

41522006


41522007

41522008

**SWL: 300kg/660lbs**

## Table of contents

<b>English</b> .....	<b>3</b>
<b>Svenska</b> .....	<b>5</b>
<b>Norsk</b> .....	<b>7</b>
<b>Dansk</b> .....	<b>9</b>
<b>Deutsch</b> .....	<b>11</b>
<b>Nederlands</b> .....	<b>13</b>
<b>Français</b> .....	<b>15</b>
<b>Italiano</b> .....	<b>17</b>

**SWL: 300kg/660lbs**

ClassicHighBackSling enables a comfortable and secure seated posture and gives the patient a feeling of security. It is designed to be left in the chair after the transfer is completed, which means that the patient sits on the sling. It provides support all the way from the hollow of the knees up over the head. ClassicHighBackSling is designed for use in most common lifting situations such as transfers from bed to wheelchair, etc. SystemRoMedic includes a wide range of functional, comfortable, high-quality slings that can be adapted for different types of lifting and for patients with different needs. The choice of model and material depends on the patient's needs and the transfer situation.

## Functional inspection

### Visual and mechanical inspection

Check the condition and function of the sling regularly. Always inspect the product after laundering. Check to ensure that seams and material are free from damage. Check to ensure that fabric is not worn or faded. Apply load to the device and check to ensure that clasps, handles, etc. withstand heavy load. If there are signs of wear, the product must be discarded.

## Always read the manual

Always read the manuals for all assistive devices used during a transfer.

Keep the manual where it is accessible to users of the product.

Do not leave the patient unattended during a transfer situation.

To prevent discomfort and the risk of the patient sliding out of the sling, trial fit the sling carefully. First, apply the sling's back loops to the slingbar, and then apply the leg support loops. Use a slingbar that has been tested.


Note that there may be a risk of the patient tipping forward when a sling with standard (not separate) leg supports is used.

The lowest allowable safe working load always determines the safe working load of the assembled system.

Always check the safe working loads for the lift and accessories before use.

Contact your dealer if you have any questions.

## ClassicHighBackSling can only be applied when the patient is fully recumbent

1. Turn the patient towards you. Place the sling so that the mid-section follows the spine and the lower edge is about 5 to 10 cm above the hollow of the patient's knees. Place the edge of the sling as far in under the patient as possible.
2. Turn the patient back over, ensuring that he or she is centred over the sling.
3. Preferably, raise the back support before lifting. Start raising the lift, when the straps become taut, stop and check that everything is ok before proceeding with the lift.
4. The sling is left in the chair.

## The ClassicHighBackSling can only be removed when the patient is lying on the bed/stretcher or floor:

Follow the same steps as with application, but in the reverse order.

### Material:

Polyester net: An airy material that can be left under the patient, if necessary. Withstands high laundry temperatures.

### Care of the product:

Read the product label.

Do not use rinsing agent.

To prolong product life, avoid tumble-drying.

### Size information

Slings can often be used for both 2 and 4-point suspension. It is essential to ensure that the width of the slingbar is appropriate for the sling size and the patient's status and capacity.

Use patient-specific slings to ensure safety and hygiene.


Size		XS/Orange	S/Red	M/Yellow	L/Green	XL/Purple
Rec. user weight	lbs	44-88	66-132	121-187	176-265	243-309
	kg	20-40	30-60	55-85	80-120	110-140
Sling overall height	inch	37,5	38,5	42	43	45
	cm	95	98	106	110	115
Sling seat measurement	inch	12	16	19	25,5	31,5
	cm	30	40	48	65	80

Contact your local distributor if you have any questions about the product and its use.

See [www.directhealthcaregroup.com](http://www.directhealthcaregroup.com) for a complete list of distributors.

Always make sure that you have the right version of the manual.

The most recent editions of manuals are available for downloading from our website, [www.directhealthcaregroup.com](http://www.directhealthcaregroup.com).

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## Simple solutions for great results

SystemRoMedic™ is the name of Direct Healthcare Group's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic™ is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

## Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic™ is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic™ offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

## Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading at/from our website;  
[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)

## For questions about the product and its use

Please contact your local Direct Healthcare Group and SystemRoMedic™ representative. A complete list of all our partners with their contact details can be found on our website; [www.directhealthcaregroup.com](http://www.directhealthcaregroup.com).

# Direct Healthcare Group

Advancing Movement & Health®

Direct Healthcare Group Sverige AB  
Torshamnsgatan 35,  
SE-164 40 Kista, Sweden

Tel: +46 (0)8-557 62 200  
[info@directhealthcaregroup.com](mailto:info@directhealthcaregroup.com)  
[www.directhealthcaregroup.com](http://www.directhealthcaregroup.com)