# MiniLift160 Classic

 handicare

CE

User manual – English

Bruksanvisning – Svenska

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REF

400641436

400641437


**SWL: 160 kg/430 lbs**

## Table of contents

<b>English</b> .....	<b>4</b>
<b>Svenska</b> .....	<b>20</b>
<b>Norsk</b> .....	<b>36</b>
<b>Dansk</b> .....	<b>52</b>
<b>Suomi</b> .....	<b>68</b>
<b>Deutsch</b> .....	<b>84</b>
<b>Nederlands</b> .....	<b>100</b>
<b>Français</b> .....	<b>116</b>
<b>Italiano</b> .....	<b>132</b>
<b>Español</b> .....	<b>148</b>
<b>Dimensions MiniLift160 Classic</b> .....	<b>164</b>
<b>Dimensions MiniLift160 Classic, with low legs</b> .....	<b>166</b>

MiniLift160 Classic is a mobile sit-to-stand lift which has been developed to, as gently as possible, assist the user when rising from a sitting to a standing position. When MiniLift160 Classic is combined with the appropriate lifting accessories, the user gets support under the feet, for the front of the lower legs and behind the back, which provides for a safe and secure

sit-to-stand procedure. The construction safely moves the user forwards and upwards in a natural movement pattern and, at the same time, the leg muscles and balance is exercised.

## Safety information

### Visual inspection

Inspect lift functions regularly. Check to ensure that material is free from damage.

### Before use

Make certain the lift is properly assembled.

Check sling bar connection and safety latch function.

Check lifting function and base-width adjustment.

Check to ensure that the actuator is correctly installed.

## Always read the user manual

Always read the user manuals for all assistive devices used during a transfer.

Keep the user manual where it is accessible to users of the product.

Always make sure that you have the right version of the user manual.

The most recent editions of user manuals are available for downloading from our website, [www.handicare.com](http://www.handicare.com).

Under no circumstances may the lift be used by persons who have not received instruction in the operation of the lift.

It is strongly prohibited to modify the original product.




# Assembly of MiniLift160 Classic

## Check to ensure that all components are included


Lift motor, control box, battery pack, base with motor, footplate, hand control and cord, manual and charging cable. Screws and screw covers should also be included (see below).

## Assembly


Place the upper section on the lower section and lower it carefully into place (1). Press it fully into place and tighten all four screws. Finally, place the plastic covers over the screws (2).

Mount the footplate on the crossbar and then lower it (picture 3).

Mount the leg support by inserting the crossbars on the back into the fixtures on the lift (4), and then tighten.

Connect cables: Hand control cord in HS, cable for lift unit on mast in outlet M1 and cable for base-width motor in outlet M2 (5).


## Final inspection

Check to ensure that no parts have been left in the packaging.


Inspect the lift for signs of wear and damage.

Check all four castor wheels and castor wheel locks.

Check all connections and fixtures including screws and bolts.


Check the emergency stop function by activating the emergency stop, and then pressing either the up or down button. If nothing happens when the up or down buttons are pressed, the emergency stop is functioning properly.

Grasp the hand control, press the up button and run the lift arm all the way up. Then, press the down button and run the lift arm all the way down.

Test base-width adjustment function. Press the button for base-width adjustment to widen the base fully, and then press the other button to narrow the base again.

Test lift function by lifting a person (not a user) using an approved sling. At the same time, check the emergency lowering function with someone on the lift. See section on Emergency lowering.

If the lift is functioning correctly, connect the charger and check to ensure that the charging lamp on the control box lights up.

### NOTE!

Before the lift is used for the first time:

- it must be charged for at least 4 hours. See section on charging batteries.
- the control box service counter needs to be reset. To reset the service counters press both lift buttons on the hand control at the same time for 5 seconds. An audio signal will indicate that the timer has been reset.

Keep the user manual where it is accessible to users of the product.

## Using the product

### Important Information

- The lift must be assembled according to the assembly instructions provided with the lift.
- The lift may only be used indoors and on a level floor.
- Lifting accessories must be properly fitted and tested in relation to the user's needs and functional ability.
- Do not leave the user unattended during a transfer situation.
- Under no circumstances may max. load be exceeded. See section on Safe working load.
- Ensure that the user's feet do not get caught between the footplate and the floor. Ensure that the user's feet do not get caught between the footplate and the base when adjusting the base width.
- MiniLift160 Classic is not intended for transporting users. It may only be used for short transfers.
- Never move the lift by pulling on the actuator!
- The lift must not be used or immersed in or under water including showering.
- The lift should not be left / stored in a high humidity environment.
- The lift must not be charged in a wet room.
- To ensure optimal function, the lift must be inspected regularly. See section on Maintenance.
- Warranty applies only if repairs or alterations are done by an authorized technician.
- Using lifting accessories other than those approved by Handicare can entail a risk.
- The lift arm must be loaded throughout the entire lowering procedure.

**To transfer with MiniLift160 Classic, the user must be able to:**

- support weight while standing
- understand instructions

## Safe working load

Different products on the same lift system (lift unit, sling bar, sling, scales and other lifting accessories) may have different allowable safe working loads. The lowest allowable safe working load always determines the safe working load of the assembled system. Always check the safe working loads for the lift and accessories before use. Contact your dealer if you have any questions.

## The control box and battery have the following features

1. **Emergency stop**
2. **Charging lamps**
3. **Display showing battery charge level**
4. **Connection for charger cable**
5. **Connection for base-width adjustment motor**
6. **Connection for lift arm actuator**
7. **Connection for hand control**
8. **Electrical emergency lowering**
9. **Handle for lifting battery**

### **Charging batteries**

A tone sounding when using the lift indicates that the battery need recharging.

Charge the lift after use to ensure that the battery is always fully charged.

Lock the castor wheels when charging the battery.

1. Connect the charging cable to a power outlet and in the connection for charger cable.

The symbol for charging is shown on the display.

2. Check to ensure that the lamps on the control box light up. The green LED lamp indicates that the charger is receiving power and the yellow LED lamp indicates that the battery is charging.
3. Charging stops automatically when the battery is fully charged.

Wall-mounted charger.


1. Remove the battery pack from the lift and place it in the wall-mounted charger.
2. Check to ensure that the LED lamp on the front of the charger lights up.

#### **NOTE!**

Before the lift is used for the first time, it must be charged for at least 4 hours.


For maximum battery life, charge batteries regularly. We recommend daily charging when the lift is used daily.

The emergency stop must be deactivated during charging.

#### **Battery information on the display**

The battery discharging will be shown in four stages:

**Battery state 1:** The battery is ok, no need for charging (100 - 50 %).

**Battery state 2:** Battery needs charging. (50 - 25 %)

**Battery state 3:** Battery needs charging. (Less than 25 %) A tone sounds when a button is pressed in this battery state.

**Battery state 4:** The battery needs charging (17V or lower). At this stage some of the functionality of the lift is lost. At this battery stage it is only possible to drive the lifting arm down. Furthermore an audio signal will sound when a control button is activated. The symbol will switch between the two pictures for 10 seconds.

The battery symbol is shown when the box is active until power down (2 minutes after use).

It is not possible to use other battery types than BAJ1/BAJ2.

The battery level is measured via voltage level. This means that it is possible to experience e.g. that the battery switches from state 1 to state 2 and back to state 1.


### Service information read-out

Basic service information can be read out on the display. To get the service information on the display please press the lifting arm up button for half a second. The information will be shown for ½ minute or until other buttons are activated.

↑	12034
A•S	1257000
1g/8h	7
Service	90/360

- Total cycles done by the actuator
- Total work done by the actuator (ampere usage times seconds in use)
- Total number of overloads
- Days since last service/Days between services

### Changing the battery

### Hand control

#### Raising/lowering the lift arm

Symbol indicate direction of travel.

Motion stops as soon as the button is released.

If the lift arm encounter an obstacle while lowering, the lift will immediately stop the motion. In order to continue, the lift arm will need to be raised slightly by using the hand control before further lowering.

#### Widening/narrowing the base

Markings on the buttons indicate function.

Motion stops as soon as the buttons are released.

STOP

## Emergency stop


### To activate the emergency stop

Press the red emergency stop button on the control box.

### Resetting

Turn the button in the direction of the arrows until the button pops out.

To prevent battery discharge, we recommend that the emergency-stop button is pressed in when the lift is not in use.

## Emergency lowering


### Manual emergency lowering

For manual emergency lowering, lift the red slider upwards.


### Electrical emergency lowering

For electrical emergency lowering, use the down button on the control box. Use a narrow object such as a pen.

Manual emergency lowering

Electrical emergency lowering

## Instructions for use

These instructions must be followed for correct use of MiniLift160 Classic.

### Hand control functions

The hand control has four functions: up and down (lift arm), out and in (base width).

Two functions cannot be used at the same time.

**NOTE:** The hand control functions do not operate if the emergency stop button is activated



**NOTE:** MiniLift160 Classic may only be used for short transfers. It is not intended for transporting users.

**Caution** Never move the lift by pulling on the actuator.


**Caution** Ensure that the user's feet do not get caught between the footplate and the floor. Ensure that the user's feet do not get caught between the footplate and the base when adjusting the base width.

## IFU for MiniLift


1. Before using the MiniLift, always plan and assess the risk before a lift and transfer:
  - a. Check that the environment is suitable and plan for the task
  - b. Check the MiniLift and ThoraxSling function
  - c. Check the MiniLift battery
  - d. Check that the MiniLift has been serviced
  - e. Visually inspect the MiniLift and ThoraxSling for signs of damage or defects
  - f. Ensure sizing of the ThoraxSling is correct in accordance with IFU
2. Requirements for use:
  - a. The client must have some standing capability to use the MiniLift – this is dependent on the individual risk assessment.
3. Attach the ThoraxSling to the client according to the ThoraxSling instructions - see IFU for ThoraxSling.
4. Buckle the safety belt. Consider not buckling the safety belt if the client uses a stoma bag, PEG, has wounds or similar. Raising a client with the MiniLift will work well even without buckling the safety belt.
5. If needed, lift the client's feet from the wheelchair's foot support and turn these away.
6. Spread the legs of the MiniLift as wide as possible and position the MiniLift in front of the client, facing the client.
7. Lock the wheels of the MiniLift and of the wheelchair/bed if applicable.
8. Ensure that the client's feet can reach the foot rest and position the client's feet wide apart on the foot rest of the MiniLift, to ensure maximum balance.
9. Adjust the knee support height and depth position so that the client's shins are supported.
10. Lift the locking handle on the standing arm and set the standing arm so that the patient can reach it comfortably with bent arms.
11. Attach suitable loops of the sling on the hooks on the standing arm intended for the loops of the slings.
12. Pull on the slings loops to make sure they are attached properly.
13. Lift the locking handle on the standing arm and set the standing arm so that the loops of the sling are under tension and are not slacking.
14. Place the client's hand on the handles of the standing arm, or let the client push off from the armrests of the wheelchair. For clients with reduced muscle function in the arms, place the arm/s in ThoraxSling armrest.

15. Consider different patterns or ways of raising the client:
  - a. Short setting of the standing arm results in a more forward leaning position that requires more muscle strength in the legs. Suitable for a shorter and more active client.
  - b. Long setting of the standing arm is more suited for a weaker patient who needs training. Suitable for a taller client.

16. Instruct the client about the lifting.
17. Caregiver positioning: The caregiver should stand on the side of the lift and press the up button on the hand control. Watch the foot placement of the client and ensure that the ThoraxSling does not slip.
18. Raise the client to a position where they can stand upright. It is important to note that some clients might not want to stand upright for fear of falling forward.
19. If needed have the client adjust their hand positions to the frame of the MiniLift.

20. Unlock the wheels of the MiniLift.
21. When transferring the patient, close the leg-spreading of the MiniLift and move the MiniLift so that the patient faces forward.
22. Perform desired activity: e.g. transfer, standing training, toileting.
23. Once client is in the desired standing position, move the MiniLift to the intended transfer position.
24. Lock the wheels of the MiniLift and the wheelchair, if the client is to be lowered there.
25. Position the client's hands on the standing arm's handles.
26. The caregiver's position should be standing to the side of the MiniLift. Use the down-button on the hand control to lower the client.
27. If possible, encourage the client to lean forward during the lowering so that it resembles standing up without aid.
28. Lower the standing arm and ensure that the client ends up positioned well in the sitting position (e.g. in a wheelchair).
29. Let the client remove their hands from the MiniLift, extend the standing arm and detach the ThoraxSling completely from the MiniLift.
30. Unlock the MiniLift wheels and carefully remove the client's feet from the MiniLift's footrest.
31. Move the MiniLift away from the client.
32. Remove the ThoraxSling from the client.

## Trouble-shooting

### If the lift or base-width adjustment cannot be activated, check the following:

- That the emergency stop function is activated.
- That all cables are properly and securely connected. Pull out the contact and plug it in again firmly.
- That battery charging is not in progress.
- That the battery is charged.

If the lift is not working properly, contact your dealer.

### If the lift makes unusual noises:

- Try to determine the source of the sound. Take the lift out of operation and contact your dealer.

## Accessories

ThoraxSling, polyester, XS-XXL, article no.: 45500003-009

ThoraxSling, wipeable, S-XXL, article no.: 45500004-009

ThoraxSling, with seat support, polyester, S-XL, article no.: 45600004-008

ThoraxSling, with seat support, disposable, non-woven, S-XL, article no.: 45690004-008

CalfStrap, article no.: 70200033

### Hand control

Hand control HB33-6, 6 buttons, with service, battery status and overload indicators, article no.: 70200089

## Maintenance

The lift must undergo thorough inspection at least once per year. Inspection must be performed by authorized personnel and in accordance with Handicare's service manual.

Repairs and maintenance may only be done by authorized personnel using original spare parts.


Used batteries are to be left at the nearest recycling station. Used batteries can also be returned to Handicare or a Handicare dealer for recycling.

### Cleaning/disinfection

Clean the lift with warm water or rubbing alcohol and ensure that the castor wheels are free of dirt and hair. Do not use cleaning agents containing phenol or chlorine, as this may damage the material.

### Storage and transportation

If the lift is not to be used for some time or e.g., during transport, we recommend that the emergency stop button be pressed in. The lift should be transported and stored in  $-10^{\circ}\text{C}$  to  $+50^{\circ}\text{C}$  and in normal humidity, 30% -75%. The air pressure should be between 700 and 1060 hPa. Let the lift reach room temperature before the batteries are charged or the lift is used.

### Service agreements

Handicare offers the possibility of service agreements for maintenance and regular testing of your mobile lift. Contact your local Handicare representative.

### Normal use environment

This equipment is intended for indoor use. It is not recommended to be used outside normal indoor temperature and humidity ranges,  $10^{\circ}\text{C}$  to  $40^{\circ}\text{C}$  and 30% to 75%.

## Symbols

	White/blue Read user manual
	May not be discarded in domestic waste
	The product complies with the requirements of the Medical Devices Directive 93/42/EEC
	Type B, according to the degree of protection against electric shock
	The device is intended for indoor use
	Class II equipment
	Important Be observant
	Red/ black Do not puch or pull the lift by the actuator

## Technical Information

<b>Lifting speed</b>	35 mm/s without load.
<b>Batteries</b>	Two 12V, 2.9 Ah valve-regulated, sealed, lead accumulator (gel-type batteries)
<b>Charger</b>	Max. 400mA
<b>Motor (mast)</b>	DC 24 V, 10.5 Ah. Operating time: 10% maximum continuous operation for 2 minutes Push: 7500N.
<b>Motor (base)</b>	DC 24 V, 12.3 Ah. Operating time: 10% maximum continuous operation for 2 minutes Push: 2000N.
<b>Sound level</b>	55.8 dB(A)
<b>Material</b>	Steel
<b>Emergency lowering</b>	Manual and electrical
<b>Castors</b>	Front 3,1", 80 mm, back 3,1", 80 mm
<b>Weight</b>	39 kg/104,5 lbs
<b>Weight loose parts (Foot plate)</b>	5 kg / 11 lbs
<b>Footplate, measurement</b>	390 x 345 mm / 15.3 x 13.5
<b>Number of lifts per charge</b>	20-25
<b>IP class</b>	IP X4
<b>Expected lifetime</b>	10 years
<b>Operating forces buttons on hand control</b>	4 N
<b>Max. load</b>	160 kg/ 430 lbs

Medical Device Class I. The product complies with the requirements of the Medical Devices Directive 93/42/EEC.

## Electromagnetic Compatibility (EMC)

Product has been tested for compliance with current regulatory standards regarding its capacity to block EMI (electromagnetic interference) from external sources. Some procedures can help reduce electromagnetic interferences:

- Use only Handicare cables and spare parts to avoid increased emissions or decreased immunity which can compromise the correct functioning of the equipment.
- Ensure that other devices in patient-monitoring and/or life-support areas comply with accepted emissions standards.

**WARNING**

**Wireless communications equipment such as wireless computer network devices, mobile phones, cordless telephones and their base stations walkie-talkies, etc. can affect this equipment and should be kept at least 1.5m away from the equipment.**

Intended Environment: Home Healthcare Environment and Professional Healthcare facility environment

Exceptions: HF Surgical Equipment and RF Shielded room of an ME SYSTEM for magnetic resonance imaging.

**WARNING**


**Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.**

The MiniLift160Classic is a device for safely raising persons from a sitting to a standing position using a moving lift arm. In case the unit stops operating due to possible EM disturbance, the unit should be switched off. Its function should be verified and the equipment should be used in another room or environment. In case of other equipment being disturbed, its function should be verified before using again.

Guidance and manufacturer's declaration - electromagnetic emissions		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	This equipment uses RF energy only for its internal functions. Therefore its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic emissions			
Immunity test	IEC60601-1-2 test level	Compliance	Electromagnetic environment - guidance
Surge IEC 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage reduction in % of UN: 100% for 1 period, 100% for 0.5 periods, 30% for 25 periods, 100% for 250 periods	Voltage reduction in % of UN: 100% for 1 period, 100% for 0.5 periods, 30% for 25 periods, 100% for 250 periods	Mains power quality should be that of a typical commercial or hospital environment. If the user of the MiniLift160Classic requires continued operation during power mains interruptions, it is recommended that the MiniLift160Classic be powered from an uninterruptible power supply or a battery.

Guidance and manufacturer's declaration - electromagnetic immunity			
Immunity test	IEC60601-1-2 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 2\text{kV}$ , $\pm 4\text{kV}$ , $\pm 8\text{kV}$ , $\pm 15\text{kV}$ air  $\pm 8\text{kV}$ contact	$\pm 2\text{kV}$ , $\pm 4\text{kV}$ , $\pm 8\text{kV}$ , $\pm 15\text{kV}$ air  $\pm 8\text{kV}$ contact	Floor should be wood, concrete or ceramic tile. If the floors are covered with synthetic material, the relative humidity level should be at least 30%.
Conducted disturbances induced by RF fields IEC 61000-4-6	3V in 0.15 MHz to 80MHz  6V in ISM and amateur radio bands between 0.15 and 80 MHz  80% AM at 1 kHz	3V in 0.15 MHz to 80MHz  6V in ISM and amateur radio bands between 0.15 and 80 MHz  80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the product, including cables, than 1.0m, if the transmitter's output power rating exceeds 1W. <sup>(a)</sup> Field strength from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. <sup>(b)</sup> Interference may occur in the vicinity of equipment marked with this symbol:  
Radiated RF electromagnetic field IEC 61000-4-3	Home Healthcare Environment 10 V/m 80 MHz to 2.7 GHz 80% AM and 1 kHz	Home Healthcare Environment 10 V/m 80 MHz to 2.7 GHz 80% AM and 1 kHz	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	385 MHz - 27 V/m 450 MHz - 28 V/m 710, 745, 780 MHz - 9 V/m 810, 870, 930 MHz - 28 V/m 1720, 1845, 1970, 2450 MHz - 28 V/m 5240, 5500, 5785 MHz - 9 V/m	385 MHz - 27 V/m 450 MHz - 28 V/m 710, 745, 780 MHz - 9 V/m 810, 870, 930 MHz - 28 V/m 1720, 1845, 1970, 2450 MHz - 28 V/m 5240, 5500, 5785 MHz - 9 V/m	
Electrical fast transient/burst IEC 61000-4-4	Signal and control lines: $\pm 1\text{kV}$ AC: $\pm 2\text{kV}$ 100 kHz repetition frequency	Signal and control lines: $\pm 1\text{kV}$ AC: $\pm 2\text{kV}$ 100 kHz repetition frequency	
Power frequency Magnetic field IEC 61000-4-8	30 A/m  50 Hz and 60 Hz	30 A/m  50 Hz and 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<p>(a) Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MiniLift160 Classic is used exceeds the applicable RF compliance level above, the MiniLift160 Classic should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary.</p> <p>(b) Over the frequency range 150 kHz to 80 MHz, field strength should be less than 1V/m.</p>			

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

# Dimensions MiniLift160 Classic

All measurement are in mm/inch  
Tolerance +/- 5 mm/ 0.2"


All measurement are in mm/inch  
Tolerance +/- 5 mm/ 0.2"


# Dimensions MiniLift160 Classic, with low legs

All measurement are in mm/inch  
Tolerance +/- 5 mm/ 0.2"

All measurement are in mm/inch  
Tolerance +/- 5 mm/ 0.2"

## Simple solutions for great results

Handicare offers a unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. We have a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

## Improved work environment, improved quality of care and cost savings

The Handicare patient transfer product range is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, Handicare offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

## Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading from our website; [www.handicare.com](http://www.handicare.com).

## For questions about the product and its use

Please contact your local Handicare representative. A complete list of all our partners with their contact details can be found on our website; [www.handicare.com](http://www.handicare.com).

Handicare offers solutions and support to increase the independence of disabled or elderly people as well as to improve the convenience of those who are caring for them.

*The Handicare Group is one of the leading healthcare companies in Europe with its own manufacturing organizations and sales companies in Norway, Sweden, Denmark, Germany, the Netherlands, Great Britain, France, China, Canada and the USA. Handicare's products are also distributed by partners in more than 40 countries worldwide. Our wide range of high-quality products includes a complete easy transfer system and other patient handling aids, stairlifts, car adaptations and bathing and toileting products.*

Handicare AB  
Torshamnsgatan 35  
SE-164 40 Kista, SWEDEN  
Tel: +46 (0)8-557 62 200  
Fax: +46 (0)8-557 62 299  
E-mail: [info@handicare.se](mailto:info@handicare.se)  
[www.handicare.com](http://www.handicare.com)