




Dyna-Form[®] Air Suresse

This is a quick reference guide for the **Dyna-Form Air Suresse**
(refer to user manual for additional information).


Product Code MAT/SURESSE/200/85/25



On/Off

Power Button

The Power Button turns the system power on and off by pressing the Power Button for at least two seconds.



CPR Function

Rapid deflation of the Mattress may be required for emergency treatment or to decommission the system. Firmly pull the Rapid Release / CPR Tag from the side of the Mattress to rapidly deflate the entire system. To re-inflate the system after the Rapid Release / CPR Tag has been removed replace as such, ensuring all sealing connectors are firmly attached and restart the Control Unit. Wait for the Mattress system to gain optimal pressure. Perform a Bottoming Out test after inflating the mattress following rapid deflation.



Soft



Medium



Firm



Dynamic



Static

Soft Mode, Medium Mode, Firm Mode

Press buttons to increase or decrease pressure setting. The **Soft**, **Medium** & **Firm** settings allow comfort to the user, without clinical compromise.

The GREEN LEDs illuminate to indicate which of the three settings is operational.

Press **Dynamic Mode** for alternative cells cyclically inflating and deflating. Static Mode will automatically revert to Alternation Mode after one hour for patient safety.

Upon power up, the system automatically reverts back to the dynamic mode operating at the previous pressure setting for patient safety. Static Mode will automatically revert to Alternation Mode after one hour for patient safety.

Press to facilitate **Static Mode** for clinical procedure / patient transfer purposes. After 20 minutes, the system automatically reverts back to the previous pressure setting for patient safety.

Press Static Mode for all cells to be fully inflated with no dynamic alternation. Static Mode will automatically revert to Alternation Mode after one hour for patient safety.

NB: Please ensure that all securing straps on the base of the mattress are secured onto the NON MOVING PARTS of the bed frame.



Audible Warning

Silences the audible warning (on / off). Audible warning will resume after 20 minutes if cause of failure not resolved.



Warning Function

The **Red Warning LED** (A,B or C) flashes, and an audible alert sounds, to indicate the control unit or mattress pressure has failed. The LED will remain illuminated until appropriate pressure is restored. The audible warning can be silenced by pressing the Warning Mute button. **The system has three different warning signals, identified by illumination sequences.**



A B C

The system cannot reach the set pressure within 8 minutes.

The pressure is too **HIGH**.



A B C

The system cannot reach the set pressure within 8 minutes.

The pressure is too **LOW**.



A B C

Power unit has **NO** power feed (Mains Failure).



Lock

Lockable Mode

Press for at least two seconds to lock the Control Unit settings – a beep sounds and the amber LED illuminates to indicate system is locked. When locked, only the Warning Mute and Lock / Unlock buttons remain operational. Press again for at least two seconds to unlock (beep sounds and amber LED turns off).

The Control Unit will automatically unlock in the event of a power failure.

Symptoms Problems / Cause Points to check

| | | |
|---------------|---|---|
| Low Pressure | The mattress is set to a mode that is too SOFT. | Change the mode button (from soft to a firmer pressure setting) as required. If the mattress is still too soft after a short period of 10 minutes, then please call an engineer. Ensure CPR valve is closed. |
| | The CPR has been knocked /moved. | Check all tubing is not kinked within the mattress. |
| | There may be a leak in the system. | Ensure that the tubing within the mattress are fully connected. |
| High Pressure | The mattress is excessively firm on a constant basis. | Set mattress to a softer setting as clinically required. |
| | | Evaluate that the mattress to be of a "less firm" state after a short period of 10 minutes. If this is not achieved, then please follow the task as below before calling an engineer for assistance. NB: Check all tubing is not kinked within the mattress. |



**Direct Healthcare
Services**
Delivering the Promise

Direct Healthcare Services Ltd
Unit 8-10 Withey Court, Western Industrial Estate
Caerphilly, CF83 1BF

Sales Office +44 (0) 845 459 9831 **F** +44 (0) 845 459 9832 **E** info@directhealthcareservices.co.uk
www.directhealthcareservices.co.uk

Although Direct Healthcare Services Ltd will endeavour to provide the exact specification as printed, we reserve the right to change specification without prior notice. All equipment & information within this brochure, should be used in conjunction with appropriate clinical judgement & nursing procedures.
All products manufactured to standards:

