



Amfibi XL

Amfibi XL is a hygiene chair designed with bariatric users in mind. With a maximum user weight of 250kg, 67cm internal width and electric heigh adjustment, Amfibi XL hygiene chair provides stability and comfort for the user while requiring minimal effort for the care staff.



Amfibi XL has a unique pillar design for quick and easy seat height adjustment. The seat can be raised by up to 81cm from its lowest position and thus requires minimised effort from care personnel.

Technical Specifications

Part N°	Product	Width, External	Max. User Weight	Seat Height	Width, Internal	Seat Depth
56-204	Amfibi XL	79cm	250kg	55-81cm	67cm	49cm