

Nightingale Bariatric Floor Bed

A bed with Endless possibilities



TRUE FLOORLINE BED WITH NO TOWERS

ONE BED - ENDLESS POSSIBILITIES

AUSTRALIAN DESIGNED & MADE



*'The only bed
of it's type in
the world'*




Nightingale Bariatric Floor Bed

*...the true floorline bed
with no towers...*






Unique dual Safe Working Load of 200kg or 300kg - simply swap handsets



Low height of 100 and high height of 725mm



No towers. Up to 400mm shorter than other floor level beds



Universal solution for Low, High & Bariatric care



Chair position at the push of a button



10 year frame warranty & 5 year Linak warranty

We are proud to introduce the first ever universal bed that suits low care, high care and bariatric care. Plus, this innovative three-in-one bed has a low height of 100mm and NO TOWERS!

From care of frail elderly or dementia patients, through to acute and bariatric settings, the patent pending Australian made Nightingale Bariatric Floor Bed has you covered. One bed, endless possibilities!

Nightingale Bariatric Floor Bed

A bed with endless possibilities

ENDLESS BED POSITIONS




High position



Medium position



Low position



Backrest function



Backrest & Kneebreak



Backrest & Kneebreak



Trendelenburg



Reverse Trendelenburg



Chair position

SPECIFICATIONS

Overall dimensions	2080mm x 905mm	
Sleep surface dimensions	2000mm x 890mm	
SWL	200kg/300kg	
NET Bed Weight	100kg	
Low Position	100mm	
High Position	725mm	
Underbed Clearance	Clearview	
Angles	Back	62°
	Knee	35°
	Trendelenburg	15°
	Reverse	15°

Warranty	Frame	10 years
	Electric	5 years
Optional Features	Underbed Light	
	Battery Backup	
	Dual Braking	
	End Board Options	
Product Codes	Single	EN9000
	King Single	EN9000KS
	Double	EN9000D
	Queen	EN9000Q
	King	EN9000K