

ECLIPSE 600

bariatric manual wheelchair



- For people of significant size
- Delivers style and performance
- Forward wheelbase for excellent maneuverability
- Transportable
- Narrow overall width



WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

WWW.PDGMOBILITY.COM

ECLIPSE 600

THE RIGHT MANUAL BARIATRIC WHEELCHAIR FOR YOUR NEXT CLINICAL APPLICATION



Nine standard frame colors to choose from and custom colors also available. Contact us for details.



Transportable.



Forward wheelbase.



Adjustable seat depth and back angle.



Narrow overall width.



Solar cushions specially designed for bariatric applications

Seat width, adjustable 20" – 48"
Seat depth, adjustable 16" – 24"
Seat height, adjustable 14" – 20"
Back angle, adjustable 85° – 100°
Back height 25" and 30"
Overall width Seat width plus 7"

Overall length
Weight capacity
Front rigging
Chair weight

31.5" - 34.5" (without front rigging)
 600 lbs or 1000 lbs upgrade
 Optional centre mount
 35 lbs (without backposts, armrests, front rigging, rear wheels, anti-tippers)

Warranty

PDG's exclusive TGIF warranty (Thank Goodness It Fits) Lifetime warranty on the frame.
 PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product. Canadian Patent # 2,198,643 US Patent # 5,957,474

USA
 PDG Product Design Group
 PMB # 173 - 1124 Fir Avenue
 Blaine, WA
 98230, USA

CANADA
 PDG Product Design Group
 Unit 103 - 318 East Kent Avenue South
 Vancouver, BC
 V5X 4N6

T: +1 604 323 9220
 F: +1 604 323 9097
 Toll Free (in N. America): 888-858-4422
 Toll Free Fax (in N. America): 877-323-9093
www.pdgmobility.com