

ELITE PRODUCT OPTIONS



FRONT & REAT
LOADING PTTY



SUBSTANTIAL
FRAMEWORK THAT
CATERS FOR SEAT
REVERSAL

X399 BARIATRIC TOILETING
/SHOWERING CHAIR



SEAT POSITIONED WITH
GAP AT FRONT



SEAT POSITIONED WITH
GAP AT BACK



REVERSE VIEW SHOWING 'PRAM' HANDLES

ELITE KEY FEATURES

- Designed to wheel over standard toilet or Clos-O-Mat or be used for assisted showering
- Manufactured from high gauge epoxy powder coated steel, this frame is robust and wide enough for bariatric comfort, but is dimensionally sized to fit through a standard door
- Seat and back are polyurethane so are soft feel for patient comfort and easy to clean
- Back frame work is reinforced and designed to fit flush with toilet cistern
- Adjustable 'pram handles' ensure carer best positioning when pushing patient
- Seat is reversible; the gap can facilitate personal hygiene with access to front or back
- Seat has PU 'plug' fitting for use as a seat only
- Pan carrier fitted to underside of the seat provides lever 'clip in' holder for the oval potty with handle
- Oval pan with handle can be inserted at front or rear of frame and so allows frame to 'park' against a wall and still be used
- Swing away arms also detach for side transfer
- Footplates clip-on and swing away
- Substantial 125mm (5") braked castors all round

ELITE CODE	X399
SEAT HEIGHT	610mm 24"
FLOOR TO UNDERSIDE	540mm 21 1/4"
WIDTH BETWEEN ARMS	550mm 21 3/4"
SEAT WIDTH	640mm 25 1/4"
SEAT DEPTH	460mm 18"
SEAT THICKNESS	46mm 1 1/4"
SEAT TO TOP OF BACK	390mm 15 1/4"
HEIGHT (HANDLES EXTENDED)	1110mm 43 1/2"
OVERALL WIDTH (INC. CASTERS)	750mm 29 1/2"
OVERALL DEPTH (INC. FOOTPLATES)	980mm 38 1/2"
ORIFICE WIDTH	270mm 10 1/2"
ORIFICE DEPTH	410mm 16"
FOOTPRINT	750mm x 980mm 29 1/2" x 38 1/2"
MAXIMUM USER WEIGHT	250kg 40st

ELITE PRODUCT OPTIONS



SWING AWAY &
DETACHABLE ARMRESTS



PRAM HANDLES
ADJUST IN HEIGHT



SWING AWAY FOOT PLATE

ELITE SPARES / ACCESSORIES



71019F/H
OVAL POTTY WITH
HANDLE



72PU500S
EASY CLEAN PU SEAT &
72PU500SP SLOT IN PLUG



72PUSB
EASY CLEAN
PU BACK