




# Bure Ortho/Ortho S/Space

Ortho is especially designed for orthopaedics departments and situations that demand extra legroom and room for manoeuvre for the user. Ortho has a powerful electric motor for powered height adjustment. Walker heights are thus adjustable with the patient in situ and the walker can be used as a stand-up aid. Ortho is also available in a small version – Ortho S (only available in the gas version) – with a max frame width of only 65 cm. In addition to Ortho's advantages, Bure Space also has extra room to get up close to wheelchairs or similar. The base frame has been made wider and the internal width extended to 67 cm. Bure Space is available with electrically operated or hydraulic height adjustment.



## Technical specifications

Part No.	Product	Width	Max user weight	Length	Height	Armrest/Inner dimension
56-317T	Bure Ortho, 125mm castors	78cm	150kg	80cm	90-130cm	25-50cm
56-318	Bure Ortho Gas, 125mm castors	78cm	150kg	80cm	90-130cm	25-50cm
56-318-S	Bure Ortho S Gas, 125mm castors	65cm	150kg	79cm	90-130cm	30-45cm
56-304	Bure Space EI, 125mm castors	78cm	150kg	79cm	90-130cm	25-50cm
56-307	Bure Space Gas, 125mm castors	78cm	150kg	79cm	90-130cm	25-50cm

All walkers may be ordered with 75, 100 or 125mm castors – specify when ordering. (Unless otherwise specified, the heights shown below refer to products with 125 mm castors). Examples of possible parts combinations are shown below.