



Bure Ortho/Ortho S/Space

Ortho is especially designed for orthopaedics departments and situations that demand extra legroom and room for manoeuvre for the user. Ortho has a powerful electric motor for powered height adjustment. Walker heights are thus adjustable with the patient in situ and the walker can be used as a stand-up aid. Ortho is also available in a small version – Ortho S (only available in the gas version) – with a max frame width of only 65 cm. In addition to Ortho's advantages, Bure Space also has extra room to get up close to wheelchairs or similar. The base frame has been made wider and the internal width extended to 67 cm. Bure Space is available with electrically operated or hydraulic height adjustment.

Ergonomic handles
Stepless adjustment in every direction.

Height-adjustable handles
See accessories.

Low minimum height
Only 85 cm (when 75 mm castors are used).

Easily reachable hand controller
The hand controller is located such that both of the carer's and user's hands can easily reach it.

Stable
The walker has anatomical pads that provide relief for shoulders, arms and the neck and shift the focus to larger muscle groups such as leg and chest muscles instead. This improves stability and makes the walker easier to steer.

Waterproofing
All Bure Double electrical components have IP 65 international protection ratings and can therefore be used in wet spaces.

Flexible stepless adjustment
This makes it easy to pre-set the walker to the desired height.

Power Rise
Bure Ortho provides excellent assistance for helping users stand up.

Brakes and directional locks
Located so that the user can reach them from the walking position.

Unrestricted front legroom
Users with contractures and patients with full leg casts can keep their legs straight and still bring the walker close to them.

Open design
Provides the same stability but with greater room for manoeuvres.

More space
Bure Space has a base frame with internal width of 67cm - more space to get close with a wheelchair.

Technical specifications

| Part No. | Product | Width | Max user weight | Length | Height | Armrest/Inner dimension |
|----------|---------------------------------|-------|-----------------|--------|----------|-------------------------|
| 56-317T | Bure Ortho, 125mm castors | 78cm | 150kg | 80cm | 90-130cm | 25-50cm |
| 56-318 | Bure Ortho Gas, 125mm castors | 78cm | 150kg | 80cm | 90-130cm | 25-50cm |
| 56-318-S | Bure Ortho S Gas, 125mm castors | 65cm | 150kg | 79cm | 90-130cm | 30-45cm |
| 56-304 | Bure Space EI, 125mm castors | 78cm | 150kg | 79cm | 90-130cm | 25-50cm |
| 56-307 | Bure Space Gas, 125mm castors | 78cm | 150kg | 79cm | 90-130cm | 25-50cm |

All walkers may be ordered with 75, 100 or 125mm castors – specify when ordering. (Unless otherwise specified, the heights shown below refer to products with 125 mm castors). Examples of possible parts combinations are shown below.