



Bure XL

Bure XL was designed for the bigger, heavier patient and therefore provides extra space and strength. Because extra powerful aids are required in the care and nursing of heavier patients, Bure XL is always supplied with electric drives.

Ergonomic handles

Ergonomically designed handles that are steplessly adjustable in every direction for optimum individual comfort. This is extremely important for stroke patients who can thus reach their “usual” hand position.

Height-adjustable handles
See accessories.

Sustainable

StandTall Walkers are designed to last. For example, they have extra reinforcement bands around the arm rest cushions to prevent them from splitting from impacts to the side.

Stable

The walker has anatomical cushions that provide relief for shoulders, arms and the neck, and shift the focus to larger muscle groups such as leg and chest muscles instead. This improves stability and makes the walker easier to steer.

Easily reachable hand controller

The hand controller is located such that both of the carer's or user's hands can easily reach it.

Flexible stepless adjustment

Electrically driven height adjustment makes it easy to pre-set the walker to the desired height. Once the patient is standing in position the height can be steplessly adjusted for comfort and the proper support for an optimal walking configuration.

Power Rise

Bure XL is a good aid for helping patients and users stand and reduces strain on employees. Some patients can even get up on their own with the aid of the walker.

Designed for weight

The Bure XL is designed to provide support for heavier patients up to 240 kg. This StandTall Walker has also been widened to provide good stability and plenty of space. The inner dimension is 44–62 cm. The extra leg room makes it especially easy for patients who need extra walking space due to e.g. leg spasms.

Open design

Provides the same stability but with greater room for manoeuvres.

Waterproofing

All Bure XL electrical components have IP 65 international protection ratings and can therefore be used in wet spaces.

Hospitals and nursing homes

Because of its size, Bure XL is best suited to hospitals and nursing homes.

Technical specifications

Part No.	Product	Width	Max user weight	Length	Height	Armrest/Internal dimensions
56-315	Bure XL EL	88cm	240kg	78cm	91-131cm	44-62cm

All walkers may be ordered with 75, 100 or 125mm castors – specify when ordering. (Unless otherwise specified, the heights shown below refer to products with 125 mm castors). Examples of possible parts combinations are shown below.