



Moving Health Forward

## Rental Solutions

### Bariatric & Standard Equipment Guide

Acute & Community Beds & Mattresses

Pressure Care Cushions

Specialist Seating

Walkers & Walking Frames

Transfer/Sit-to-stand Aids

Mobile Hoists

Slings & Sheets

Commodes/Shower Chairs

Gantry Systems & Hoisting Accessories

Bariatric Wheelchairs

Scale Systems

**Emergency sales & rental (24/7/365)**

**0800 879 9289**



Moving Health Forward

## Vision

To advance the health of individuals with affected movement, utilising innovative and proven ways to enhance the certainty of care needed today and the care we plan to provide for tomorrow.

## Mission

DHG constantly support our healthcare partners, and the people whose healthcare needs they serve, by providing solutions that improve the vital connections between movement and health. All leading to the best recovery outcomes and healthy independence in daily living.

## Our Story

DHG is a global medical device company that supports efforts to move health forward by delivering innovative and leading solutions for individuals whose lives are impacted by reduced movement.

We strive to constantly support our healthcare partners, and the people whose healthcare needs they serve, working to keep improving the vital connections between movement and health, enabling accelerated recovery and healthy independence in daily living.

At DHG, we are passionate to work in close collaboration with healthcare professionals to enhance outcomes, lower the cost of healthcare, deliver efficiencies, and expand access to our solutions.

## Our Brands



### Pressure Ulcer Prevention

DHG specialises in the manufacturing of pressure ulcer prevention solutions with a comprehensive portfolio of innovative, award-winning products including mattresses, cushions and overlays.

With significant experience in the provision of clinically proven, value for money solutions, customers benefit from a combination of innovative pressure care technologies, designed in partnership with leading clinical and healthcare establishments.

In 2016, Direct Healthcare Services' Intelligent Pressure Care Management won the prestigious Queen's Award for Enterprise in Innovation following the introduction of the Dyna-Form® Mercury Advance; a keystone hybrid mattress.



### Safe Moving & Handling

DHG's SystemRoMedic range has extensive heritage in the safe movement and handling of users across the spectrum of care environments. Providing solutions for lifting systems and manual transfer aids, the SystemRoMedic portfolio enables safer working ergonomics for the carer, whilst providing clinical effectiveness for the user. The clinical focus of the SystemRoMedic looks to release more time back to care by providing a highly-experienced, trained moving & handling team who can assess for the most appropriate moving, handling, specialist seating, in-bed or out-of-bed solution to provide effective care.



### Specialist Therapies

DHG's Specialist Therapies portfolio offers advanced technologies for use with avoidable harms.

Our Negative Pressure Wound Therapy (NPWT) offers intuitive, simple-to-use, cost effective solutions for the accelerated healing of chronic wounds.

With the knowledge that Venous Thromboembolism (VTE) is a global, life-threatening issue, DHG offer Intermittent Pneumatic Compression (IPC) solutions for use in the prevention of Pulmonary Embolism.



### Bathroom Safety

Renowned for more than 80 years, DHG's Linido bathroom safety solutions provide increased independence, safety and comfort in bathrooms, showers and toilet areas across private and public environments, as well as in hospitals. DHG offers a wide range of safe, practical, durable and ergonomic aids with modern design; creating accessible, safe and aesthetic bathrooms.

Our bathroom safety products are extensively tested for their functionality, including safety testing as well as ergonomic use and comfort, to ensure to comfort of safety every day. The complete range provide users in any bathroom or toilet situation with as much autonomy as possible.



### Rental & Service Solutions

Through DHG's Nightingale brand, we have over 20 years' experience in ad-hoc or contracted rental and servicing solutions that enable a focus on patient care across any operational budget. We offer a 24/7/365 rental service providing a range of pressure area care solutions, beds for acute and community environments, safe moving & handling equipment, specialist seating, and bariatric care equipment; ensuring that your patient receives the right care at the right time. To facilitate easier choice, we have designed a number of packages focusing on early discharge, single-handed care, rehabilitation, or capacity pressures.

With a network of nationwide, ISO 13485 accredited rental and service centres, we work closely with healthcare professionals to ensure our solutions are reliable, safe and effective. Our solutions are delivered by a highly experienced Technical Service team, who will provide on-site support and training during the delivery process.

Our ProServe® Service Solutions provide ISO 13485 compliant servicing options for the wide range of products your hospital or equipment service have on site. We offer a range of technical service packages designed to suit the needs of our customers in today's ever changing healthcare environment.

As a global supplier of equipment that is designed to facilitate better health outcomes, accelerate patient recovery and provide easier daily living for the elderly, disabled or those providing care and support to them, DHG is acutely aware that for equipment to be used to its full potential, those using it must have an understanding of everything a piece of equipment can do and how it should be used.



The Ethos Education Programme underpins this theory by providing practical tips and techniques supported by theory-based learning for a wide range of people including those responsible for specifying, installing or using DHG equipment. The Ethos Education Programme provides access to 8 CPD accredited workshops that will encourage you to review current practices to ensure client and carer safety when involved with moving and handling tasks or assessing for a specialist seat. The Ethos Education Programme is designed to deliver a comprehensive overview of how equipment can be used to assist in single and multi-carer environments whilst giving you the opportunity to consider your individual circumstances and following on from the workshop, how you could reduce your risk of injury whilst simultaneously improving the moving and handling experience and dignity for those within your care.







Plus-size Care Management Webcast Series:

## “Applying Innovative Systems In The Management of Plus-size Individuals”

Presented by Anita Rush

### Learning Aim:

Delegates will have the opportunity to develop their knowledge base to enable effective management and implement safer systems of work for Plus-Size individuals.

### During this Webcast, delegates will gain an understanding of:

- The global incidence of obesity and obesity related mortality
- Different terminologies (meadtric and bariatric), their influence on plus-size individuals body dynamics and their impact on normal body movement and subsequent moving and handling tasks
- The challenges of managing a plus-size individual and implementation of suitable systems to enable benefits both for the individual and the care providers
- The key elements involved in the risk and mobility assessments required for a plus-size individual
- How using an ergonomic system approach plus size interventions can enable best practice and harm free care
- How pivotal the discharge planning process and a cohesive communication channel across all stakeholders involved is to enable a safe discharge
- Complex cases that can be solved and safe systems of work implemented (case study)
- The need for a strategic approach to the management of plus size individuals

### The Webcasts are split into three parts:

• Part I – Evidence Based Practice

• Part II – Risk Management of Plus-size Individuals

• Part III – Complexities of Plus-size Management

### Anita Rush

M.Sc Health Ergonomics, Dip Health Care Studies, SRN Clinical Lead and Clinical Nurse Specialist Equipment

Anita is a highly experienced Clinical Lead Nurse Specialist for equipment provision across Berkshire. The role generates many complex patient assessments: including tissue viability, postural management, seating and manual handling for the population of Berkshire and beyond, especially in the field of Plus Size Management.

Anita has developed a competency-based training and education programme that enables a proactive approach to risk management. Sourcing solutions for patients with complex needs is difficult, but with good outcomes and you see an improvement in the patient's quality of life then it is all worthwhile.

Anita has lectured for the Disabled Living Foundation, published a number of articles, has presented at many national and international conferences on the management Plus Size Individuals. Anita also delivers bespoke Plus Size Training sessions for Health and Social Services within the UK including Ireland.

Anita has co-authored two publications: Chapter 12, in 'The Guide to the Handling of People, a Systems Approach' – 6th edition reprinted May 2013, published by Backcare the charity for healthier backs and in collaboration with the National Back Exchange. The second publication 'Moving and Handling of Plus Size People – an illustrated guide' was published by the National Back Exchange in 2013. Anita is a member of the National Back Exchange.

## BARIATRIC BED RANGE

### Aurum®+ Bariatric Bed (500kg SWL)

#### Acute Bariatric Adjustable Bed Frame With Integrated Weighing Scales



The Aurum® + Bariatric width-adjustable, profiling bed frame is designed for size, strength and practicality. Fully compliant to the latest hospital bed standards, the bed features width and length adjustment.

Featuring the most precise integrated scales on the market, ensuring accurate weight monitoring throughout the patients care journey.

Complete with stowable safety rails to ensure caregivers can conduct manual handling activities without obstruction.

**Order Code:** T100PBASP

#### Specifications:

- 78st/500kg SWL (69st/440kg max patient weight)
- Width adjustable from 900mm to 1200mm
- Extendable length up to 2200mm

#### Profiling Positions



### Olympia Bariatric Bed (500kg SWL)

#### Acute Bariatric Width-Adjustable Bed Frame

The Olympia Bariatric width-adjustable, profiling bed frame is fully compliant to the latest hospital bed standards. It can be used as a standard bed at 900mm and is expandable to 1100mm, 1200mm and 1300mm, with an extendable length to 2200mm. The Olympia Bariatric has a user friendly, quiet operation to assist with the provision of dignified care.

**Order Code:** 0111

#### Specifications:

- 78st/500kg SWL (68st/435kg max patient weight)
- Width adjustable from 900mm to 1300mm (with 1100 and 1200 options)
- Extendable length up to 2200mm

#### Profiling Positions



## Versatech 1100 Ultra Low Bariatric Bed (500kg SWL) Expandable Bariatric Bed with Integrated Weigh Scale

A groundbreaking innovation in bariatric care, the VersaTech 1100 has an ultra low height of just 8 inches from the ground to enhance patient safety and support ingress and egress.

Easily extendable in width from 36" to 48" the bed caters to a wide range of patients and can be used as both a standard and bariatric bed. The adjustable frame also benefits from length extension at the touch of a button, extending the bed from 80-88"

Featuring integrated weighing scales with a bed exit alarm for continuous patient monitoring and ultra quiet motors for optimum patient comfort.

With the option of VersaDrive technology, the bed can be easily maneuvered by caregivers, reducing strain and improving workflow efficiency.

**Order Code:** 0117

### Specifications:

- 78st/500kg SWL
- Width adjustable from 890mm to 970mm to 1040mm to 1190mm (35" to 38" to 41" to 47")
- Extendable length from 2003mm to 2130mm to 2230mm (80" to 84" to 88")

### Profiling Positions



## XXL Xpress Bariatric Care Bed (270kg SWL)

The XXL Xpress medical bariatric bed combines both the form of a modern design with the functionality needed to assist carers in giving superior levels of care to patients. It helps to reduce the risks associated with the manual handling of this patient group. The extra width provides a patient surface that is 4 feet in width or more (1200mm or 1400mm) in a bed that can be easily transported in and out of community environments.

**Order Code:** 0110

Plus mattress option if required.

### Specifications:

- 43 stone/270kg SWL
- 2000mm (L) extendable to 2200mm (L) x 1200mm (W)

### Profiling Positions



## STANDARD BED RANGE

Compatible with foam, hybrid/ dynamic, turn assist and low air loss mattress options.



### Olympia Hi/Lo Acute Care Bed (230kg SWL)

The Lago Hi/Lo Acute Care Bed is a perfect symbiosis of safety and comfort. A low bed position of  $\pm 25\text{cm}$  reduces the risk of falling to an absolute minimum, whereas the ergonomic working height of  $\pm 80\text{cm}$  guarantees efficient care.

**Order Code:** 0112

**Specifications:**

- 36 stone/230kg SWL
- Height adjustment: 250mm–800mm

**Profiling Positions**



Compatible with foam, hybrid/ dynamic, turn assist and low air loss mattress options.



### Aerys Hi/Lo Community Bed (230kg SWL)

The Aerys Hi-Lo bariatric bed combines safety and comfort. A low bed position of  $\pm 24.5\text{ cm}$  reduces the risk of falling to an absolute minimum, whereas the ergonomic working height of  $\pm 80\text{ cm}$  guarantees efficient care. The Aerys bariatric bed has been especially designed to meet the needs and the requirements of Acute hospitals, Long-term care units and Nursing homes/Care homes and being easily transferable through standard doorways it can be installed where other beds could prove difficult.

**Order Code:** 0113

**Specifications:**

- 36 stone/230kg SWL
- Height adjustment: 250mm–800mm
- Mattress platform dimensions: 2000mm x 900mm

**Profiling Positions**





## Deprimo Low Bed (250kg SWL)

The deprimo® floor level bed has been developed to reduce the risk of vulnerable patients injuring themselves by falling out of bed, increasing patient safety and reducing the costs of in-hospital injury. The fully articulating bed has a 120mm (4.7") ultra-low height and is available with crash mats that are placed either side of the bed to enable a patient to simply slide on and off the bed.

Innovative patented stowable safety rails\* can be easily and safely stored under the mattress platform.

**Order Code:** T100DLB

### Specifications:


- 39 stone/250kg SWL
- Height adjustment: 120mm–820mm (4.7"-32.3")
- Mattress platform width 900mm (35.4")



### Profiling Positions



## BARIATRIC SUPPORT SURFACES RANGE



### Dyna-Form® Maximus (413kg Weight Limit)

Specifically designed to meet the needs of bariatric patients up to 413kg/ 60 stone. Its finest quality Combustion Modified Ether (CME) foam type and construction is designed to support the weight of heavier patients, prevent 'bottoming out' and offer effective pressure reduction for users.

This product provides not only the best pressure re-distribution properties, but great patient comfort & product longevity.


Complete with durable, multi-stretch, vapour permeable cover and toughened polyurethane base.

**Mattress Order Code:** 0409

**Risk Category:** High Risk

**Dimensions:** 1980mm (L) x 1200mm (W) x 150mm (D)

**Intelligent®**  
Pressure Care  
Management



### Dyna-Form® Mercury Advance Bari WA

(380kg Weight Limit)

A width-adjustable mattress replacement system that can deliver both pressure redistribution in static mode and pressure relief in dynamic mode, the Dyna-Form® Mercury Advance Bari WA is designed for use with patients at Very High Risk of developing pressure damage.

Offering high levels of patient comfort, this unique system has the facility to 'step up' to that of a dynamic mattress when clinically required. Similarly the mattress function can be 'stepped down' as the patient's condition improves. Can be used at 900mm, 1100mm, 1200mm and 1300mm widths.

**Mattress Order Code:** 0415

**Pump Order Code:** 1209

**Risk Category:** Very High Risk

**Dimensions:** 1980mm (L) x 900–1100/1200–1300mm (Width Adjustable) x 150mm (D)



## Verto Lateral Turning Mattress

(450kg Weight Limit)

Designed for use with patients suffering from the inability to move or turn themselves and have therefore developed pressure sores on those areas of the body where insufficient blood is supplied.

The lateral rotation feature enables a gentle, gradual, 30-40 degree turning angle of the patient allowing capillaries to fill, preventing blood starvation and the formation of pressure sores. In this way, existing pressure sores can be treated successfully. The mattress is a low air loss system which keeps the patient cool and dry without blowing air directly onto the patient.



**Mattress Order Code:** T100TM

**Risk Category:** Very High Risk

**Dimensions:** 2000mm (L) x 1200mm (W) x 200mm (D)

**Turn Angles:** 1/4, 1/2, 3/4 and full turn (40°)

**Turn Time Options:** 10/20/30/60 minutes

## STANDARD SUPPORT SURFACES RANGE

### Dyna-Form® Mercury Advance

(254kg Weight Limit)



The Dyna-Form® Mercury mattress is a dynamic mattress replacement system combined with the benefits of modern foam technology. It is specifically designed for patients considered to be at 'Very High Risk' of pressure ulcer development. Offering high levels of patient comfort, this unique system has the facility to 'step up' to that of a dynamic mattress when clinically required. Similarly, the mattress function can be 'stepped down' as the patient's condition improves.

These features make it particularly beneficial for use within the patient's home or palliative care environment and help reduce logistic and decontamination costs. The clinical benefits of a single system are equally applicable to those of a modern hospital setting.

**Order Code:** 0420

**Risk Category:** Very High Risk

**Dimensions:** 1980mm (L) x 880mm (W) x 150mm (H)

### Quattro® Plus

Suitable for all patient types, the Quattro® Plus mattress system combines an active 1-in-4 cell cycle, with TISSUEgard, Ortho-Differential Support and a multistretch cover to provide patients with pressure relief to prevent and manage patients at risk of or with existing pressure ulcers.



**Order code:** 0418

**Risk category:** High-Very High Risk

**Dimensions:** 1950mm (L) x 880mm (W) x 180mm (D)

## STANDARD SUPPORT SURFACES RANGE

### Quattro® Acute (250kg Weight Limit)

The Quattro® Acute mattress system uniquely offers DEEP CELL THERAPY™, combining 1-in-4 alternation with continuous low pressure to provide pressure ulcer prevention and treatment in the most dependent of patients.

**Order Code:** 0419

**Risk Category:** High-Very High Risk

**Dimensions:** 1950mm (L) x 900mm (W) x 180mm (H)



### Carital® Optima (300kg Weight Limit)

The Carital® Optima is a dynamic, anti-deformation mattress replacement system. Offering high levels of patient comfort, the mattress utilises continuous low pressure to provide pressure redistribution through high levels of immersion and envelopment.

**Order Code:** 0303

**Risk Category:** Very High Risk

**Dimensions:** 2000/2100/2200/2300mm (L) x 700/750/800/850/900mm (W) x 130mm (H)



## BARIATRIC CUSHION RANGE

**Intelligent®**  
Pressure Care  
Management



### Dyna-Form® Mercury Advance Bari Alternating Air Cushion (318kg Weight Limit)

The Dyna-Form® Mercury Advance Bari Alternating Cushion is a dynamic replacement seating system, combined with the benefits of modern foam technology. It is specifically designed for patients considered to be at Very High Risk of pressure ulcer development.

Offering high levels of patient comfort, this unique system has the facility to 'step up' to that of a dynamic cushion when clinically required. Similarly, the cushion function can be 'stepped down' as the patient's condition improves.

**25" seat width**  
**Order Code: 1101**

**28" seat width**  
**Order Code: 1101**

**30" seat width**  
**Order Code: 1101**

**32" seat width**  
**Order Code: 1101**

**34" seat width**  
**Order Code: 1101**

**30" seat width**  
**Order Code: 1101**

## STANDARD CUSHION RANGE

**Intelligent®**  
Pressure Care  
Management



### Dyna-Form® Mercury Advance (152kg Weight Limit)

The Dyna-Form® Mercury Advance Alternating Cushion is a dynamic replacement seating system, combined with the benefits of modern foam technology. It is specifically designed for patients considered to be at "Very High Risk" of pressure ulcer development.

Offering high levels of patient comfort, this unique system has the facility to "step up" to that of a dynamic cushion when clinically required. Similarly, the cushion function can be 'stepped down' as the patient's condition improves. A higher maximum weight capacity, up to 24 stone /152kg, allows the product to meet the challenges of heavier patients.

**19" seat width**  
**Order Code: 1101**

**22" seat width**  
**Order Code: 1101**

## BARIATRIC SEATING RANGE



### Static Bedside Chair

Complete with drop arms and height adjustment, this bariatric bedside chair has pressure relieving memory foam in the seat and back cushion to provide added comfort for the user.

**28" seat width** (318kg SWL)  
**Order Code:** 0502

**32" seat width** (318kg SWL)  
**Order Code:** 0501



### Fero Obesitas Bedside Mobile Chair

A versatile bedside chair with electric backrest recline operation and drop arms. This chair can also be used as a patient transfer chair as it comes with larger castors, rear handle and a footplate.

**30" seat width** (318kg)  
**Order Code:** 0507



### Pro-Axis Riser Recliner Chairs

Hoist compatible – Option to hire integral Mercury Advance with riser recliner

These heavy duty chairs are ergonomically designed and offer pressure relieving seats, footrests and waterfall backs with full electric rise and recline functionality. All chairs are also compatible with mobile hoists and stand aids. Adjustable seat depth up to 5" with dynamic back rest adjustment to accommodate different body shapes. The new 34" seat width has a SWL of 70st/444kg.

**25" seat width** (318kg SWL)  
**Order Code:** 0504

**31" seat width** (318kg SWL)  
**Order Code:** 0503

**28" seat width** (318kg SWL)  
**Order Code:** 0504

**34" seat width** (444kg SWL)  
**Order Code:** 0503



### Hi-Lo Tilt-In Space Bariatric Day Chair (300kg SWL)

Offering the capability and functionality of a bed with full electric operation and a high dependency castellated foam surface. This chair is designed to be used as the first step from bed to chair for high dependency patients and allows for easier bed to chair lateral transfers. Can be used to assist with sit to stand exercises during early stages of a patients rehabilitation. Large castors and a rear handle also allow for patient transportation.

**Order Code:** 0508

## STANDARD SEATING RANGE

Kirton specialist seating is manufactured to include clinically proven pressure relieving seat modules – Reflexion foam or Dyna-Tek Intelligent Air.



### Pro-Axis Riser Recliner Chairs

Hoist compatible – Option to hire integral Mercury Advance with riser recliner

These heavy duty chairs are ergonomically designed and offer pressure relieving seats, footrests and waterfall backs with full electric rise and recline functionality. All chairs are also compatible with mobile hoists and stand aids. Adjustable seat depth up to 5" with dynamic back rest adjustment to accommodate different body shapes.



**19" seat width** (120/140kg SWL)  
**Order Code:** 0505

**22" seat width** (120/140kg SWL)  
**Order Code:** 0505



### G-2 (180kg SWL)

For the semi- and non-ambulant user with high posture and pressure care needs, the G-2 features electric tilt-in-space action and built-in, helical-motion lateral supports.

**Order Code:** 0708



Manual



Powered



Regulated Motion



Tilt-in-Space



Articulating Legrest



### Florien II (180kg SWL)

The Florien II is a dynamic and highly mobile seating system that enables greater independence by promoting correct posture and facilitating a supportive, comfortable seated position. Featuring manual tilt-in-space, leg rest elevation and removable side arms for easier hoist access and side transfers.

Available with lateral support.

**Order Code:** 0707



Manual



Tilt-in-Space



Articulating Legrest



### Florien Elite (160kg SWL)

Florien Elite is ideal for ambulant, semi-ambulant and non-ambulant users as well as in multi-user environments. The unique upholstered arms float across the seat to achieve appropriate seat width for each user. The arms also can swivel backward to allow a side transfer or easier hoist access.

**Order Code:** 0717



Manual



Tilt-in-Space



Back Angle Recline



Articulating Legrest





## Duo Range (160kg SWL)

The Duo family of multi-positioning seating systems comprises the Duo Major, Duo Minor and Duo Mini. Their dynamic features provide the ultimate in positioning and postural management for those who are at risk of developing an unconventional body shape, unable to maintain a midline position and require additional support.

**Order Code:** 0712



Manual



Powered



Regulated Motion



Tilt-in-Space



Back Angle Recline



Articulating Legrest



Seat Height Adjustment



## Cura Care Chair (180kg SWL)

The Cura Care Chair is a high dependency chair, designed to be a vital part of a patient's care plan. The simple operation of the seat allows easy positioning of the patient and when combined with the pressure management cushioning built into the seat and backrest, patient comfort is cared for. The seat width measures at 18".

**Order Code:** T100CC



Manual



Tilt-in-Space



Back Angle Recline



Articulating Legrest



Seat Height Adjustment

## BARIATRIC WHEELCHAIR RANGE



### Goliath Wheelchair (300kg SWL)

The Goliath heavy-duty wheelchair is built for strength and durability, constructed from high-quality steel tube. The wheelchair is fully folding.

The Goliath has elevating, height adjustable and detachable armrests to help caregivers with patient transfers. It is also equipped with a rigid, upholstered seat cushion which is fully detachable.

#### Specifications

- 28" 30" and 31" Seat Widths available
- 21.5" Seat Depth

**Order Code:** T500GW



### Phoenix Wheelchair (185kg SWL)

The Phoenix heavy duty wheelchair is built for durability and versatility. Constructed from high-strength aluminium, the wheelchair is fully folding.

The Phoenix wheelchair has elevating, height adjustable and detachable arm rests to help care givers with patient transfers. It comes complete with a rigid, upholstered seat cushion with webbing in seat and backrest allowing complete adjustment.

#### Specifications

- 590mm (23") Seat Width

**Order Code:** T300PW

## COMMODES / SHOWER CHAIRS – BARIATRIC RANGE



### Freeway T100 Bariatric Shower Chair

The Freeway T100 Bariatric Shower Chair is suitable for a variety of applications for bariatric clients. The T100 benefits from increased diameter tubing, special heavy duty castors, heavy duty footrests with height adjustable heel supports which are removable and swing away.

Other commodes/showering accessories available.

**28" seat width (350kg)**  
**Order Code:** 0706

**30" seat width (350kg)**  
**Order Code:** 0706

**32" seat width (350kg)**  
**Order Code:** 0706

### Extra Wide Commode (453kg SWL)



Designed to accommodate the widest of patients comfortably, the commode features padded, drop arms to assist with side transfers and four strong braked castors provide a stable seat for the patient.

#### Specifications:

- 453kg SWL
- 36" Seat Width
- 20.5" Seat Height

**Order code:** T100C/70

### Standard Bariatric Commode (337kg SWL)



Our Standard Bariatric Commode can also be used as a shower chair and will fit over a standard sized toilet, making the commode extremely versatile.

Made with an aluminium frame, the commode is high strength yet still lightweight and it has four strong braked castors to ensure stability whilst in use. The swing away arms enable easy side transfers.

#### Specifications:

- 337kg SWL
- 28" Seat Width
- 21" Seat Height

**Order code:** T100C/53

## WALKING FRAMES AND WALKERS – BARIATRIC RANGE



### Walking Frame (318kg SWL)

Our height adjustable Bariatric Walking Frame is lightweight but strong enough to cater for patients weighing up to 318kg (50st). The walking frame has an extra wide frame top to accommodate larger patients and is constructed with moulded hand grips to provide comfort for the user.

**Order Code:** 0902



### Stand Tall XXL Walking Aid (325kg SWL)

The Stand Tall is a multi-functional aid; It can be used for walking, mobilising from seated to standing, it is safe to use, gives independence and has an IV option. The Stand Tall has 2 built-in handgrips on the side which allow the user to stand up in small steps from the bed or chair being in full control of the situation. This minimises the risk of falling and therefore also minimises the risk of injury to caregiver or patient himself. Once upright the user can walk around.

**Order Code:** 0901



### Rota Stand Solo (200kg SWL)

The Rota Stand Solo revolutionises sit to stand and pivot transfers. The patented T bar mechanism in the base allows a sit to stand and pivot transfer to be carried out by a single nurse or carer.

Not only does this device minimise costs incurred with patient care, it reduces the time it takes to carry out the task.

**Order Code:** T100RSS



### Rota Stand Compact (200kg SWL)

Suitable for sit-to-stand and pivot transfers, the Rota Stand Compact combines all the benefits of the Rota Stand Solo with the added value of being lightweight and foldable.

The compact and lightweight design means this standing aid is ideal for community assessments by therapists who often have to lift equipment in and out of their vehicles. Equally the Rota Stand Compact encourages more independence for users who wish to travel to friends and relatives.

**Order Code:** T100RSC

## TRANSFER / SIT-TO-STAND AIDS – BARIATRIC RANGE



### **ReTurn7600** (205kg SWL)

Assists with sit-to-stand and patient rehabilitation. The ReTurn7600 is an easy-to-use and incredibly effective assistive device providing convenient and natural support for sit-to-stand and transfer, but it is also a rehabilitation device. The ReTurn7600 enables safe and active sit-to-stand and transfer to or from bed, wheelchair or toilet. The ReTurn7600 can also be used to facilitate repositioning farther back in a chair or wheelchair.

**Order Code:** 0904

---

## – STANDARD RANGE



### **ReTurn7500i** (150kg SWL)

Enables safe, natural support for a sit-to-stand and transfer. Can also be used to facilitate repositioning further back in a chair or wheelchair or used as a rehabilitation device.

**Order Code:** 0905



### **ReTurn7400** (150kg SWL)

ReTurn7400 has a low rising ladder that is ideal for shorter adults and children. It also works well in special care situations, such as in an x-ray unit.

**Order Code:** 0906



### **QuickMove (170kg SWL)**

For users with impaired balance and reduced strength in their legs. The QuickMove supports and activates the user when rising to a standing position. The user must have some strength in their arms to actively participate in the transfer.

**Order Code:** 0511



### **MiniLift160 Classic & ThoraxSling (160kg SWL)**

A compact, mobile sit-to-stand lift to gently assist the user when rising from a sitting to a standing position. Combined with re-useable ThoraxSling, the user is supported under the feet, at the front of the lower legs and across their back.

**Order Code:** 0512



### **MiniLift200 & ThoraxSling (200kg SWL)**

A mobile sit-to-stand lift to gently assist the user when rising from a sitting to a standing position. Combined with re-useable ThoraxSling, the user is supported under the feet, at the front of the lower legs and across their back.

**Order Code:** 0513

## MOBILE HOISTS – BARIATRIC RANGE



### Eva600EE Mobile Hoist (270kg SWL)

A versatile mobile hoist, complete with a SlingBar and with electrical base widening for users weighing up to 270kg/ 600lbs. Eva600EE is strong and stable, yet lightweight.

**Order Code:** 0515

---



### Liko Viking XL Mobile Hoist (300kg SWL)

Includes re-usable sling

An increasing need for lifting heavy patients in different care situations led to the development of a mobile lift with high lifting capacity, made of aluminium to be as easy to manoeuvre as possible. This hoist comes with patient stability bars allowing for use as a rehabilitation frame with the Lift Pants as an accessory (please see accessories).

**Order Code:** 0601

---

## MOBILE HOISTS – STANDARD RANGE



### Eva450EE Mobile Hoist (205kg SWL)

Includes re-usable sling

A versatile mobile hoist, complete with a SlingBar and with electrical base widening for users weighing up to 205kg/450lbs. Eva450EE is ideal lifting in connection with transfers to and from the floor, chair or bed, as well as for horizontal lifting and gait training.

**Order Code:** 0514



## GANTRY SYSTEMS – BARIATRIC RANGE



### Liko Ultra Twin Gantry Hoist (400kg SWL)

Includes re-usable sling

With a very high lifting capacity this gantry hoist has two lift motors that work together to meet very high demands. The Ultra Twin not only lifts and transfers, it also helps position the patient in sitting, semi-reclining or horizontal positions. Adjustable width span of 2m, 2.5m and 3m available.

**Order Code:** 0602

**Specifications:**

Total Height 2270mm – 2570mm

## HOISTING ACCESSORIES – BARIATRIC RANGE



### Liko Ultra Twist Plate (500kg SWL)

For use with Gantry Hoist only

Allowing for safe patient rehabilitation between the span of the gantry hoist system with 360 degree rotation.

**Order Code:** 1115

## GANTRY SYSTEMS – STANDARD RANGE



### Castor Gantry and RiseAtlas 450 Ceiling Lift (205kg SWL)

Castor gantry is a safe, high-performance and portable freestanding gantry which, in combination with a 3m rail and RiseAtlas450 lift, functions as a stationary lift system. Castor gantry is an ideal solution for occasional lifting needs in a specific location, or where it is not possible to install a permanent rail system on the ceiling.

**Order Code:** 0516

**Specifications:**

Total Height 2000mm – 2400mm

Lifting Range 1920mm – 2270mm

## SLINGS – BARIATRIC RANGE



### **BariVest** (300kg/570kg SWL)

For use with a Gantry or Mobile Hoist

Can be used with a gantry system or mobile lifts with a four-point spreader bar for rising, standing training and gait training. Available in XXXL (recommended user weight: 180-300kg) or XXXXL (recommended user weight: 250-570kg).

**Order Code:** XXXL: 0541 XXXXL: 0542



### **Liko Ultra Lift Pants** (500kg SWL)

For use with a Gantry or Viking XL Mobile Hoist

Lift Pants can be used with an overhead lift or with a mobile lift with sufficient lifting height allowing for safe patient rehabilitation. Available in XXL or XXXL.

**Order Code:** 0604

## SLINGS – STANDARD RANGE



### **HighBackSling** (300kg SWL)

Polyester. Available in sizes: XXS, XS, S, M, L, XL & XXL.

**Order Code:** 0570



### **WalkingVest** (300kg SWL)

For rising and gait training with a stationary or mobile lift. Available in sizes: XXS, XS, S, M, L, XL & XXL.

**Order Code:**

**XXS:** 0530 **XS:** 0531 **S:** 0532 **M:** 0533 **L:** 0534 **XL:** 0535 **XXL:** 0536



## SHEETS

### Safe Handling Sheet

For horizontal lifting and positioning with a two-point spreader bar.

For positioning users in bed, without subjecting the caregiver to unnecessary strain.

Available in sizes L polyester (300kg SWL), XL Polyester (400kg SWL) and XXL polyester reinforced (570kg SWL)

**Order Code:** L: 0521 XL: 0522 XXL: 0523



## EVACUATION SHEET

### S-CAPEPOD (300kg SWL)

The easy to use S-CAPEPOD allows for quick evacuation of immobile bariatric patients up to 300kg (47.3st).

This simple to deploy evacuation device is stored between the frame of the bed and the mattress and when required can be unfolded and fitted around the patient, taking around 60 seconds to get the patient ready for evacuation.\*

Using the hand-hold straps at the head and foot-end allow the mattress and occupant to be easily manoeuvred along the floor and down stairs to a point of safety.

\*training may be required in order to achieve swift evacuation

**Order Code:** T100SCAPE3

## PATIENT WEIGHING



### M200 Chair Scales (300kg SWL)

Class III Compliant

Heavy duty wheeled chair scales with a high capacity wide seat, hinged armrests and foot plates. This chair allows for dignified patient weighs with a simple to use display unit which includes a BMI function.


**Order Code:** 0803

# Bariatric Equipment Selection Guide


Our professional teams are always on hand to help you choose the best solution for your needs. Here are some patient details you might want to consider when selecting bariatric equipment:

## Body Shape


Apple



Pear Abducted



Pear Adducted



Bulbous Gluteal



| Age Group | Child   | Young Adult   | Adult         | Mature Adult   | Elderly     |
|-----------|---------|---------------|---------------|----------------|-------------|
| Height    | to 5ft  | 5ft to 5ft 6" | 5ft 6" to 6ft | 6ft to 6 ft 6" | over 6ft 6" |
| Weight    | to 20st | 20–30st       | 30–40st       | 40–50st        | 50–65st     |

| Mobility Assessment                   | Pressure Ulcer Risk Assessment |
|---------------------------------------|--------------------------------|
| No Limitation (Can stand and walk)    | Low                            |
| Slightly Limited (Walks with support) | Low to Medium                  |
| Limited (Can bear weight)             | Medium                         |
| Very Limited (Chair bound)            | Medium to High                 |
| Completely Immobile (Bed bound)       | High                           |
| Not Conscious                         | High to Very High Risk         |

**Simple Patient Assessments make equipment choice easier and more cost effective**

|                               |                 |
|-------------------------------|-----------------|
| Order Maker: (Name and Title) | Tel:            |
| Hospital/Ward:                | Email:          |
| Patient Name/ID:              | Weight: Height: |
| Body Shape:                   | Mobility:       |
| PAC Needs:                    | Date and Time:  |
| Ward Contact:                 | Tel:            |
| Trust:                        | ORDER NO:       |

**Always order equipment in line with local protocol and clinical assessment**

# Weight Conversion

Metric conversions given to the nearest KG

| KG        | Stone | KG         | Stone | KG         | Stone | KG         | Stone |
|-----------|-------|------------|-------|------------|-------|------------|-------|
| <b>1</b>  | 0.15  | <b>60</b>  | 9.44  | <b>200</b> | 31.49 | <b>340</b> | 53.54 |
| <b>2</b>  | 0.31  | <b>70</b>  | 11.02 | <b>210</b> | 33.06 | <b>350</b> | 55.11 |
| <b>3</b>  | 0.47  | <b>80</b>  | 12.59 | <b>220</b> | 34.64 | <b>360</b> | 56.69 |
| <b>4</b>  | 0.62  | <b>90</b>  | 14.17 | <b>230</b> | 36.21 | <b>370</b> | 58.26 |
| <b>5</b>  | 0.78  | <b>100</b> | 15.74 | <b>240</b> | 37.79 | <b>380</b> | 59.83 |
| <b>6</b>  | 0.94  | <b>110</b> | 17.32 | <b>250</b> | 39.36 | <b>390</b> | 61.41 |
| <b>7</b>  | 1.10  | <b>120</b> | 18.89 | <b>260</b> | 40.94 | <b>400</b> | 62.98 |
| <b>8</b>  | 1.25  | <b>130</b> | 20.47 | <b>270</b> | 42.51 | <b>410</b> | 64.56 |
| <b>9</b>  | 1.41  | <b>140</b> | 22.04 | <b>280</b> | 44.09 | <b>420</b> | 66.13 |
| <b>10</b> | 1.57  | <b>150</b> | 23.62 | <b>290</b> | 45.66 | <b>430</b> | 67.71 |
| <b>20</b> | 3.14  | <b>160</b> | 25.19 | <b>300</b> | 47.24 | <b>440</b> | 69.28 |
| <b>30</b> | 4.72  | <b>170</b> | 26.77 | <b>310</b> | 48.81 | <b>450</b> | 70.86 |
| <b>40</b> | 6.29  | <b>180</b> | 28.34 | <b>320</b> | 50.39 | <b>460</b> | 72.43 |
| <b>50</b> | 7.87  | <b>190</b> | 29.91 | <b>330</b> | 51.96 | <b>470</b> | 74.01 |

Notes:



Moving Health Forward

PRESSURE ULCER PREVENTION

SAFE MOVING & HANDLING

SPECIALIST THERAPIES

BATHROOM SAFETY

SPECIALIST SEATING

RENTAL & SERVICE SOLUTIONS



DHG  
Withey Court, Western Industrial Estate  
Caerphilly, United Kingdom  
CF83 1BF

T: +44 (0) 845 459 9831  
F: +44 (0) 845 459 9832  
E: [info@directhealthcaregroup.com](mailto:info@directhealthcaregroup.com)



Date: JANUARY 2025

**DHG - HEALTHCARE.COM**