




Dyna-Form™ Mattress Overlay

The **Dyna-Form Mattress Overlay** has been created for users considered to be **“At Risk”** of pressure ulcer development. To be used in conjunction with any standard mattress, ‘Dyna-Form’ design quality comprises of a specialist “castellated cut” CME foam surface, complete with a sewn multi-stretch, vapour permeable cover and non-slip base. This product is ideal for use in the homecare environment.

Also available in viscoelastic foam.

Features

Benefits

1a 1b	CME castellated foam (with tear duct cut)	1a smaller castellated cells for heels and 1b larger cells sacally
2	Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)	Reduces shear and friction forces, whilst assisting infection control policies
3	Sewn or welded cover (NB: welded cover option carries a price premium)	Prevents fluid ingress & aids decontamination processes (can be laundered up to 95 °C)
4	Non-slip base (NB: carries a price premium)	Prevents slippage on top of mattress replacement, maintaining patient safety
5	Fully replaceable parts	Ensures great value for money

Technical Specifications

Product Code	MAT0410001
Risk Category	At Risk
Weight Limit	26 Stone / 165kg
Warranty	2 Years Foam / 2 Years Cover
Dimensions	(Length) 198cm x (Width) 88cm x (Height) 6-9cm, or made to measure Available in Double. (Variable sizes available on request)
Foam Density	34–36kg/m ³ Nominal Hardness – 125 Newtons
Fatigue Class	Very Severe
Fire Retardancy	BS7177; 2008 (Crib 5 Composite)
Product Weight	6kg
Colour	Light Blue
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)