




Dyna-Form™ Maximus

The **Dyna-Form Maximus** mattress is specifically designed to meet the needs of bariatric patients up to 65 stone / 413kg. Its finest quality Combustion Modified High Resilience (CME) foam type and construction is designed to support the weight of heavier patients, prevent “bottoming out” and offer effective pressure reduction for users considered to be at “**High Risk**” of pressure ulcer development. This product provides not only the best pressure re-distribution properties, but great patient comfort & product longevity. Supplied at 120 cm wide as standard it can however be made to the dimensions of any bariatric bedframe. It comes complete with a durable, multi-stretch, vapour permeable cover and a toughened polyurethane base.

Available in viscoelastic foam.

Features

- ① Finest quality combustion modified high resilience (CME) foam construction
- ② Multi-stretch, vapour permeable cover fabric (impregnated with an anti-fungal agent)
- ③ Toughened PU / non-slip base

Benefits

- Provides excellent comfort and pressure re-distribution whilst preventing bottoming out & aiding patient transfer
- Reduces shear & friction forces, whilst assisting in the control of infection (can be laundered up to 80 °C)
- Provides a robust protective surface against hydrolysis / strikethrough and prevents slippage

Technical Specifications

Product Code	MAT0910001
Risk Category	High Risk
Weight Limit	65 stone / 413kg
Warranty	4 Years
Dimensions	(Length) 198cm x (Width) 120cm x (Height) 15cm
Foam Density	Blue 38–40kg/m ³ Nominal Hardness – 200 Newtons
Fatigue Class	Very Severe
Fire Retardancy	BS7177 (Crib 5 Composite)
Mattress Weight	17kg
Colour	Blue
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)