



# EASY SoftBack

EASY™ is a well-designed hygiene stool with a great many advanced functions. The hygiene chair is designed with telescoping pillars to create as much space and flexibility as possible for care personnel. The chair is suitable for use in shower and toilet environments in care institutions and at home.

EASY™ can be raised, lowered and tilted forwards and backwards even with the user in the chair. It's functions make the EASY™ Unequalled aid for personnel while also providing great comfort for the user (for example, its ability to tilt really helps take the strain). EASY™ is available with our new, soft Comfort seat with an opening at the rear for a more natural, hygienic working method. EASY™ SoftBack is available in 2 versions -with a conventional backstrap/necksupport (SoftBack Comfort) or with Ribbons/PUR (SoftBack Flexi).

## EASY SOFTBACK

### Adjustable back-strap

The soft, comfortable back-strap relieves pressure on bony areas. It is easily adjustable for best support and seat depth. Machine washable.

### Individual fit

EASY has a neck support that is adjustable in height and from side to side for comfortable support and individual fit.



EASY SoftBack Flexi

### Integrated seat with opening to rear

(Accessory)

The seat's ergonomic design and flexible structure provide the best possible comfort to the user, while the opening to the rear provides the carer with a natural, more hygienic working method.

### Backrest Flexi

Our new backrest is actually three ribbons made of PUR which are designed to adapt to the back shape of the user. The ribbons are firmly attached to the frame and they can easily be cleaned/wiped off.

### Lockable armrests

The armrests can be flipped up or locked in their lower position. This is for safety when the chair is tilted.

### Plastic or metal

A plastic bedpan and guides are available as accessories for EASY. Note that guides are available in different versions to suit plastic or metal bedpans. See accessories

### Tilt function

EASY can be inclined approx 5 degrees forward and 30 degrees backward. This makes it easy to use when sitting down, getting up and for personal hygiene. The tilt function also improves the working environment and reduces strain on employees.

### Swing-away, detachable footrests

The footrests can be swung away and detached and there is also a heel strap provided as standard to prevent the patient's foot from slipping from the footrest during tilting.

### Adjustable seat height

An infinitely variable seat height from 53CM all the way to 93cm means the chair will fit every patient and reduce staff workload to a minimum.

### Directionally locking castor

The right front castor is directionally lockable for easy manoeuvrability. The image shows a double castor (accessory).

### Lockable castors

All castors are lockable except the directionally lockable castor. The image shows a double castor (accessory).

### Pillars

The telescoping aluminium pillars have an open design that allows water and dirt to simply run through.



### Easy control

Electrical height adjustment is via hard-wired controller with a high ingress protection rating (IP65) for excellent safety.

## Technical specifications

Part No.	Product	Width, external	Max. user weight	Seat height	Width between armrests	Seat depth
56-270	EASY SoftBack Comfort	65cm	150kg	53-93cm	49cm	45cm
56-270/LOW	EASY SoftBack Comfort Low	65cm	150kg	50-81cm	49cm	45cm
56-272	EASY SoftBack Flexi	65cm	150kg	53-93cm	49cm	45cm

Examples of possible parts combinations are shown below.