

Strong, safe and easy to manoeuvre – with electrical base widening and low legs

Eva600EE is a versatile mobile patient lift for all lifting situations, especially for heavier users weighing up to 270 kg/600 lbs. When used with Handicare's functional and easy-to-use SystemRoMedic™ lifting accessories, Eva600EE is the right choice for lifting in connection with transfers to and from the floor, chair or bed, as well as for horizontal lifting and gait training. This variant of Eva600EE features low legs and is especially suitable for use with lower beds.





Eva600EE, with low legs

- Extra strong and stable, SWL 270 kg/600 lbs
- Aluminium for low weight easy to manoeuvre
- Mast adjustable to three height positions large lifting range
- Height-adjustable handle bars comfortable and ergonomic working positions
- Stable aluminium sling bar with reliable safety latches
- Wider sling bar providing a more open lifting sling ideal for larger users
- Special version with low legs for use with lower beds
- Encapsulated twin castors in the rear easy rolling and maintenance-free
- Double front castors with a diameter of 35 mm/ 1.4"
- Electrical base widening
- Electrical emergency stop button readily available on the control box
- Electrical and manual emergency lowering
- · High capacity battery with long service life
- Built-in charging
- Quick guide on the lift
- Many accessories for safe, easy and individually adapted lifting of users

PRODUCT SPECIFICATION

Eva600EE. Art. No. 60100010

Mobile lift with electrical base widening and low legs, complete with SlingBar L and a hand control.

TECHNICAL INFORMATION

Safe working load (SWL)	270 kg /600 lbs
Material —	Aluminium
Total height	1330-2020 mm /52.4-79.5"
Base length	1400 mm /55.1"
Base width, outer dimensions	750-1150 mm /29.5-45"
Base width, inner dimensions	670-1180 mm /26.4-46.5"
Base widening	Electrical
Total weight	39 kg /86 lbs
Rear castors, diameter	100 mm /3.9"
Front castors, diameter	2 x 35 mm /1.4"
Lifting range	500-1810 mm /19.7-71.2"
Lifting speed without load	34 mm/s /1.3"/s
Emergency lowering	Manual and electrical
Battery capacity	2.9 Ah
Material, sling bar	Aluminium
Sling bar width	600 mm /23.6"
Protection class	IPX4



Eva600EE complies with requirements of the Medical Devices Directive for Class 1 products (MDD93/42/EEC).

The armrests accessory, AmbulationArm, can be employed by the user during gait training but they can also be used by the caregiver for manoeuvring the lift. Eva mobile lift fitted with AmbulationArm. New control box with diagnostics, e.g. lift



counter and overload

ACCESSORIES

Ergonomic and easy-to-use hand control with 4 buttons

To create optimal solutions for users with different needs, Eva600EE can be used in combination with lifting slings and many other clever and functional lifting accessories in Handicare's SystemRoMedic $^{\text{TM}}$ lifting range.

Art. No.	Product name	Description
80100010	AmbulationArm	_Armrests, for all Eva mobile lifts
70200001	SlingBar S	_Two-point aluminium sling bar, 35 cm/13.8"
70200003	SlingBar M	_Two-point aluminium sling bar, 45 cm/17.7"
70200042	SlingBarSpreader M	Side bars for SlingBar S/M/L
70200006	StretcherBar	For horizontal lifting
70200017	Sling bar RFL X4	_Four-point steel sling crossbar
70100002	Charder MHS2500	Scale, SWL 300 kg/660 lbs
70100003	Charder MHS2500	_Scale, SWL 400 kg/880 lbs
70200008	SwiftHook	For easy changing of sling bars
70200089	Hand control HB33-6 —	- 6 buttons, with service, battery status and

Other models and variants of Eva mobile lifts

Eva mobile lifts are available in several models and variants with different lifting capacity and with electrical or manual base-width adjustment. They are also available in special versions with low legs for use with lower beds.

Art. No.	Product name	Description
60100001	Eva400EE	Mobile lift with electrical base widening, SWL 180 kg/400 lbs
60100011	Eva400EM	Mobile lift with manual base widening, SWL 180 kg/400 lbs
61000002	Eva450EE	Mobile lift with electrical base widening, SWL 205 kg/450 lbs
60100006	Eva450EE	Mobile lift with electrical base widening, low legs, SWL 205 kg/450 lbs
60100007	Eva450EM	Mobile lift with manual base widening, low legs, SWL 205 kg/450 lbs
60100003	Eva600EE	Mobile lift with electrical base widening, SWL 270 kg/600 lbs
61000002 60100006 60100007	Eva450EE Eva450EE Eva450EM	Mobile lift with manual base widening, SWL 180 kg/400 lbs Mobile lift with electrical base widening, SWL 205 kg/450 lbs Mobile lift with electrical base widening, low legs, SWL 205 kg/450 lbs Mobile lift with manual base widening, low legs, SWL 205 kg/450 lbs

Please visit www.handicare.com for more information.



Maskinvägen 17, 972 54 Luleå, Sweden Tel: +46 (0)8-557 62 200 Fax: +46 (0)8-557 62 299 www.handicare.com info@handicare.se