



# Dyna-Form® PowerTurn Elite

Designed for challenging bariatric patients that have very little or no mobility, the Dyna-Form® Power-Turn Mattress System combines two clinically effective therapies: Lateral Rotation and True Low Air Loss in a portable system. A full, Left or Right Turn Mode is offered with either a 25° Small Turn or a 40° Big Turn.

## Features

Electrically operated back rest

Epoxy powder coat finish

Drop arms

## Benefits

Back rest recline from 85° to 60° to ensure optimal patient comfort

Easy to clean and favourable with Infection Control

Easy patient side transfer into the chair

## Technical Specifications

<b>SWL</b>	453kg/71 stone
<b>Risk Category</b>	Very High Risk
<b>Warranty</b>	1 year
<b>Dimensions</b>	91cm (W) x 203cm (L) x 22cm (H)
<b>Product Weight</b>	Mattress: 10kg Pump: 5kg