

# Vega505EE

SystemRoMedic

SWL 230 kg / 505 lbs



## Sturdy and safe, with an extremely large lifting range

Vega is a mobile patient lift adapted for safe and easy lifting of users weighing up to 230 kg/505 lbs. The lift is sturdy, stable and secure and, at the same time, light-weight for a steel lift. It is also relatively narrow and compact, which means that it can fit into and reach practically anywhere, even in tight spaces. In spite of the handy size, Vega offers an extremely large lifting range and, therefore, the lift can manage most lifting situations and requirements; both low and high lifting, seated or supine, to and from the floor, bed or chair, and gait training.

**Direct Healthcare  
Group**

Advancing Movement & Health®



## Unique design: the V-concept

Vega features a unique design with a lift arm that can be both raised to an almost upright position and lowered all the way down to the floor. In addition, the mast is split in two and the handle bars are integrated in the mast. This design, which we call the V-concept, offers a number of advantages, both for the caregiver and the user.

Vega features unique, ergonomic handle bars integrated in the mast. The grip-friendly handle bars offer multiple hand placement options and heights, allowing the caregiver to work more conveniently and comfortably and in an ergonomically correct manner. The design also provides intuitive properties, making the manoeuvring of the lift natural and easy.

Furthermore, the split mast, which signifies the V-concept, offers a new location, *in* the mast, for the actuator and control-box, which, in turn, provides more space and improved comfort, both for the user and the caregiver.

### Simplicity and safety

The ergonomic and easy-to-use hand control has clearly marked buttons for lifting and lowering and for electric base widening. All motion stops as soon as the buttons are released.

The emergency stop button is readily available on the control box. Vega also features both manual and electrical emergency lowering.

### Reliable power supply

The battery has a high capacity allowing for performance of many lifts per charge. This means that charging is needed less frequently. A sound signal can be heard when the battery level is low and an indicator on the control box shows when charging is in progress.



Ergonomic and easy-to-use hand control with 4 buttons.



New control box with diagnostics, e.g. lift counter and overload indicator.

### Extra equipment as standard

Vega is supplied complete and ready to use, fitted with a hand control and with the following accessories and extra equipment as standard:

- SlingBar M, two-point aluminium sling bar, standard width 450 mm/17.7"
- SwiftHook, accessory for SlingBar
- Built-in charging
- Encapsulated twin castors, Ø100 mm/3.9", with bumpers



for Vega



**Unique design with many advantages:**

- Stable and secure, still comparatively light-weight
- For users of up to 230 kg/505 lbs
- Convenient size, fits in and reaches practically anywhere
- Extremely large lifting range
- Excellent for both very low and very high lifting and gait training
- Unique design with intuitive handle bars integrated in the mast
- The design provides more space and comfort for both user and caregiver
- Easy and ergonomic to manoeuvre, prevents work-related injuries
- Quick guide on the lift
- Encapsulated twin castors - easy rolling and maintenance-free
- Electrical base widening
- Emergency stop button readily available on the control box
- Electrical and manual emergency lowering
- High capacity battery with long service life
- Built-in charging
- Many accessories for safe, easy and individually adapted lifting of users

*Pictured here is an earlier version of Vega506EE fitted with the old hand control and control box. This version is no longer available.*

