



QUATTRO[®] plus

Technical Specifications

POWER UNIT

Width	346mm
Height	258mm
Depth	156mm
Weight	4.0kg
Cycle time	16 mins

MATTRESS

Length*	1950mm
Width*	880mm
Depth	180mm
Weight	9.4kg
Inner construction materials	PU air cells
Cover material	PU coated stretch nylon
No. of cells	30
Max. user weight	200kg / 31stone
Min. user weight	None

Suitable for all patient types, the **Quattro[®] Plus** mattress system combines an active 1-in-4 cell cycle, with TISSUEgard, Ortho-Differential Support and a multi-stretch cover to provide patients with pressure relief to prevent and manage patients at risk of or with existing pressure ulcers.

The active 1-in-4 cell cycle supports 75% of the patient's body at all times, on groups of three inflated cells, whilst the fourth cell deflates sufficiently to redistribute the pressure and encourage tissue reperfusion. Additional benefits of the 1-in-4 system may include increased patient comfort, and less awareness of the support surface movement.

TISSUEgard™ enables the partial immersion and envelopment of the patient into the support surface, reducing the pressure on the patient's skin and decreasing the pressure differential between inflated and deflated cells, helping to reduce associated shearing forces.

Ortho-Differential Support provides firmer outer edges of the mattress to facilitate patient transfers and provide extra support, whilst creating a softer central area of the mattress.

*UK standard size