

Quick Selection Guide - Pressure Care Cushion Range

Direct Healthcare Group's cushion range offers a variety of seating solutions suitable for use from 'At Risk' to 'Very High Risk' patient groups. Utilising the highest quality materials the range offers optimum pressure redistribution through the use of innovative designs.

At Risk

Product Code	Cushion Name	Risk Category	Foam	Depth	Weight Guide	Postural Management	Contoured	Gel Inserts	Castellated	Low Profile
CUS0510004	Dyna-Tek® Profile	At Risk	RX39	2"	114kg/18 stone	✓	✓			✓
CUS0110001	Dyna-Tek® Pad Cushion	At Risk	RX36	4"	114kg/18 stone				✓	
CUS0110007	Dyna-Tek® Pad Cushion (Low Profile)	At Risk	RX36	2"	114kg/18 stone				✓	✓
CUS0210001	Dyna-Tek® Superior	At Risk	RX36 top, RX39 base	4"	127kg/20 stone				✓	
CUS0210000	Dyna-Tek® Superior Bari	At Risk	RX36 top, RX39 base	4"	318kg/50 stone				✓	
S43	Castle-Top Plus (Firm)	At Risk	CM60	3.5"	127kg/20 stone				✓	
S44	Castle-Top Plus (Medium)	At Risk	RX39	3.5"	102kg/16 stone				✓	
S45	Castle-Top (Firm)	At Risk	CM60	3"	127kg/20 stone				✓	
S46	Castle-Top (Medium)	At Risk	RX39	3"	102kg/16 stone				✓	

High Risk

Product Code	Cushion Name	Risk Category	Foam	Depth	Weight Guide	Postural Management	Contoured	Gel Inserts	Castellated	Low Profile
CUS0410001	Dyna-Flex® Cushion	High Risk	RX39 base, Viscoelastic top	3"	222kg/35 stone	✓				
CUS0610001	Dyna-Tek® Posture Visco	High Risk	RX39 base, Viscoelastic top	3"	178kg/28 stone	✓	✓			
CUS0310001	Dyna-Tek® Owl	High Risk	Viscoelastic	3"	152kg/24 stone	✓	✓			
C1	Community One	High Risk	Viscoelastic top & bottom, CM60 core	3"	95kg/15 stone	✓				
C2	Community Two	High Risk	Sunmate top & bottom, CM60 core	3"	114kg/18 stone	✓				
PC	Mercury Junior	High Risk	D5 base, RX39 overlay	2"	64kg/10 stone	✓	✓			✓
QM50	Mercury 50	High Risk	D5	2"	95kg/15 stone	✓	✓			✓
QM100	Mercury 100	High Risk	D5 base, RX39 overlay	3"	114kg/18 stone	✓	✓			
QM150	Mercury 150	High Risk	Reflex base, Viscoelastic insert	3"	102kg/16 stone	✓	✓			
QM175	Mercury 175	High Risk	D5 base, RX39 overlay, Visco insert	3"	114kg/18 stone	✓	✓			
PCX	Mercury Junior Xtreme	High Risk	Reflex core, Chip foam sides	2.25"	64kg/10 stone	✓	✓			✓
PCX50	Mercury 50 Xtreme	High Risk	D5	2"	95kg/15 stone	✓	✓			✓

Product Code	Cushion Name	Risk Category	Foam	Depth	Weight Guide	Postural Management	Contoured	Gel Inserts	Castellated	Low Profile
PC100X	Mercury 100 Xtreme	High Risk	D5 base, RX39 overlay	2"	114kg/18 stone	✓	✓			
PC150X	Mercury 150 Xtreme	High Risk	D5 base, Visco insert	3"	102kg/16 stone	✓	✓			
PC175X	Mercury 175 Xtreme	High Risk	D5 base, RX39 overlay, Visco insert	3.5"	114kg/18 stone	✓	✓			
QC1	Q-Care 1	High Risk	CM60 core, RX39 overlay	3"	102kg/16 stone	✓	✓			
S6	Cold Mould 2"	High Risk	Cold Mould	2"	133kg/21 stone	✓	✓			✓
S7	Cold Mould 3"	High Risk	Cold Mould	3"	146kg/23 stone	✓	✓			
S8	Cold Mould 4"	High Risk	Cold Mould	4"	152kg/24 stone	✓	✓			
S11	Qbi-Flex	High Risk	CM60 base, Pink visco core, Yellow visco base	3"	95kg/15 stone	✓	✓			

Very High Risk

Product Code	Cushion Name	Risk Category	Foam	Depth	Weight Guide	Postural Management	Contoured	Gel Inserts	Castellated	Low Profile
CUS1210001	Dyna-Form® Mercury Advance Cushion	Very High Risk	RX39 in air cells	4"	152kg/24 stone	✓			✓	
CUS0810001	Dyna-Tek® Intelligent Air	Very High Risk	RX39 in air cells	3"	152kg/24 stone	✓			✓	
CUS0810007	Dyna-Tek® Intelligent Air (Low Profile)	Very High Risk	RX39 in air cells	2"	152kg/24 stone	✓			✓	✓
CUS0710001	Dyna-Tek® Gel	Very High Risk	Viscoelastic	3"	159kg/25 stone	✓	✓	✓		
CUS0910000	Transflo	Very High Risk	RG50	3.5"	178kg/28 stone	✓		✓		
CUS0920000	Transflo Leisure	Very High Risk	RG50	2.5"	140kg/22 stone	✓		✓		✓
CUS0930000	Transflo Heavy Duty	Very High Risk	RG50	3.5"	203kg/32 stone	✓		✓		
CUS0970000	Transflo Contour	Very High Risk	RG50	3.5"	178kg/28 stone	✓	✓	✓		
CUS0940000	Transflo Bari	Very High Risk	RG50	3.5"	203kg/32 stone	✓		✓		
C3	Community Three	Very High Risk	Sunmate top & bottom, Pudgee core	3"	127kg/20 stone	✓				
CG	Community Gel	Very High Risk	Visco top & bottom, CM60 core	3"	95kg/15 stone	✓		✓		
CG2	Community Gel Two	Very High Risk	Sunmate top & bottom, CM60 core	3"	127kg/20 stone	✓		✓		

Product Code	Cushion Name	Risk Category	Foam	Depth	Weight Guide	Postural Management	Contoured	Gel Inserts	Castellated	Low Profile
QM200	Mercury 200	Very High Risk	D5 base, Visco overlay	3.5"	114kg/18 stone	✓	✓			
QM300	Mercury 300	Very High Risk	D5 base, Visco overlay, Gel insert	3.5"	114kg/18 stone	✓	✓	✓		
QM200X	Mercury 200 Xtreme	Very High Risk	D5 base, Visco overlay	3.5"	114kg/18 stone	✓	✓			
QM300X	Mercury 300 Xtreme	Very High Risk	D5 base, Visco overlay, Gel insert	3.5"	114kg/18 stone	✓	✓	✓		
S30	Q-Care 2	Very High Risk	CM60 core, RX39 overlay, Pudgee insert	3"	102kg/16 stone	✓	✓	✓		
QCINT1	Q-Care Integra 1	Very High Risk	CM65 core, RX39 overlay, Chipfoam sides	2"	102kg/16 stone	✓	✓			✓
QCINT2	Q-Care Integra 2	Very High Risk	CM60 base, RX39 overlay, Sunmate insert	2"	114kg/18 stone	✓	✓			✓
S9B	Sunmate 2"	Very High Risk	Sunmate	2"	95kg/15 stone	✓				✓
S9C	Sunmate 3"	Very High Risk	Sunmate	3"	114kg/18 stone	✓				
S10B	Pudgee 2"	Very High Risk	Pudgee	2"	76kg/12 stone	✓				✓
S10C	Pudgee 3"	Very High Risk	Pudgee	3"	95kg/15 stone	✓				
S1C	Laminar	Very High Risk	Sunmate base, Pudgee top	2"	127kg/20 stone	✓				✓

Product Code	Cushion Name	Risk Category	Foam	Depth	Weight Guide	Postural Management	Contoured	Gel Inserts	Castellated	Low Profile
S1B	Contour	Very High Risk	Sunmate base, Pudgee core, RX39 top	3"	191kg/30 stone	✓				
S1A	Super Contour	Very High Risk	Sunmate base, Pudgee core, RX39 top	3"	127kg/20 stone	✓			✓	
S78	Interum	Very High Risk	D5 base, Sunmate overlay, Gel insert	4.75"	114kg/18 stone	✓	✓	✓		
S79	Optimus	Very High Risk	D5 base, Pudgee overlay, Gel insert	4.75"	114kg/18 stone	✓	✓	✓		
S77	Qbi-Gel	Very High Risk	Gel top & bottom, Sunmate core	1.5"	191kg/30 stone	✓				✓

Direct Healthcare
Group

Advancing Movement & Health®



LIT-00045 Issue 1
Date: August 2018

Sales Offices

UK & Europe

Direct Healthcare Group Ltd.
Withey Court, Western Industrial Estate,
Cerphilly, UK,
CF83 1BF

T: +44 (0) 800 043 0881
info@directhealthcaregroup.com
www.directhealthcaregroup.com

Asia Pacific

Direct Healthcare Group PTY Ltd.
PO Box 562, Wembley,
Western Australia, 6913

T: +61 (0) 423 852 810
info@directhealthcaregroup.com.au
www.directhealthcaregroup.com