



Stool

Code **G12JDS02** Series **HUG**

Net weight kg 5 Dimensions [WxDxH] mm 410x370x445

In use weight capacity Vertical kg 150

Use and installation info

Use a soft, moist cloth soaked with liquid soap and/or other liquid sanitizing substances. DO NOT USE abrasive liquid detergents, chlorine bleach or strong acids that does not improve the aseptic aspect of the equipment and could lead to opacity and deterioration of the product.

Fasteners

Packaging information

Dimensions [WxDxH] mm 435x145x485 Weight kg 5,9

Available colours

Primary: WE Signal White, DD Graphite, PB Coral

Secondary: G9 Lime, P5 Aubergine, R2 Copper, WE Signal White, N2 Bright

