



Folding grab bar

Code **G18JCS09** Series **HUG**

Net weight kg 3,4 Dimensions [WxDxH] mm 100x850x120

In use weight capacity Vertical kg 150 Horizontal Kg 100

Use and installation info

Fasteners

Ponte Giulio suggest to use fasteners suitable to the wall characteristics on which the product will be installed.

Packaging information

Dimensions [WxDxH] mm 890x145x110 Weight kg 3,73

Available colours

Primary: C1 Chrome, PB Coral, DD Graphite, WE Signal White

Secondary: WE Signal White, C1 Chrome, PB Coral, D4 Dark Grey

