



# Tergo®Glide

**Tergo®Glide** is an ergonomic transfer aid that can be used in the transfer and repositioning of dependant patients, with the additional benefit of an integral, absorbent incontinence pad. The patented **Tergo®Glide** design provides a smooth, effective lateral movement, whilst preventing unwanted migration down the bed.

The practical and robust grips enable the caregiver to establish a firm, ergonomic working position; providing a reduction in physical loads associated with patient moving and handling. The large size and high safe working load of the **Tergo®Glide** provides a transfer solution that is suitable in a wide range of patient demographics. Whilst constant, in situ placement reduces disturbance of the patient and enables independent transfer if required.

## Features

- ① Low-friction base material with directional properties
- ② Top material is an absorbent incontinence pad
- ③ Practical and robust grips

## Benefits

- Provides easy lateral transfer and repositioning movements whilst preventing migration down the bed as the backrest is raised
- Provides the additional benefit of incontinence management within a transfer sheet, preventing the need for additional material *in situ*
- Provides a firm and secure hold of the pad to facilitate use of the product

## Technical Specifications

<b>Product Code</b>	<b>TRG0010001</b> (S-M) <b>TRG0010002</b> (L-XL)
<b>Weight Limit</b>	250kg/39 stone
<b>Standard Sizes</b>	88cm x 120cm (S-M) 105cm x 120cm (L-XL)
<b>Structure</b>	Integrated, multi-layered with handles
<b>Absorbency</b>	30dl/m <sup>2</sup>
<b>Warranty</b>	1 year