

VISKOA MATTRESS FIRE RETARDANT



The VISKOA fire retardant mattress has been developed in accordance with the following British Standards:

- BS EN 597-1:1995
- BS EN 597-2:1995
- BS EN 6807: 2006, section 9

The VISKOA fire retardant mattress is recommended for:

- The prevention of pressure ulcers, patients (from 30 to 150 Kg) presenting a medium to high risk of developing pressure ulcers in lying position.
- The treatment of existing pressure ulcers from stage 1 and 2 – according to clinical evaluation.

It is ideal for long-term institutions, nursing/care homes environments.

MAIN FEATURES

- Single piece mattress
- Single support
- Made from visco-elastic foam and HR foam
- Dimensions: 195 x 88 x 14 cm
- For patient up to 150 Kg
- Visco-elastic foam : 6 cm high – density: 60 Kg/m³
- High Resilient (HR) foam: 8 cm high – density: 30 Kg/m³

- Delivered with a Promust PU HD cover
- Packaging: Box of 1 or 10 units
- Warranty 3 years and 1 year for the **cover (only for communities)**

ADVANTAGES & BENEFITS

<p>Made from visco-elastic foam to improve the immersion of the body into the material thus allowing to reduce the pressure exerted on the skin</p>	<p>Safety and comfort for the patient → prevention of pressure ulcers for patients at high risk ; treatment of existing pressure ulcers from stage 1 to 2 according clinical evaluation</p>
<p>Fire Retardant Treatment to protect the patient against ignition sources</p>	<p>Safety for the patient when he is contact with fire or electrical sources</p>
<p>Promust PU HD cover (properties: waterproof, vapour resistant and stretchable) to ensure an optimal hygiene for the patient</p>	<p>Comfort for the patient</p>
<p>Packaging: box of 10 units available (pallet box)</p>	<p>Practical: mattresses are ready to be used Environmentally friendly: less packaging used</p>

WINDCARE
S O U P

WINDCARE
S O U P