

 **Dyna-Form®**
Air Suresse Alternating Mattress System



Dyna-Form Air Suresse - Alternating Mattress System

The box contains a completely assembled Alternating Mattress System consisting of:

A. Dyna-Form Air Suresse Mattress including Mattress Base, CPR Handle and Transportation Mode, Air Hoses, Quick Connectors.

B. Digital Control Unit

C. Power Cord

D. Carry Bag

E. User Manual


C.




D.



E.



Control Unit Panel (see page 6)



Transport Function - Rapid Release Handle (see page 9)

