



Dyna-Form™ Birthing Mattress

The **Dyna-Form Birthing Mattress** has been created for heavy use within the maternity environment. Its high frequency welded cover provides an excellent barrier from fluid ingress in an environment which could be an infection control risk for new mothers. As epidural procedures can often be lengthy, the foam design is most sympathetic to the entire birthing procedure, with pressure care at the forefront of its design. This product can be made to measure in order to fit most birthing bed designs. **Dyna-Form's** design quality comprises of a specialist high density CME foam surface, complete with a high frequency welded, multi-stretch, vapour permeable cover and non-slip base.

Features

| Features | Benefits |
|--|---|
| ① CME pressure reducing foam insert | Provides excellent comfort and pressure re-distribution properties |
| ② Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent) | Reduces shear & friction forces, whilst assisting infection control policies |
| ③ High frequency welded cover | Prevents fluid ingress & aids decontamination processes (can be laundered up to 95°C) |
| ④ Non-slip base | Prevents slippage, maintaining patient safety |
| ⑤ Fully replaceable parts | Ensures great value for money |

Technical Specifications

| | |
|---|--|
| Product Code | MAT/BIRTH/SPECIAL |
| Risk Category | High Risk |
| Weight Limit | 40 Stone / 254kg |
| Warranty | 2 Years |
| Dimensions | Made to measure |
| Foam Density | Blue – 39–40kg/m ³ Nominal Hardness – 200 Newtons |
| Fatigue Class | Very Severe |
| Fire Retardancy | BS7177; 2008 (Crib 5 Composite) |
| Weight | Approximately 10kg & 5kg respectively |
| Colour | Wine |
| Hypochlorite Solution compatible | 1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines) |