



The **Dyna-Form® Community** mattress can be used with patients considered to be up to that of “**High Risk**” of pressure ulcer development. Using modern foam designs, castellated on one side to aid pressure ulcer management, whilst providing excellent patient comfort; this mattress can be used without a turning programme, aiding moving & handling procedures. The **Dyna-Form® Community** mattress can work comfortably with any electric profiling bed frame. Its durable design uses only the very best CME foam & fabrics to provide a modern mattress complete with infection control enhancing cover.

Also Available in Pro-Formance® Strikethrough Resistant Cover.

Features

Benefits

| | | |
|-------|--|---|
| 1a 1b | Castellated foam (with tear duct cut) | 1a smaller castellated cells for heels and 1b larger cells sacrally |
| 2 | Single sided cut surface | Works without a turning programme |
| 3 | Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent) | Reduces shear & friction forces, whilst assisting infection control policies |
| 4 | Fully sewn or high frequency welded cover. (NB: welded cover option carries a small price premium) | Prevents fluid ingress & aids decontamination processes (can be laundered up to 95°C) |
| 5 | Fully replaceable parts | Ensures great value for money |

Technical Specifications

| | |
|---|---|
| Product Code | MAT/COM/198/88/15 |
| Risk Category | High Risk |
| Weight Limit | 22 Stone / 140kg |
| Warranty | 1 Year |
| Dimensions | (Length) 198cm x (Width) 88cm x (Height) 15cm – However Variable Sizes Available Upon Request |
| Foam Density | 33kg/m ³ Nominal Hardness – 120 Newtons |
| Fatigue Class | Very Severe |
| Fire Retardancy | BS7177; 2008 (Crib 5 Composite) |
| Weight | 10kg |
| Colour | Blue |
| Hypochlorite Solution compatible | 1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines) |