



The **Dyna-Form Jupiter** mattress can be used with patients considered to be up to that of “**High Risk**” of pressure ulcer development. Using modern foam designs, castellated on one side to aid pressure ulcer management, whilst providing excellent patient comfort; this mattress can be used without a turning programme, aiding moving & handling procedures. **Jupiter** can work comfortably with any electric profiling bed frame. Its durable design uses only the very best CME foam & fabrics to provide a modern mattress complete with infection control enhancing cover.

Also Available in Pro-Formance™ Strikethrough Resistant Cover.

Features

Benefits

1a 1b	Castellated foam (with tear duct cut)	1a smaller castellated cells for heels and 1b larger cells sacrally
2	Single sided cut surface	Works without a turning programme
3	Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)	Reduces shear & friction forces, whilst assisting infection control policies
4	Fully sewn or high frequency welded cover. (NB: welded cover option carries a small price premium)	Prevents fluid ingress & aids decontamination processes (can be laundered up to 95 °C)
5	Fully replaceable parts	Ensures great value for money

Technical Specifications

Product Code	MAT/JUP/198/88/15
Risk Category	High Risk
Weight Limit	26 Stone / 165kg
Warranty	3 Years
Dimensions	(Length) 198cm x (Width) 88cm x (Height) 15cm – However Variable Sizes Available Upon Request
Foam Density	34–36kg/m ³ Nominal Hardness – 125 Newtons
Fatigue Class	Very Severe
Fire Retardancy	BS7177; 2008 (Crib 5 Composite)
Weight	10kg
Colour	Light Blue
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)