

The **Dyna-Form Mattress Overlay** has been created for users considered to be **"At Risk"** of pressure ulcer development. To be used in conjunction with any standard mattress, **'Dyna-Form'** design quality comprises of a specialist "castellated cut" CME foam surface, complete with a sewn multi-stretch, vapour permeable cover and non-slip base. This product is ideal for use in the homecare environment.

	Features	Benefits	Technical Specifications	
1a 1b	CME castellated foam (with tear duct cut)	1a smaller castellated cells for heels and 1b larger cells sacrally	Product Code Risk Category	OLAY/PAD/198/88/9 At Risk
2	Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)	Reduces shear and friction forces, whilst assisting infection control policies	Weight Limit Warranty	26 Stone / 165kg 2 Years Foam / 2 Years Cover
3	Sewn or welded cover (NB: welded cover option carries a price premium)	Prevents fluid ingress & aids decontamination processes (can be laundered up to 95°C)	Dimensions	(Length) 198cm x (Width) 88cm x (Height) 9cm Available in Double.
4	Non-slip base (NB: carries a price premium)	Prevents slippage on top of mattress replacement, maintaining patient safety	Foam Density	(Variable sizes available on request) 34–36kg/m³ Nominal Hardness – 125 Newtons
5	Fully replaceable parts	Ensures great value for money	Fatigue Class	Very Severe
			Fire Retardancy	BS7177; 2008 (Crib 5 Composite)
			Product Weight	6kg
			Colour	Light Blue
			Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)