

The **Dyna-Tek Owl** Cushion is a uniquely designed postural management product, providing exceptional pressure re-distribution via its "Owl" contoured design. This product is specifically made for users considered to be at **"High Risk"** of pressure ulcer development, with consideration towards pelvic tilt aspects amongst other key areas. Designed predominantly to house the ischial tuberosities within it's circular surface, this product shares the patient's weight equally, in areas normally vulnerable to skin breakdown. The **Dyna-Tek Owl** provides the ultimate combination of pressure care and patient comfort.

	Features	Benefits	Technical Specifications	
1	Unique "Owl" surface design	Creates exceptional pelvic stability and ischial housing	Product Code Risk Category	CUSH/OWL/46/43/9 High Risk
2	Profiled foam surface	Provides excellent pressure re- distribution & important patient comfort	Weight Limit Warranty	24 Stone / 152kg 2 Years
3	Finest quality CME foam	Prevents patient bottoming out and aids patient transfer	Dimensions	(Length) 46cm x (Width) 43cm x (Height) 9cm (18" x 17" x 3.5") Variable Sizes Available Upon Request
4	Multi-stretch, vapour permeable cover (impregnated with an anti fungal agent)	Reduces shear & friction forces, whilst assisting infection control policies (can be laundered up to	Foam Density	CME Base – 38–40kg/m³ Nominal Hardness – 200 Newtons
5	Non-slip base	95°C)  Aids moving & handling procedures, providing ultimate patient safety	Fatigue Class Fire Retardancy Weight Colour	Severe BS7176; 2007 1.3kg Blue
6	Available with sag compensator	Suitable for wheelchair users	Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)