

CHNICALSPECIFICATION











General Use

- Specialist care units
- · General and care of the elderly wards
- Care home and hospice environments
- Domestic home use
- For those seated for lengthy periods
- For those at risk from pressure related problems
- For the semi and ambulant user

Features and Benefits

- Gas action or Electrically powered Tilt-in-Space, Back angle recline and Elevating legrest
- Three pressure relieving seat module
- Advanced Comfort Backrest
- Innovative pressure relieving Reflexion™ Foam on user contact surfaces
- Multi-positioning, integral headrest
- Adjustable seat height and seat length
- Integral, adjustable lateral supports
- Multi-directional, breathable Dartex material on all user contact surfaces

Pressure Relieving Seat Options (integral)

- Intelli-Gel® Integral seat module
- Transflo Gel seat module
- Reflexion™ Foam seat module

Accessories

- Height adjustable footboard
- Laminated tray with straps
- Cushioned or Universal Headrest
- Lap-Belts and Pelvic Positioners

Weight Limit

160kg/25 stone



Tilt -in-Space

















