

Coral

Commode/shower stool with back and arm supports

Stable support

Height-adjustable

Compact and light-weight

Ergonomic design

3002010030
Coral, commode/shower stool with back and arm supports



A compact and height-adjustable commode/shower stool for independent comfort and convenience

Coral commode/shower stool with back and arm supports is an assistive device allowing the user to sit down safely and conveniently during toileting or showering. It is supplied complete with toilet rails and can be used freestanding with a bedpan or bucket or standing over a toilet. The commode/shower stool provides a stable and height-adjustable support for active users with weak legs, low stamina and/or poor balance and who require a little more support during sit-to-stand. For maximum stability, the legs of the stool are fitted with non-slip end caps with drainage holes. The seat height is easily adjusted to ensure the perfect fit for every user and situation.

Coral

Molded toilet seat provides a good sitting position and a comfortable seating angle

The rust-free frame is fitted with a practical and comfortable molded toilet seat made from durable plastic with a smooth surface. The seat provides a good sitting position and the comfortable seating angle of one degree creates a safe situation and a pleasant experience during toileting and showering as the user will not slide off the seat. The lower parts in the back and front of the seat make it easier for the user to reach intimate body parts without standing up.

Back and arm supports provide extra support and security

For users with reduced strength and balance, sit-to-stand and stand-to-sit procedures can be difficult. This model of Coral features a firm, comfortable back support and fixed arm supports that provide the extra support needed during these movements. The back and arm supports also offer the user an extra sense of security while seated without hindering a more active and independent user.

Light-weight and compact, easy to fit

The commode/shower stool is light-weight and compact and fits easily into any wet room, also when space is limited. The design, with an opening in the back support, and the shape of the arm supports makes it easy for users with reduced strength or hand function to grip and to lift the stool.



Also available with shower seat

This model of Coral is also available with a solid shower seat (item no. 3002010130). Coral stools are also available without back and arm supports.



For more information, please visit our website, www.handicare.com, or contact your local Handicare partner.

PRODUCT SPECIFICATIONS

Art.no	Description
3002010030	Coral, commode/shower stool with back and arm supports, complete with toilet rails

FEATURES AND BENEFITS

- Stable and secure, ergonomic design
- Safe and very easy to use
- Compact, hygienic and attractive design
- Fits into any wet room, also when space is limited
- Very easy to assemble – no tools required
- Stainless steel and plastic, suitable for wet environments
- Easy maintenance
- Molded toilet seat provides a good sitting position and a comfortable seating angle
- Lower parts in the back and front facilitate intimate hygiene
- The design of the lower parts prevents the user from slipping through the seat
- Easily adjustable seat height; 420-570 mm, adjustable in 7 steps of 25 mm
- Firm yet flexible and comfortable back support
- Fixed arm supports for support during transfers and for security while sitting
- Back and arm supports are designed not to hinder a more active and independent user

TECHNICAL INFORMATION

Material, frame and legs: _____ Stainless steel
 Material, seat, back and arm supports: _____ GF 30 % PP plastic
 Seat height: _____ 420 - 570 mm / 16.5" - 22.4", adjustable
 Seat measurements, width x depth: _____ 490 mm x 410 mm / 19.3" x 16.1"
 Toilet opening, width x depth: _____ 210 mm x 280 mm / 8.3" x 11"
 Distance between arm supports: _____ 450 mm / 17.7"
 Height, back support: _____ 330 mm / 13"
 Width, back support: _____ 349 mm / 13.7"
 Height, arm supports (measured from seat): _____ 210 mm / 8.3"
 Length, arm supports: _____ 249 mm / 9.8"
 Total height: _____ 750 - 900 mm / 29.5" - 35.4", adjustable
 Outer measurements, width x depth: _____ 550 mm x 525 mm / 21.7" x 20.7"
 Free height under seat (for placing over a toilet): _____ 365 - 515 mm / 14.4" - 20.3"
 Free width under seat (for placing over a toilet): _____ 410 mm / 16.1"
 Weight: _____ 7 kg / 15 lbs
 Max. load on arm support: _____ 65 kg / 143.3 lbs per arm support
 Max. user weight: _____ 150 kg / 330 lbs

ACCESSORIES

Item No.	Product
LI02979.411	Bedpan, with lid and handle, Ø 300 mm, polished stainless steel
LI02979.430	Bucket, with lid and handle, Ø 300 mm, polished stainless steel
3002920402	Soft toilet seat, PUR, oval opening, with cut-outs for fixed arm supports
3002920403	Soft toilet seat, PUR, keyhole opening, with cut-outs for fixed arm supports
3002920501	Transportation seat, PUR, with cut-outs for fixed arm supports
3002920200	Calf support, textile, with Velcro



handicare

Maskinvägen 17,
 972 54 Luleå, Sweden
 Tel: +46 (0)8-557 62 200
 Fax: +46 (0)8-557 62 299
www.handicare.com