

Coral

Shower chair

Height adjustable

Stable support

Compact and light-weight

Ergonomic design



3002010140 Coral, shower chair

A comfortable and height-adjustable shower chair for increased independence and convenience

Coral shower chair is an assistive device allowing the user to sit down safely and conveniently during showering or at the washbasin. The shower chair provides a stable and height-adjustable support for less active users with weak legs, low stamina and/or poor balance who sit on the chair for longer periods of time and who require more support during sit-to-stand. For maximum stability, the legs of the chair are fitted with non-slip end caps with drainage holes. The seat height is easily adjusted to ensure the perfect fit for every user and situation.



handicare

Coral

Molded shower seat provides a good sitting position and a comfortable seating angle

The rust-free frame is fitted with a practical and comfortable molded shower seat made from durable plastic with a smooth surface. The seat provides a good sitting position and the comfortable seating angle of one degree creates a safe situation and a pleasant experience during showering as the user will not slide off the seat.

Backrest with folding armrests provides extra support and security

For users with reduced strength and balance, sit-to-stand and stand-to-sit procedures can be difficult. This model of Coral is fitted with a stable, comfortable backrest with folding armrests that provide maximum support during these movements. The back- and armrests also offer the user an extra sense of security while seated without hindering a more active and independent user.

The backrest frame is fitted with a comfortable and easy-to-care-for mesh textile cover offering soft and pliable support. It also dries very quickly. The depth of the backrest is adjustable with Velcro at the lower end. The armrests can be individually folded up and down, thereby both facilitating sideways transfers and providing more space for the user while also permitting easier access for the caregiver. The length of the armrests allows for maximum support during forward or backward transfers.

Light-weight and compact, easy to fit

The shower chair is light-weight and compact and fits easily into any wet room.

For more information, please visit our website, www.handicare.com, or contact your local Handicare partner.



FEATURES AND BENEFITS

- Stable and secure, ergonomic design
- Safe and very easy to use
- Compact, hygienic and attractive design, fits into any wet room
- Very easy to assemble – no tools required
- Stainless steel and plastic, suitable for wet environments
- Easy maintenance
- Molded shower seat provides a good sitting position and a comfortable seating angle
- Easily adjustable seat height; 420-570 mm, adjustable in 7 steps of 25 mm
- Comfortable and ergonomically shaped backrest provides soft and pliable support
- The textile cover for the backrest is easy to remove and washable at 60-85 degrees for optimal hygiene
- Folding armrests provide ergonomic support during transfers and security while sitting
- Back- and armrests are designed not to hinder a more active and independent user

PRODUCT SPECIFICATIONS

| Art.no | Description |
|------------|---------------------|
| 3002010140 | Coral, shower chair |

TECHNICAL INFORMATION

| | |
|--|--|
| Material, frame and legs: | Stainless steel |
| Material, seat and armrests: | GF 30 % PP plastic |
| Material, backrest, frame/cover: | GF 30 % PP plastic / textile, polyester mesh |
| Seat height: | 420 - 570 mm / 16.5" - 22.4", adjustable |
| Seat measurements, width x depth: | 490 mm x 460 mm / 19.3" x 18.1" |
| Backrest, height x, width: | 520 mm x 490 mm / 20.5" x 19.3" |
| Height, armrests (measured from seat): | 230 mm / 9" |
| Length, armrest: | 400 mm / 15.7" |
| Total height: | 940 - 1090 mm / 37" - 42.9", adjustable |
| Outer measurements, width x depth: | 620 mm x 630 mm / 24.4" x 24.8" |
| Weight: | 10 kg / 22 lbs |
| Max. load on armrest: | 75 kg / 165.3 lbs |
| Max. user weight: | 150 kg / 330 lbs |

ACCESSORIES

| Item No. | Product |
|------------|---|
| 3002920500 | Transportation seat, PUR |
| 3002920200 | Calf support, textile, with Velcro |
| 3002921300 | Safety belt, for hip or chest |
| 3002920900 | Front support, aluminium/PUR |
| 3002921100 | Armrest lock, prevents folding function |



handicare

Maskinvägen 17,
972 54 Luleå, Sweden
Tel: +46 (0)8-557 62 200
Fax: +46 (0)8-557 62 299
www.handicare.com