

# Coral

Stools

*Stable support*

*Height-adjustable*

*Compact and light-weight*

*Ergonomic design*

3002010100  
Coral, shower stool



3002010000  
Coral, commode/shower stool



## **Compact and height-adjustable stools for independent comfort and convenience**

The Coral stool is an assistive device allowing the user to sit down safely and conveniently during showering or at the washbasin. The stool is available in two variants; with shower seat and with toilet seat. Coral commode/shower stools provide stable and height-adjustable support for more active users with weak legs, low stamina and/or poor balance and who require a little bit of support and an opportunity to sit down, for example in the shower. For maximum stability, the legs of the stools are fitted with non-slip end caps with drainage holes. The seat height is easily adjusted to ensure the perfect fit for every user and situation.

# Coral

## Two variants; with shower seat and with toilet seat

The Coral stool is available in two variants. Both variants feature a rust-free frame fitted with a practical and comfortable molded seat made from durable plastic with a smooth surface. The shower stool has a shower seat and the commode/shower stool has a toilet seat. Both seats provide a good sitting position and a comfortable seating angle of one degree which creates a safe situation and a pleasant experience during showering and toileting as the user will not slide off the seat. The lower parts in the back and front of the toilet seat make it easier for the user to reach intimate body parts without standing up.

The commode/shower stool is supplied complete with toilet rails and can be used freestanding, with a bedpan or bucket, or standing over a toilet.

## Light-weight and compact, easy to fit

Coral stools are light-weight and compact and fit easily into any wet room, also when space is limited.

## Also available with back and arm supports

Coral stools are also available with back and arm supports.



3002010130

3002010030

For more information, please visit our website, [www.handicare.com](http://www.handicare.com), or contact your local Handicare partner.

## PRODUCT SPECIFICATIONS

Art.no	Description
3002010100	Coral, shower stool
3002010000	Coral, commode/shower stool, complete with toilet rails

## FEATURES AND BENEFITS

- Stable and secure, ergonomic design
- Safe and very easy to use
- Compact, hygienic and attractive design
- Fits into any wet room, also when space is limited
- Very easy to assemble – no tools required
- Stainless steel and plastic, suitable for wet environments
- Easy maintenance
- Molded seat provides a good sitting position and a comfortable seating angle
- Lower parts in the back and front of the toilet seat facilitate intimate hygiene
- The design of the lower parts prevents the user from slipping through the seat
- Easily adjustable seat height; 420-570 mm, adjustable in 7 steps of 25 mm

## TECHNICAL INFORMATION

Material, frame and legs: \_\_\_\_\_ Stainless steel  
 Material, seat: \_\_\_\_\_ GF 30 % PP plastic  
 Seat height: \_\_\_\_\_ 420 - 570 mm / 16.5" - 22.4", adjustable  
 Seat measurements, width x depth: \_\_\_\_\_ 490 mm x 410 mm / 19.3" x 16.1"  
 Toilet opening, width x depth: \_\_\_\_\_ 210 mm x 280 mm / 8.3" x 11"  
 Total height: \_\_\_\_\_ 420 - 570 mm / 16.5" - 22.4", adjustable  
 Outer measurements, width x depth: \_\_\_\_\_ 490 mm x 525 mm / 19.3" x 20.7"  
 Free height under seat (for placing over a toilet): \_\_\_\_\_ 365 - 515 mm / 14.4" - 20.3"  
 Free width under seat (for placing over a toilet): \_\_\_\_\_ 410 mm / 16.1"  
 Weight: \_\_\_\_\_ 6 kg / 13 lbs  
 Max. user weight: \_\_\_\_\_ 150 kg / 330 lbs

## ACCESSORIES

Item No.	Product
LI02979.411	Bedpan, with lid and handle, Ø 300 mm, polished stainless steel
LI02979.430	Bucket, with lid and handle, Ø 300 mm, polished stainless steel
3002920400	Soft toilet seat, PUR, oval opening
3002920401	Soft toilet seat, PUR, keyhole opening
3002920500	Transportation seat, PUR
3002920200	Calf support, textile, with Velcro



**handicare**

Maskinvägen 17,  
 972 54 Luleå, Sweden  
 Tel: +46 (0)8-557 62 200  
 Fax: +46 (0)8-557 62 299  
[www.handicare.com](http://www.handicare.com)