

**StandUp**  
SystemRoMedic™

*The perfect  
solution for  
safe and quick  
lifts from the  
floor*



**StandUp** is a lift chair and a new kind of assistive device in the SystemRoMedic™ range. Using StandUp, a single caregiver can quickly, safely and easily assist a person lying on the floor, unable to get up on their own. In a secure and comfortable way, StandUp lifts the user up into a seated position, ready for rising or transferring to a wheelchair. StandUp is both compact and light-weight. It weighs only 11.5 kg/25 lbs, which is less than half compared to a portable mobile lift. At the same time, it is strong and sturdy enough to manage lifting users weighing up to 160 kg/352 lbs.



# StandUp

## - the perfect solution for safe and quick lifts from the floor

StandUp is very convenient to always have close at hands. It fits into all cars and is easy to carry with the attached carrying strap over the shoulder. StandUp is suitable for use in all environments and settings, even in small spaces where a mobile lift is hard to fit. StandUp operates silently, no assembly is required and the user instructions are always available right on the product.

StandUp is an assistive device that actually gets used, thereby providing both safe working conditions for health care staff and good care for the users. It has been specially developed for use in homecare, nursing homes and special institutions but it is equally useful also in hospitals.



Easy to use, instructions on the seat.

Compact and light-weight.

StandUp unfolds effortlessly in three simple steps and is immediately ready for use.



Quick, secure and comfortable lifts from the floor to a seated position.

Fold out the legs and fold down the seat.

Carrying strap for easy transport.

### PRODUCT SPECIFICATION:

StandUp, lift chair with power cable \_\_\_\_\_ Art. No. 60800001  
 StandUp, lift chair with battery \_\_\_\_\_ Art. No. 60800002

### ACCESSORIES:

Spacious Bag, large bag with wheels, for convenient transport and storage of StandUp and other assistive devices \_\_\_\_\_ Art. No. 00951

### TECHNICAL INFORMATION

Measurements, folded: \_\_\_\_\_ 45 x 74 x 19 cm/17,7 x 29 x 7,4"  
 Weight: \_\_\_\_\_ 11,5 kg/25 lbs  
 Seat, height: \_\_\_\_\_ 5-44 cm/2-17,3"  
 Seat, width x depth: \_\_\_\_\_ 43 x 39,5 cm/16,9" x 15,6"  
 Lifting speed: \_\_\_\_\_ 7,6 mm/0,3" per sec without load  
 Input: \_\_\_\_\_ 230V AC/50Hz/180VA int. 10 % max 2/18 min  
 Motor: \_\_\_\_\_ DC 24V, 4,5 A. IP51  
 Duty: \_\_\_\_\_ max operating time = 10 % at max 2 min continuous operation  
 Push: \_\_\_\_\_ 6000N  
 Noise level: \_\_\_\_\_ 48dB(A); measuring method DS/EN ISO 3746, without load  
 Material: \_\_\_\_\_ Aluminium  
 Protection class: \_\_\_\_\_ IP X4  
 Expected lifetime: \_\_\_\_\_ 10 years  
 Safe working load: \_\_\_\_\_ 160 kg/352 lbs



- Compact and light-weight compared to portable mobile lifts and inflatable lift-chairs and cushions
- Incredibly convenient; all-in-one, always at hands, always ready to use
- Quick, secure and comfortable assistance for users who have fallen to the floor
- Lifts from the floor to a seated position
- For all environments and settings
- Safe and easy to use, user instructions available on the product
- Prevents work related injuries
- An economical and cost-saving solution
- For users weighing up to 160 kg/352 lbs
- Available in two variants; with power cable and with battery
- Both variants are supplied complete with carrying strap and an extra long version of SitWalk (6085) for helping the user get up on the seat

SystemRoMedic™ is a complete solution for easy transfers. The concept encompasses four product areas that meet most needs in different transfer situations.

**Transfer**, between two locations. **Positioning**, within the same location.

**Support**, for mobility. **Lifting**, both manual and mechanical.



**handicare**

82 First avenue  
 Pensnett Estate,  
 Kingswinford  
 West Midlands DY6 7FJ  
 T: 01384 408700  
 www.handicare.co.uk  
 enquiries@handicare.co.uk