

# RiseBasic300M™

SystemRoMedic™



## Safe and easy lifting in all situations and for users weighing up to 138 kg/300 lbs

- SWL 138 kg/300 lbs
- Manual transfer in rail
- Modern, elegant and hygienic design
- Unusually small measurements and extremely low weight
- For all rail systems, in all settings and in all types of rooms
- Large lifting area and possibility of both very low and very high lifting
- In all situations and for all user needs
- Very easy to use, providing safe working conditions and allowing the caregiver to remain close to the user throughout the entire lifting procedure
- Moves silently and gently, soft start and stop
- Ergonomic and easy-to-use hand control
- On/Off-button on the hand control prevents unauthorized use of lift functions
- Switch manoeuvring, for manoeuvring of turntables using the hand control, is available as an accessory
- Manual emergency stop
- Electrical and manual emergency lowering
- Overload protection
- Over speed protection
- Slack strap switch
- Efficient strap control preventing accidental twisting of the lift strap
- High capacity, long lasting batteries and an energy-saving electric motor
- Easy and convenient hand control charging
- A unique charging solution, with one trolley for two charging options (in-rail/end point charging) is available as an accessory
- Discreet information panel
- Easy to service
- Diagnostics, a special software available as an accessory allows a PC to access further information about the ceiling lift unit
- Many accessories for safe, easy and individually adapted lifting of users in most transfer situations

# RiseBasic300M™

For professional and highly functional care environments with a harmonic and human atmosphere.



### Easy and convenient hand control charging is standard

The charging station for the ceiling lift unit is integrated in a small and discreet holder and parking place for the hand control. The charging station is easily installed on the wall, in convenient proximity of the ceiling lift unit's parking place.



## PRODUCT SPECIFICATION

| Item No. | Product       |
|----------|---------------|
| 50100039 | RiseBasic300M |

RiseBasic300M is supplied with hand control and charger. Sling bar ordered separately.

## ACCESSORIES

| Item No.  | Product                    | Description  |
|-----------|----------------------------|--|
| 70200075  | QuickTrolleyUniversal 360  | Quick release trolley, turns 360 degrees                         |
| 70200062  | In-rail/end point charging | Two options for charging   |
| 70200066  | Switch manoeuvring         | For manoeuvring of turntables using the hand control             |
| 70200048  | Software for diagnostics   | Software   |
| 50400048  | Parking Placard            | For convenient storing of sling bars, manuals and lifting slings |
| 50400332  | Parking Placard BariBar    | For convenient storing of BariBar and lifting slings             |
| 70200071  | SlingBar XS                | 250 mm / 9.8"  |
| 70200001  | SlingBar S                 | 350 mm / 13.8"   |
| 70200002  | SlingBar M                 | 450 mm / 17.7"   |
| 70200003  | SlingBar L                 | 600 mm / 23.6"   |
| 70200042  | SlingBarSpreader M         | Side bars for SlingBar S/M/L                                     |
| 70200072  | BariBar                    | Four-point sling bar for ceiling lift units, SWL 454 kg/1000 lbs |
| 70200017  | Sling bar RFL X4           | Four-point sling crossbar, SWL 300 kg/660 lbs                    |
| 70200053- |                            |  |
| 70200058  | Lift strap extension       | Extension straps, 6 lengths: 370-970 mm / 14.6"-38.2"            |

## TECHNICAL INFORMATION

|  |  |
|--|--|
| Safe working load, RiseBasic300M:      | 138 kg / 300 lbs                       |
| Material, chassis:                     | Aluminium                              |
| Material, cover:                       | Non flammable ABS plastic              |
| Measurements, length x width x height: | 266 x 151 x 156 mm / 10.5 x 5.9 x 6.1" |
| Weight:                                | 6.8 kg / 15 lbs                        |
| Transfer in rail:                      | Manual                                 |
| Hand control:                          | Electrical, 4 buttons                  |
| Lifting range:                         | 2400 mm / 94.5"                        |
| Lifting speed, without load:           | 30 mm / 1.2" per sec                   |
| Lifting speed, with load:              | 20 mm / 0.8" per sec                   |
| Sound level, with load:                | 51 dB(A)                               |
| Efficient strap control:               | Standard                               |
| Emergency lowering:                    | Electrical and manual                  |
| Emergency stop:                        | Manual                                 |
| Batteries:                             | 2 x 12V 2.9Ah                          |
| Hand control charging:                 | Standard                               |
| In-rail charging:                      | Optional, accessory                    |
| End point charging:                    | Optional, accessory                    |
| Protection class, lift unit:           | IP20                                   |
| Protection class, hand control:        | IPX4                                   |



**CE** RiseBasic330M meets the requirements of the Medical Devices Directive for Class 1 products 93/42/EEC.

Please visit [www.handicare.co.uk](http://www.handicare.co.uk) for more information.



# handicare

82 First avenue  
Pensnett Estate,  
Kingswinford  
West Midlands DY6 7FJ  
T: 01384 408700  
[www.handicare.co.uk](http://www.handicare.co.uk)  
[enquiries@handicare.co.uk](mailto:enquiries@handicare.co.uk)

SystemRoMedic™ is a complete solution for easy transfers. The concept encompasses four product areas that meet most needs in different transfer situations.

**Transfer**, between two locations. **Positioning**, within the same location.

**Support**, for mobility. **Lifting**, both manual and mechanical.