

# ReTurn7600™

SystemRoMedic™



*Patented\* design  
For larger and heavier  
users weighing up to  
205 kg / 450 lbs  
Activates the user  
Eliminates heavy lifting  
for the caregiver*

## **Safe and active sit-to-stand and transfer – without any heavy lifting**

ReTurn7600 is an easy-to-use and incredibly effective assistive device for sit-to-stand and transfer, but it is also a rehabilitation device. ReTurn7600 enables safe and active sit-to-stand and transfer to or from bed, wheelchair or toilet. ReTurn7600 can also be used to facilitate repositioning farther back in a chair or wheelchair. For the caregiver, ReTurn7600 enables easy transfers without any heavy lifting. For users who still retain some degree of strength, ReTurn7600 stimulates the natural pattern of movement during sit-to-stand, while at the same time strengthening the muscles and functional ability. This makes it the ideal assistive device in a rehab unit, and a user who uses ReTurn7600 at home can normally live at home longer with the help of only one caregiver.



# ReTurn7600



## Product specification

Item No.	Product	Description
7600	ReTurn7600	Transfer platform, complete with wide rising ladder



## Different models of ReTurn provide solutions for all situations and users

The basic principle behind ReTurn is that the user's own strength and ability must be utilized and improved. All ReTurn models therefore help to activate the user and encourage the user's natural patterns of movement.



More information about the different models of ReTurn and all accessories can be found at [www.handicare.co.uk](http://www.handicare.co.uk)

## Easy-to-use and very versatile

ReTurn7600 is a manual sit-to-stand and transfer aid that is very easy to understand and use, both for the user and the caregiver. ReTurn7600 is also a versatile product that can be used in several different ways to safely and efficiently assist a user during sit-to-stand and positioning in a seated position, as well as for short transfers indoors. In addition, ReTurn7600 is also an excellent training device, making it an ideal aid to the rehabilitation process.

ReTurn7600 has been developed for users weighing max. 205 kg/450 lbs. This model therefore features a wider base plate, providing more space for larger and heavier users. ReTurn7600 also has lower-leg supports that are adapted for these users, and a rising ladder that is both higher and wider than the rising ladders on ReTurn7500 and ReTurn7400.

## Incredibly easy to manoeuvre

ReTurn7600 has a unique, patented\* design which makes the transfer platform incredibly easy to manoeuvre. ReTurn7600 rolls, turns and pivots smoothly and easily on all surfaces, even in confined spaces. The device is easily positioned close to a bed, wheelchair or toilet, and moving ReTurn7600 requires minimal exertion, even with a very large and heavy user.

## Stable and secure base plate

A low centre of gravity and a broad underlying support surface, with wheels at all corners, make the base plate stable and safe. Two large, fixed and lockable centre wheels make ReTurn7600 easy to roll and directionally stable, while the castors at the corners make it very easy to manoeuvre.

## The rising ladder allows space for a natural pattern of movement and closeness

ReTurn7600 is equipped with a rising ladder that is fully open at the bottom, while the upper section allows many possibilities for grasping the ladder, including a central 'power grip'. ReTurn's rising ladder leans forward, away from the user, allowing the user plenty of space to get really close to the assistive device, so as to be able to 'climb' and/or pull himself up while rising to a standing position. The user is also supported by the ladder during standing training and transfer. The design allows plenty of work space also for the caregiver, who can remain very close and help to support the user from different directions.

The lower-leg supports can be adjusted both vertically and laterally. They can also be angled, making it easy to individually adapt the settings according to the user's height and other factors. This function also makes it possible for several individuals to use the same ReTurn7600.

ReTurn7600 is easy to assemble and disassemble for transport and storage.

## Technical information

SWL	205 kg / 450 lbs
Overall height	1210 mm / 47.64"
Base plate, length x width x height	795 mm x 672 mm x 45 mm / 31.30" x 26.46" x 1.77"
Weight	26 kg / 57.32 lbs

\* Patent: SE-0300528-7

Registered Community Design: RCD-529631-1

## Accessories

By combining ReTurn7600 with different accessories the user can gain additional support and a greater sense of security during sit-to-stand, standing training and transfer. The choice of accessories depends on the situation and the user's needs and functional ability.

Item No.	Product	Description
6031--6038	ReTurnBelt**	Size XXS--XXXL
6233--6236	ReTurnBelt, disposable, non-woven	Size S--XL
6133--6136	ReTurnBelt, wipeable	Size S--XL
6061	SupportStraps	Accessory for belts
7216	Padding for leg support	For ReTurn7600

\*\* Patent: SE-0401517-8



**handicare**

82 First avenue  
Pensnett Estate,  
Kingswinford  
West Midlands DY6 7FJ  
T: 01384 408700  
[www.handicare.co.uk](http://www.handicare.co.uk)  
[enquiries@handicare.co.uk](mailto:enquiries@handicare.co.uk)

SystemRoMedic™ is a complete solution for easy transfers. The concept encompasses four product areas that meet most needs in different transfer situations.

**Transfer**, between two locations. **Positioning**, within the same location.

**Support**, for mobility. **Lifting**, both manual and mechanical.