

Solutions for infection control

SystemRoMedic™



Disposable and wipeable assistive devices facilitate infection control efforts

Handicare's SystemRoMedic™ range includes an increasing selection of disposable and wipeable transfer aids made from various, strong, safe and convenient materials with many practical features. The main common nominator for all of these devices is that they provide for optimal hygiene, thereby effectively contributing to the prevention of the spread of multi-resistant and other types of infectious bacteria.

Personal transfer aids for short-term use provide optimal hygiene and safety

Handicare's SystemRoMedic™ range includes an increasing selection of personal transfer aids for short-term use. Most of these are made from a sturdy and durable non-woven (polypropylene) material, but we are also using a soft and breathable but still durable cotton/polyester blend. One of our most recent transfer aids in this category, the disposable EasySlide sliding mat, is made from a thin, pliable and environmentally friendly polyethylene (PE) plastic. This is a strong and durable material which, when incinerated, gives off only

carbon dioxide and water. Used sliding mats can also be recycled as soft plastic.

Common features and benefits:

- Intended for a single user and for short-term use
- Instead of laundering, the product is discarded after use, when it becomes soiled or when the user no longer needs it
- Prevents the spread of multi-resistant and other types of infectious bacteria

DisposableBelt™

SystemRoMedic™



For short-term use. Support for sitting transfers, sit-to-stand, standing and walking.

Item No.	Product	Size	Unit
6054	DisposableBelt, non-woven, one size	Medium/Large	5 pcs



ReTurnBelt™

SystemRoMedic™



ReTurn7500 and ReTurnBelt in polyamid/nylon

The patented* SystemRoMedic™ ReTurnBelt is an excellent complement to Handicare's ReTurn transfer platforms. For users with weak legs and requiring extra support.

EasySlide, disposable™

SystemRoMedic™



EasySlide on a roll – a personal sliding mat for short-term use. For turning and repositioning in bed, for supine and seated transfers, for repositioning in a seated position and for application of lifting slings. Available on a roll and in a practical single-pack.

Item No.	Product	Size	Unit
2040	EasySlide, disposable, on a roll	50 x 70 cm	2 x 60 pcs
2045	EasySlide, disposable, on a roll, incl practical wall- and tabletop dispenser in steel	50 x 70 cm	2 x 60 pcs
2041	EasySlide, disposable, single-packed	80 x 70 cm	140 pcs
2043	EasySlide, disposable, single-packed	120 x 70 cm	100 pcs
2042	EasySlide, disposable, single-packed	200 x 70 cm	70 pcs



Item No.	Product	Size	Unit
6233	ReTurnBelt, disposable, non-woven	S	5 pcs
6234	ReTurnBelt, disposable, non-woven	M	5 pcs
6235	ReTurnBelt, disposable, non-woven	L	5 pcs
6236	ReTurnBelt, disposable, non-woven	XL	5 pcs



*) Patent: SE-0401517-8



HighBackSling™

SystemRoMedic™

Standard model for users who need a lot of support

Handicare's disposable HighBackSling is designed for short-term use in most common lifting situations, e.g. sitting transfers between bed and wheelchair, and when the user needs extra head support. HighBackSling can also be used for lifting from a recumbent position on the floor.

BariSling™

SystemRoMedic™

Special model for extremely large and heavy users

Our disposable BariSling is a lifting sling that is designed for seated lifting of users weighing 200 kg/440 lbs or more.

Common features and benefits:

- Multiple uses during one user's hospital stay
- Sitting transfers (from bed to chair, chair to chair, chair to bed, etc.)
- Infection control (use with patients on contact, droplet and airborne precautions)
- Compatible with lifts and sling bars from other manufacturers
- Bariatric option with up to 450 kg/1000 lbs capacity is available



Item No.	Product	Size	Unit
40890004	HighBackSling, disposable, non-woven	S	5 pcs
40890006	HighBackSling, disposable, non-woven	M	5 pcs
40890007	HighBackSling, disposable, non-woven	L	5 pcs
40890008	HighBackSling, disposable, non-woven	XL	5 pcs
40890009	HighBackSling, disposable, non-woven	XXL	5 pcs



Size S / M / L / XL: SWL 200 kg / 440 lbs
Size XXL: SWL 300 kg / 660 lbs

47091010	BariSling, disposable, non-woven SWL 450 kg / 1000 lbs	XXXL	2 pcs
----------	---	------	-------



ThoraxSling™

SystemRoMedic™

For sit-to-stand and standing training in combination with MiniLift200, MiniLift160 and MiniLift125

ThoraxSling is a specially adapted sit-to-stand vest that is used in combination with Handicare's SystemRoMedic™ sit-to-stand lifts MiniLift200, MiniLift160 and MiniLift125. It is easy for the caregiver to apply, while giving the user an added sense of security during rising and standing training. ThoraxSling also facilitates toileting. ThoraxSling is available in two models; one which provides support behind the back (ThoraxSling) and one which also provides support under the seat (ThoraxSling, with seat support).

Item No.	Product	Size	Unit
45590004	ThoraxSling, disposable, non-woven	S	5 pcs
45590006	ThoraxSling, disposable, non-woven	M	5 pcs
45590007	ThoraxSling, disposable, non-woven	L	5 pcs
45590008	ThoraxSling, disposable, non-woven SWL 200 kg / 440 lbs	XL	5 pcs
45690004	ThoraxSling, with seat support, disposable, non-woven	S	5 pcs
45690006	ThoraxSling, with seat support, disposable, non-woven	M	5 pcs
45690007	ThoraxSling, with seat support, disposable, non-woven	L	5 pcs
45690008	ThoraxSling, with seat support, disposable, non-woven SWL 200 kg / 440 lbs	XL	5 pcs



MiniLift200 and ThoraxSling in polyester

ThoraxSling, with seat support, 45690006






StretcherSling

SystemRoMedic™

For horizontal lifting with StretcherBar

StretcherSling is a lifting sling that is used in combination with mobile lifts or stationary lift systems and Handicare's SystemRoMedic™ StretcherBar, a specially designed sling bar for horizontal lifting and transferring of users in a lying position.

StretcherSling has ten suspension points; five on either side of the user, providing evenly distributed support under the entire body.



Item No.	Product	Size	Unit
46592007	StretcherSling, disposable, non-woven SWL 300 kg / 660 lbs	One size	5 pcs



SafeHandlingSheet

SystemRoMedic™

For horizontal lifting and positioning with a two-point sling bar

Handicare's disposable SafeHandlingSheet is used for horizontal lifting, for moving from a lying position to sitting on the bedside and for positioning of a user in bed, without subjecting the caregiver to unnecessary strain.

SafeHandlingSheet is used in combination with mobile lifts or stationary lift systems and a two-point sling bar.

The disposable SafeHandlingSheet is made from a soft and breathable cotton/polyester material. The two largest sizes are designed without straps running across the sheet which can increase pressure.

- For supine transfers, turning and repositioning
- Supports prevention of pressure sores
- Skin safe, breathable cotton blend material provides improved comfort for the user
- Infection control (use with patients on contact, droplet and airborne precautions)
- Bariatric options are available



Item No.	Product	Size	Unit
47512006	SafeHandlingSheet, disposable, cotton	M	5 pcs
47512007	SafeHandlingSheet, disposable, cotton	L	5 pcs
47512008*	SafeHandlingSheet, disposable, cotton	XL	5 pcs
47512009*	SafeHandlingSheet, disposable, cotton	XXL	5 pcs

* Without straps



Size M / L: SWL 300 kg / 660 lbs

Size XL: SWL 400 kg / 900 lbs

Size XXL: SWL 570 kg / 1250 lbs



TurningSling[™]

SystemRoMedic

For turning of users in bed

TurningSling is used, in combination with a two-point sling bar mounted on a stationary or mobile lift, for turning of users in bed. TurningSling also facilitates hygiene, treatment of sores, changing of bed linen and dressing/undressing.

- Ideal for coccyx, sacral, and spinal wound care and dressing changes
- Ideal for toileting patients who use a bedpan and for brief changes and patient hygiene including bed baths and linen changes
- Infection control (use with patients on contact, droplet and airborne precautions)

Item No.	Product	Size	Unit
49592000	TurningSling, two loops, disposable, non-woven	One size	5 pcs
49592005	TurningSling, one loop, disposable, non-woven	One size	5 pcs

SWL 250 kg / 550 lbs



LimbSling

SystemRoMedic

For easy lifting of an arm or leg

LimbSling is designed to comfortably and conveniently assist lifting of a user's leg or arm, for example, to gain better access for treatment and dressing sores. LimbSling is very practical for large and heavy users. LimbSling is applied by placing it under the limb to be lifted and connecting the strap loops to the sling bar.

- Ideal for wound care, limb holding, lymphedema, ROM therapy and turning
- Infection control (use with patients on contact, droplet and airborne precautions)

Item No.	Product	Size	Unit
49292000	LimbSling, disposable, non-woven	One size	5 pcs

SWL 150 kg / 330 lbs



Strong, safe and comfortable wipeable material

Handicare also offers an increasing selection of wipeable transfer aids made from a strong, safe and comfortable material (nylon/thermoplastic polyurethane, TPU) with a smooth and soft surface and many practical features.

The surface of the entire product can be cleaned using disinfectants. In fact, we call this material "wipeable" but these products can even be submerged in water or used in the shower without getting wet.

EasyBelt™ SystemRoMedic



Support for sitting transfers, sit-to-stand, standing and walking.

Item No.	Product	Size	Unit
6113	EasyBelt, wipeable	S	1 pce
6114	EasyBelt, wipeable	M	1 pce
6115	EasyBelt, wipeable	L	1 pce
6116	EasyBelt, wipeable	XL	1 pce



ThoraxSling™ SystemRoMedic



Item No.	Product	Size	Unit
45550004	ThoraxSling, wipeable	S	1 pce
45550006	ThoraxSling, wipeable	M	1 pce
45550007	ThoraxSling, wipeable	L	1 pce
45550008	ThoraxSling, wipeable,	XL	1 pce
45550009	ThoraxSling, wipeable SWL 205 kg / 450 lbs	XXL	1 pce



ReTurnBelt™ SystemRoMedic




*) Patent: SE-0401517-8

Item No.	Product	Size	Unit
6133	ReTurnBelt, wipeable	S	1 pce
6134	ReTurnBelt, wipeable	M	1 pce
6135	ReTurnBelt, wipeable	L	1 pce
6136	ReTurnBelt, wipeable	XL	1 pce



EasyGlide
SystemRoMedic™



For supine transfers and for placing a lifting sling

The oval EasyGlide is a sliding board which is used for supine transfers from one bed or stretcher to another but also for positioning higher up in the bed. EasyGlide is particularly suitable for use on soft surfaces such as pressure-relieving mattresses. EasyGlide can also be used to assist placement of a lifting sling both in the recumbent and in seated positions.

- Very easy to use, quick-guide on the underside
- Black polyethylene with a smooth surface
- Hygienic; can be disinfected with 70 % ethanol, 45 % isopropanol or similar after use

Item No.	Product	Size	Unit
5040	EasyGlide, oval	85 x 39 cm	1 pce
5042	EasyGlide, oval	85 x 39 cm	2 pcs



EasyGlide
SystemRoMedic™



Facilitates positioning further back in chair and application of lifting slings

The smaller oval EasyGlide sliding boards are used in pairs to facilitate positioning further back in chair. Using EasyGlide, a more active user can push or shuffle themselves backwards in the chair, into a correct and comfortable position. The sliding boards are small and convenient, with soft, round shapes and smooth surfaces; features which make them very easy to place and remove and to always have readily available and close at hands..

- Very easy to use, quick-guide on the underside
- Black polyethylene with a smooth surface
- Hygienic; can be disinfected with 70 % ethanol, 45 % isopropanol or similar after use

Item No.	Product	Size	Unit
5045	EasyGlide, oval (mini)	45 x 20 cm	2 pcs



MiniStretcher
SystemRoMedic™

For safe and convenient lifting of patients onto ambulance stretchers

For ambulance staff no one situation is like the other. They are constantly subjected to new challenges when lifting patients onto ambulance stretchers; stairs, toilets, narrow spaces and heavy lifting in awkward positions are just a few examples. This is where Handicare's SystemRoMedic™ MiniStretcher comes to the rescue. It is easy to slide in place underneath the patient and has a solid construction with grip-friendly handles to secure safe heavy lifting. The flexible MiniStretcher adapts to the shapes of ambulance stretchers, so they can remain under patients for the entire journey. Naturally, the MiniStretcher fulfills all medical hygiene requirements. Since MiniStretcher can be exposed to both dirt and blood, it is waterproof and washable.

- Very easy to use, clear instructions are available on the stretcher
- The neck support provides good support under the patient's head
- Light padding on each side of a flexible plastic board makes MiniStretcher comfortable enough to be left under the patient during the ambulance transport
- The handles provide good gripping options and are placed near the friction points
- The durable surface is very easy to clean

Item No.	Product	Size	Unit
4017	MiniStretcher SWL 200 kg / 440 lbs	135 x 43 cm	1 pce



Disposable - but not for single use

All of our disposable transfer aids may be used for one user throughout their hospital stay. Handicare's disposable transfer aids are safe, cost effective and convenient options facilitating infection control efforts.

At the moment, the SystemRoMedic™ disposable range includes the following items, and others are on their way!



Wipeable transfer aids in several materials

At the moment, the SystemRoMedic™ wipeable range includes the following items, and others are on their way!



For updates and more information, please visit our website, www.handicare.com.



Maskinvägen 17,
972 54 Luleå, Sweden
Tel: +46 (0)8-557 62 200
Fax: +46 (0)8-557 62 299
www.handicare.com
info@handicare.se

SystemRoMedic™ is a complete solution for easy transfers. The concept encompasses four product areas that meet most needs in different transfer situations.

Transfer, between two locations. **Positioning**, within the same location.

Support, for mobility. **Lifting**, both manual and mechanical.