

BASIC TECHNICAL SPECIFICATION

DESCRIPTION	VALUE
internal dimension	200x86/90cm
external dimension	104x221cm
clearance of undercarriage	13cm with 125mm Tente Linea castors
clearance of undercarriage	15,5cm with 155mm Tente Linea castors
hi-lo adjustment	28-80cm with 125mm Tente Linea castors
hi-lo adjustment	30-82cm with 155mm Tente Linea castors
bed extension	10, 22, 32cm
height of siderails	40cm (split siderails); 39cm (3/4 siderails)
max. mattress height	18cm with split siderails
max. mattress height	17cm/27cm with extender with 3/4 siderails
safe working load	250kg
max. patient weight	185kg
weight of bed	145kg according to accessories
backrest angle	70°
thigh rest angle	34°
calf rest angle	35°
TR/ATR angle	14°
Power supply	230 V; +/- 10%
Maximum input	max 1,6 A, 370 VA
Water and dust protection level	IP X4
Linet is a company certified under European norms and quality systems:	ENISO 9001:2008, EN ISO 13485:2009, EN ISO 14001:2005
The products are produced under European norms for medical technology and hospital beds:	EN 60601-1-1; EN 60601-1-2; EN 60601-1-4; EN-60601-1-6; EN 601-2-52; 93/42/EEC

BASIC TECHNICAL DRAWINGS

