


TECHNICAL SPECIFICATION

kirton

ULTIMARANGE




Mini Ultima



Seat heights are taken with the seat set at 0° tilt (horizontal).

Midi Ultima



Seat heights are taken with the seat set at 0° tilt (horizontal).

Dimensions in mm

General Use

- For use at home
- For the semi and non-ambulant user
- For those with moderate to severe seating problems
- For children and small adults

Features and Benefits

- Two sizes: Mini and Midi
- Manual Tilt-in-Space and Legrest elevation
- Three seat widths: 250mm, 300mm, 350mm
- Three pressure relieving seat options
- Adjustable seat length
- Cushioned Headrest
- Height adjustable, flip-up footrest
- Built-in lateral supports
- Contoured back with headrest
- Durable 75mm twin wheels, rear braked

Pressure Relieving Seat Options (integral)

- Intelli-Gel® Integral seat module
- Transflo Gel seat module (Midi Ultima only)
- Reflexion™ Foam seat module

Back Style Options

- Advanced Comfort Backrest
- Reflexion™ Foam Backrest

Accessories

- Laminated tray
- Foot Sandals
- Adjustable headrest
- Lap-Belt
- Midline Vest
- Pelvic Positioners
- Flat seat board

Weight Limit

- 57kg/9 stone



Tilt-in-Space



Elevating Legrest




Manual



5 Year Frame Warranty



Made in Britain



Intelli-Gel
Integral



REFLEXION™
SEATING SYSTEM

Contact Us

Call Freephone 0800 212709 (UK only) and speak to our friendly Customer Service Team. E: info@kirtonhealthcare.co.uk

About Us

The Kirton Healthcare Group
23 Rookwood Way, Haverhill, Suffolk, CB9 8PB

T: 01440 705352 F: 01440 706199

E: info@kirtonhealthcare.co.uk

www.kirton-healthcare.co.uk

@Kirton_Seating